



ErgoTouch

Perfect combination of ergonomics and style.

The optimised recessed grips ensure comfortable, precise and effortless working. The special plastic surface and the additional handle waves offer a perfect grip in different holding positions.

PeriLux




Use of the colour codes.

Example:

ET588/4 - PK = pink
ET588/4 - WE = white
ET588/4 - SC = black

Place the colour code after the main article number for the desired colour variant.

Dental instrument to sever tooth-fixed ligaments. The extremely flat and sharp edged form of the working ends allows a very gentle penetration into the periodontal gap.



carl martin
 Solingen
 Finest Dental Instruments

NEUENKAMPER STR. 80 - 86 · D-42657 SOLINGEN · GERMANY
 TEL. +49(0)212/81 00 44 · FAX +49(0)212/8 73 45
 E-MAIL: INFO@CARLMARTIN.DE · INTERNET: WWW.CARLMARTIN.DE