

# Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **PlastiSept eco Wipes**  
Issue/Revision: 23.08.2023

Version: 3.1  
Replaces version: 3.0

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Trade name: **PlastiSept eco Wipes**  
UFI: not applicable

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses: Cleaning and disinfecting agent  
Intended purpose: Alcohol-free wet wipes for the cleaning and disinfection of sensitive surfaces of medical devices such as medical and dental units (incl. upholstery and suspension tables) and X-ray machines particularly with acrylic glass or metal surfaces.  
Uses advised against: None at intended use.  
Note: The product is intended for professional users.

### 1.3. Details of the supplier of the safety data sheet

Manufacturer/Supplier: ALPRO MEDICAL GMBH  
Mooswiesenstraße 9  
D-78112 St. Georgen (Germany)  
Telephone: +49 7725 9392-0  
Telefax: +49 7725 9392-91  
E-Mail: [info@alpro-medical.de](mailto:info@alpro-medical.de)  
Internet: [www.alpro-medical.de](http://www.alpro-medical.de)

E-mail address for the competent person responsible for the safety data sheet: [doku@alpro-medical.de](mailto:doku@alpro-medical.de)

### 1.4. Emergency telephone number

In-house emergency telephone number: +49 7725 9392-0  
Monday – Friday from 08:00 am to 04:30 pm (UTC+1);  
for chemical information and legal information on hazardous substances only  
Poison centre: +49 761 19240  
Poisoning information centre, Freiburg, Germany  
(24 h / 7 d), English is spoken  
National Poisons Information Service (UK): +44 344 892 0111  
National Poisons Information Service (NPIS)  
(24 h / 7 d), Medical Professionals Only

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

Classification in accordance with Regulation (EC) No 1272/2008 [CLP]

Classification	Classification procedure
Aquatic Chronic 3; H412	Calculation method

Full text of hazard classes as well as H-phrases: see under SECTION 16.1.

### 2.2. Label elements

Label elements in accordance with Regulation (EC) No 1272/2008 [CLP]

Hazard pictograms: -

# Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **PlastiSept eco Wipes**  
Issue/Revision: 23.08.2023

Version: 3.1  
Replaces version: 3.0

Signal word: -  
Hazard components for labelling: -  
H-phrases: H412 Harmful to aquatic life with long lasting effects.  
P-phrases: -

## 2.3. Other hazards

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII (see SECTION 12.5.).

The substances in the mixture have no endocrine disrupting properties according to Regulation (EC) No. 1907/2006, Annex XIV (see SECTION 11 and SECTION 12.6.).

The substances in the mixture are below the declaration limit for substances on the list of substances of very high concern according to Regulation (EC) No 1907/2006, Article 59, paragraph 10)

No further hazards known.

## SECTION 3: Composition/information on ingredients

### 3.1. Substances

Not applicable.

### 3.2. Mixtures

Chemical characterisation: Nonwoven wipes soaked with disinfecting solution. The disinfecting solution is a mixture of substances listed below with non-hazardous additions in aqueous solution.

#### Hazardous ingredients

Chemical name	Identification numbers	Classification in accordance with Regulation (EC) No 1272/2008	Weight %
N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine	CAS No: 2372-82-9 EC No: 219-145-8 REACH Registration No: 01-2119980592-29-XXXX	Acute Tox. 3; H301 Skin Corr. 1B; H314 STOT RE 2; H373 Aquatic Acute 1; H400 Aquatic Chronic 1; H410 <i>M-Factor acute: 10</i>	≥ 0.25 - < 1
Amines, N-C12-14-alkyltrimethylene-di-	CAS No: 90640-43-0 EC No: 292-562-0 REACH Registration No: 01-2119957843-25-XXXX	Acute Tox. 3; H301 Skin Corr. 1B; H314 STOT RE 1; H372 Aquatic Acute 1; H400 Aquatic Chronic 1; H410 <i>M-Factor acute: 100</i> <i>M-Factor chronic: 1</i>	≥ 0.1 - < 0.25

Full text of hazard classes and H-phrases: see SECTION 16.1.

Occupational exposure limits: see SECTION 8.1

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

General information: First aider: Pay attention to self-protection!

Following inhalation: Supply fresh air. In case of complaints consult a physician.

# Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **PlastiSept eco Wipes**  
Issue/Revision: 23.08.2023

Version: 3.1  
Replaces version: 3.0

- Following skin contact: By way of precaution, wash skin with water and soap. In case of skin reactions, consult a physician.
- Following eye contact: Flush eyes immediately with flowing water for 10 to 15 minutes holding eyelids apart. Remove contact lenses, if present and easy to do. Consult an ophthalmologist.
- Following ingestion: Rinse mouth with water. Let drink only a few sips of water (foaming product). Do not induce vomiting (risk of aspiration). Consult a physician immediately.

#### 4.2. Most important symptoms and effects, both acute and delayed

None known.

#### 4.3. Indication of any immediate medical attention and special treatment needed

No information available.

### SECTION 5: Firefighting measures

#### 5.1. Extinguishing media

Suitable extinguishing media: Water spray jet, alcohol resistant foam, extinguishing powder, carbon dioxide (CO<sub>2</sub>)

Unsuitable extinguishing media: Full water jet

#### 5.2. Special hazards arising from the substance or mixture

Hazardous combustion products: Carbon monoxide (CO), carbon dioxide (CO<sub>2</sub>), nitrogen oxides (NO<sub>x</sub>)

#### 5.3. Advice for firefighters

Special protective equipment: Not necessary

Further information: Cool endangered containers with water spray jet.

### SECTION 6: Accidental release measures

#### 6.1. Personal precautions, protective equipment and emergency procedures

##### For non-emergency personnel

Use personal protective equipment. See SECTION 8.2.

Avoid skin and eye contact.

##### For emergency responders

Use personal protective equipment. See SECTION 8.2.

#### 6.2. Environmental precautions

Do not discharge into drains or surface and ground water.

#### 6.3. Methods and material for containment and cleaning up

##### Containment

Not necessary

##### Cleaning up

Collect in suitable, closed containers for disposal.

##### Other information

Inappropriate containment and cleaning methods are not known.

**Safety data sheet**  
**according to Regulation (EC) No 1907/2006**



Trade name: **PlastiSept eco Wipes**  
 Issue/Revision: 23.08.2023

Version: 3.1  
 Replaces version: 3.0

**6.4. Reference to other sections**

Information on safe handling see SECTION 7.1.  
 Information on personal protective equipment see SECTION 8.2.  
 Information on disposal see SECTION 13.

**SECTION 7: Handling and storage**

**7.1. Precautions for safe handling**

**Precautions**

Avoid contact with skin and eyes.  
 Do not mix with other cleaners or disinfectants.

**Advice on general occupational hygiene**

When using do not eat, drink or smoke. Wash hands before breaks and at end of work. Keep away from medicines, food, feed, cosmetics and stimulants.

**7.2. Conditions for safe storage, including any incompatibilities**

Requirements for storage rooms and vessels: Keep only in the original container.  
 Advice on common storage: Not necessary  
 Further information on storage conditions: Not necessary  
 Storage class ([DE] TRGS 510): Do not store above 30 °C

**7.3. Specific end use(s)**

Apart from the uses mentioned in section 1.2 no other specific end uses are stipulated.

**Industry and sector specific guidance**

- [DE] TRGS 525 – Hazardous substances in medical care facilities (Section 7 Activities with disinfectants); Issue: September 2014;  
 Source: GMBI 2014 page 1294-1307 of 13.10.2014 [No 63]; 10.07.2015 [No 27]; www.baua.de
- [DE] DGUV rules 207-206 – Prevention of chemical risks when handling disinfectants in health service, Issue: 2016.12; Source: www.dguv.de/publikationen

**SECTION 8: Exposure controls/personal protection**

**8.1. Control parameters**

**Occupational exposure limits**

Country	Limit values				Legal Basis	Remarks
	Long term (8 hours)		Short term (15 minutes)			
	ppm	mg/m <sup>3</sup>	ppm	mg/m <sup>3</sup>		
N-(3-Aminopropyl)-N-dodecylpropan-1,3-diamin (CAS-Nr.: 2372-82-9)						
EU						no limit value specified
UK						no limit value specified

Used abbreviations, symbols, numerals and explanations in column „Remarks“

**Biological limit values**

Does not contain substances above concentration limits fixing a biological limit value.

# Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **PlastiSept eco Wipes**  
Issue/Revision: 23.08.2023

Version: 3.1  
Replaces version: 3.0

## Information on monitoring procedures

BS EN 482:2021-04; Title: Workplace exposure - Procedures for the determination of the concentration of chemical agents;  
British version of EN 482:2021

BS EN 689:2018; Title: Workplace atmospheres - Guidance for the assessment of exposure by inhalation to chemical agents for comparison with limit values and measurement strategy;  
British version of EN 689:2018

BS EN 14042:2003-04-24; Title: Workplace atmospheres - Guide for the application and use of procedures for the assessment of exposure to chemical and biological agents;  
British version of EN 14042:2003)

CEN/TR 17055:2017; Title: Workplace exposure. Measurement of chemical agents complying with the requirements given in EN 482 and either one of EN 838, EN 1076, EN 13205, EN 13890 and EN 13936. Choice of procedures

ISO TR 14294:2011; Title: Workplace atmospheres. Measurement of dermal Exposure. Principles and methods

## 8.2. Exposure controls

### Appropriate engineering controls

#### Technical and organisational protective measures

No special protective measures necessary.

### Personal protective equipment

Eye/face protection: Safety glasses with side protection according to BS EN 166.

Skin protection:

Hand protection: Protective gloves according to BS EN ISO 374-1 and BS EN 21420  
Splash guard:  
Protective gloves: type C; permeation-resistant at least 10 minutes  
Permanent contact (> 480 min):  
Protective gloves: type A or B; code letters: G, O, A;  
permeation-resistant at least 30 minutes

Other skin protection: Not necessary when used as intended.

Respiratory protection: Not necessary when used as intended.

Thermal hazards: No special protective measures necessary.

### Environmental exposure controls

Do not discharge into drains.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

(Information refer to the disinfecting solution)

Appearance/physical state: clear, colourless liquid

Odour: characteristic

Odour threshold: no data available

Melting point/freezing point: no data available

Initial boiling point and boiling range: no data available

# Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **PlastiSept eco Wipes**  
Issue/Revision: 23.08.2023

Version: 3.1  
Replaces version: 3.0

Flammability:	not applicable	
Lower explosive limit:	not applicable	
Upper explosive limit:	not applicable	
Flash point:	not applicable	
Auto-ignition temperature:	no data available	
pH (undiluted):	10.7 – 11.3	(20 °C)
Kinematic viscosity:	no data available	
Solubility in water:	completely soluble	
Partition coefficient: n-octanol/water	not applicable	
Vapour pressure:	no data available	(... °C)
Density:	0.99 – 1.01 g/cm <sup>3</sup>	(20 °C)
Relative vapour density:	no data available	
Particle characteristics:	not applicable	

## 9.2. Other information

### 9.2.1. Information with regard to physical hazard classes

Explosive substances/mixtures and products containing explosives:	not applicable.
Flammable gases:	not applicable
Aerosols:	not applicable
Oxidising gases:	not applicable
Gases under pressure:	not applicable
Flammable liquids:	not applicable
Flammable solids:	not applicable
Self-reactive substances and mixtures:	not applicable
Pyrophoric liquids:	not applicable
Pyrophoric solids:	not applicable
Self-heating substances and mixtures:	not applicable
Substances and mixtures, which emit flammable gases in contact with water:	not applicable
Oxidising liquids:	not applicable
Oxidising solids:	not applicable
Organic peroxides:	not applicable
Substances and mixtures corrosive to metals:	not applicable
Desensitised explosives:	not applicable

# Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **PlastiSept eco Wipes**  
Issue/Revision: 23.08.2023

Version: 3.1  
Replaces version: 3.0

## 9.2.2. Other safety characteristics

Electrical conductivity (undiluted): 130-170  $\mu\text{S}/\text{cm}$  (20 °C)

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

No hazardous reactions when handled and stored as intended.  
Do not mix with other cleaners or disinfectants.

### 10.2. Chemical stability

The product is stable when handled and stored as intended.

### 10.3. Possibility of hazardous reactions

None known

### 10.4. Conditions to avoid

None known

### 10.5. Incompatible materials

No material incompatibilities when used as directed.  
Prolonged exposure may cause flash rust on unalloyed steels. Upholstery surfaces may become dull.

### 10.6. Hazardous decomposition products

Does not decompose when used as intended.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

#### Acute toxicity

##### Product

Acute toxicity - oral:	Acute Toxicity Estimate $\text{ATE}_{\text{mix}} > 2000 \text{ mg/kg}$ => no classification
Acute toxicity - dermal:	Acute Toxicity Estimate $\text{ATE}_{\text{mix}} > 2000 \text{ mg/kg}$ => no classification
Acute toxicity - inhalation:	Acute Toxicity Estimate $\text{ATE}_{\text{mix}} > 20 \text{ mg/l}$ => no classification

##### Ingredients

#### Amines, N-C12-14-alkyltrimethylenedi- (CAS No: 90640-43-0):

Acute toxicity - oral:  $\text{LD}_{50}$ : 200 mg/kg; species: rat; method: OECD 423

#### N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine (CAS No: 2372-82-9):

Acute toxicity - oral:  $\text{LD}_{50}$ : 261 mg/kg; species: rat; method: OECD 401

Acute toxicity - dermal:  $\text{LD}_{50}$ : > 600 mg/kg; species: rat; method: OECD 402

#### Skin corrosion/irritation

##### Product

No classification. [calculation method]

**Safety data sheet**  
according to Regulation (EC) No 1907/2006



Trade name: **PlastiSept eco Wipes**  
Issue/Revision: 23.08.2023

Version: 3.1  
Replaces version: 3.0

**Serious eye damage/irritation**

Product

No classification. [calculation method]

**Respiratory or skin sensitisation**

Product

No data available.

**Germ cell mutagenicity**

Product

No data available.

**Carcinogenicity**

Product

No data available.

**Reproductive toxicity**

Product

No data available.

**STOT-single exposure**

Product

No classification. [calculation method]

**STOT-repeated exposure**

Product

No classification. [calculation method]

**Aspiration hazard**

Product

No data available.

**Information on other hazards**

Endocrine disrupting properties:

No substances are contained that have endocrine disrupting properties for humans.

**SECTION 12: Ecological information**

**12.1. Toxicity**

Harmful to aquatic life with long lasting effects. [calculation method]

**12.2. Persistence and degradability**

Biodegradability:

The disinfecting solution is biodegradable according to OECD criteria. The statement has been derived from the properties of the ingredients.

**12.3. Bioaccumulative potential**

No data available.

**12.4. Mobility in soil**

No data available.



# Safety data sheet

## according to Regulation (EC) No 1907/2006



Trade name: **PlastiSept eco Wipes**  
Issue/Revision: 23.08.2023

Version: 3.1  
Replaces version: 3.0

### 12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII.

### 12.6. Endocrine disrupting properties

No substances are contained that have endocrine disrupting properties for non-target organisms.

### 12.7. Other adverse effects

No data available.

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

#### Disposal of the product

Product residues must be disposed of as hazardous waste in compliance with the Directive 2008/98/EC on waste as well as national and regional regulations. Do not dispose of via the waste water. Leave product in the original container as possible. Do not mix with other waste materials.

Waste codes / waste designations according to EWC

Product residues: 15 02 02\* absorbents, filter materials (including oil filters not otherwise specified), wiping cloths, protective clothing contaminated by dangerous substances

#### Disposal of the packaging

Packaging contaminated with product is considered as hazardous waste and must be disposed of accordingly.

Waste codes / waste designations according to EWC

Contaminated packaging: 15 01 10\* packaging containing residues of or contaminated by hazardous substances

Recommendation

Contaminated packaging must be emptied optimally and can be recycled after appropriate cleaning (rinse with water).

## SECTION 14: Transport information

### 14.0. Transport classification

No dangerous good in sense of the transport regulations in road traffic (ADR), railway traffic (RID), inland waterway traffic (ADN), maritime traffic (IMDG-Code) and air traffic (ICAO-TI/IATA-DGR).

### 14.1. UN number

-

### 14.2. UN proper shipping name

-

### 14.3. Transport hazard class(es)

-

### 14.4. Packing group

-

**Safety data sheet**  
according to Regulation (EC) No 1907/2006



Trade name: **PlastiSept eco Wipes**  
Issue/Revision: 23.08.2023

Version: 3.1  
Replaces version: 3.0

**14.5. Environmental hazards**

-

**14.6. Special precautions for user**

-

**14.7. Maritime transport in bulk according to IMO instruments**

-

**14.8. Further information**

-

**SECTION 15: Regulatory information**

**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

**EU-Regulations**

REGULATION (EC) No 1005/2009 on substances that deplete the ozone layer  
not applicable

REGULATION (EC) No 2019/1021 on persistent organic pollutants and amending.  
not applicable

REGULATION (EU) No 649/2012 concerning the export and import of hazardous chemicals  
not applicable

REGULATION (EU) No 648/2004 on detergents  
Non-ionic surfactants: < 5%  
Disinfectants  
Perfumes

DIRECTIVE 2012/18/EU (Seveso III Directive) on the control of major-accident hazards involving dangerous substances, amending and subsequently repealing Council Directive 96/82/EC  
not applicable

DIRECTIVE 2010/75/EU on industrial emissions (integrated pollution prevention and control)  
VOC content: < 1 %

REACH – List of substances subject to authorisation (Annex XIV)  
not applicable

REACH – Restrictions on the manufacture, placing on the market and use of certain dangerous substances, mixtures and articles (Annex XVII)  
not applicable

COUNCIL DIRECTIVE 94/33/EC on the protection of young people at work  
not applicable

COUNCIL DIRECTIVE 92/85/EEC on the introduction of measures to encourage improvements in the safety and health at work of pregnant workers and workers who have recently given birth or are breastfeeding  
not applicable

**15.2. Chemical safety assessment**

For this mixture no chemical safety assessment has been carried out.

# Safety data sheet

according to Regulation (EC) No 1907/2006



Trade name: **PlastiSept eco Wipes**  
Issue/Revision: 23.08.2023

Version: 3.1  
Replaces version: 3.0

## SECTION 16: Other information

### 16.1. Full text of hazard classes and H-phrases

#### Hazard classes

Acute Tox.	Acute toxicity
Aquatic Acute	Acute aquatic hazard
Aquatic Chronic	Long-term aquatic hazard
Skin Corr.	Skin corrosion
STOT RE	Specific target organ toxicity (repeated exposure)

#### H-phrases (Hazard statements)

H301	Toxic if swallowed.
H314	Causes severe skin burns and eye damage.
H372	Causes damage to organs <or state all organs affected, if known> through prolonged or repeated exposure <state route of exposure if it is conclusively proven that no other routes of exposure cause the hazard>.
H373	May cause damage to organs <or state all organs affected, if known> through prolonged or repeated exposure <state route of exposure if it is conclusively proven that no other routes of exposure cause the hazard>.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects.
H412	Harmful to aquatic life with long lasting effects.

### 16.2. Abbreviations and acronyms

ADN	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par voie de <u>n</u> avigation intérieure (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par <u>r</u> oute (European Agreement concerning the International Carriage of Dangerous Goods by Road)
BS	<u>B</u> ritish <u>S</u> tandards
CAS	<u>C</u> hemical <u>A</u> bstracts <u>S</u> ervice
CLP	Regulation on <u>C</u> lassification, <u>L</u> abelling and <u>P</u> ackaging of Substances and Mixtures
[DE]	National German regulations
DGUV	<u>D</u> eutsche <u>G</u> esetzliche <u>U</u> nfall <u>v</u> ersicherung (English: German statutory accident insurance)
EC	<u>E</u> uropean <u>C</u> ommunity
EEC	<u>E</u> uropean <u>E</u> conomic <u>C</u> ommunity
EN	European Standard
EU	<u>E</u> uropean <u>U</u> nion
EWC	<u>E</u> uropean <u>W</u> aste <u>C</u> atalogue
GMBI	<u>G</u> emeinsames <u>M</u> inisterial <u>b</u> latt (English: Joint Ministerial Gazette)
IATA-DGR	<u>I</u> nternational <u>A</u> ir <u>T</u> ransport <u>A</u> ssociation - <u>D</u> angerous <u>G</u> oods <u>R</u> egulations
ICAO-TI	Technical Instructions For The Safe Transport of Dangerous Goods by Air
IMDG-Code	<u>I</u> nternational <u>M</u> aritime <u>C</u> ode for <u>D</u> angerous <u>G</u> oods
LD <sub>50</sub>	Median lethal dose
LGK	<u>L</u> ager <u>k</u> lasse (English: Storage class)
OECD	<u>O</u> rganization for <u>E</u> conomic <u>C</u> o-operation and <u>D</u> evelopment
PBT	<u>P</u> ersistent, <u>b</u> ioaccumulative and <u>t</u> oxic
REACH	<u>R</u> egistration, <u>E</u> valuation, <u>A</u> uthorisation and Restriction of <u>C</u> hemicals
RID	<u>R</u> èglement concernant le transport <u>I</u> nternational ferroviaire de marchandises <u>D</u> angereuses (Regulations Concerning the International Carriage of Dangerous Goods by Rail)
TRGS	<u>T</u> echnische <u>R</u> egeln für <u>G</u> efahr <u>s</u> toffe (English: Technical Rules for Hazardous Substances)
UN	<u>U</u> nited <u>N</u> ations
UTC	Coordinated Universal Time (French: Temps Universel Coordonné)
VOC	Volatile Organic Compounds

# Safety data sheet

## according to Regulation (EC) No 1907/2006



Trade name: **PlastiSept eco Wipes**  
Issue/Revision: 23.08.2023

Version: 3.1  
Replaces version: 3.0

vPvB      Very persistent and very bioaccumulative

### 16.3. Key literature references and sources for data

- Regulation (EC) No 1907/2006 (REACH), Annex II
- European Chemicals Agency (ECHA) – Guidance on the compilation of safety data sheets; Version 4.0 (December 2022); <https://echa.europa.eu/documents>
- Regulation (EC) No 1272/2008 (CLP regulation)
- European Chemicals Agency (ECHA) – Guidance on Labelling and Packaging in accordance with Regulation (EC) No 1272/2008; Version 4.2 (03/2021); <https://echa.europa.eu/documents>
- European Chemicals Agency (ECHA), Registered substances; <https://echa.europa.eu/information-on-chemicals/registered-substances>
- European Chemicals Agency (ECHA), C&L Classification and Labelling Inventory; <https://echa.europa.eu/information-on-chemicals/cl-inventory-database>
- Institute for Occupational Safety and Health of the German Social Accident Insurance (IFA): GESTIS database on hazardous substances and GESTIS - International limit values for chemical agents; <https://www.dguv.de/dguv/ifa/index.jsp>
- German Environmental Agency (Umweltbundesamt), Section IV 2.4: Office of Documentation and Information on Substances Hazardous to Waters RIGOLETTO (catalogue of Substances Hazardous to Waters); <https://webrigoletto.uba.de/rigoletto>

### 16.4. Methods according to Article 9 of Regulation (EC) No 1272/2008 for the evaluation of information for classification purposes

Calculation method according to the criteria in Annex I 1272/2008.

pH value measurement.

Material compatibility and corrosiveness in practical tests

### 16.5. Training advice

Provide adequate information, instructions and training for users.

### 16.6. Indication of changes

A dash in the left hand margin indicates an amendment from the previous version.

---

The information given in the safety data sheet only applies to the described product in connection with its intended use. This information is based on the latest state of our knowledge at the time of revision. In particular, it describes our product under the aspect of its hazards and safety measures to be taken. It does not constitute any guarantee of product properties and quality features.

---