

ivoclar

Small.
But mighty.



PrograMill®
Dry

The masterpiece
in milling excellence



Making People Smile

Compact milling. Big impact.

Simple. Fast. Efficient. PrograMill® Dry is a new 5-axis dry milling machine for zirconia and PMMA, ready for the laboratory of tomorrow.

Combining maximum performance with a compact footprint, it delivers exceptional ease of use and precision, making high-quality dry milling more accessible than ever.



PrograMill® Dry

The masterpiece in milling excellence

Every aspect of PrograMill Dry has been designed to support seamless workflows and intuitive operation. The result is greater comfort, efficiency and performance.

With its compact footprint, PrograMill Dry is the ideal solution for dental laboratories and practice laboratories with limited space – ready for a new era in milling excellence.



Simply digital. Smart by design.

PrograMill Dry makes it easy and cost-effective to get started with digital milling, while delivering reliable, precise results.

Its smart features, powered by the latest version of PrograMill CAM software, support efficient workflows and intuitive operation.

5-axis system

Ionizer

PrograMill Dry

ivoclar

9,2-inch touchscreen with
step-by-step guidance

RFID-based material management
with programmable tags

Advanced jet stream
extraction system

Zero point clamping system
for high-accuracy milling

Milling chamber with a smart,
user-friendly design and an
extra-large door

Everything made easy

PrograMill Dry has been designed with one guiding principle in mind: making every step as simple as possible. Easy installation, intuitive operation and self-repair guides help you to work more independently, efficiently and quickly.



GO TO IVOCALAR
3D-SHOWROOM!

Open space. Excellent access.

PrograMill Dry is a true space-saving solution with a compact footprint that fits easily into laboratories with limited space. At the same time, its extra-large door provides excellent visibility and easy access to the milling chamber, simplifying everyday operation and maintenance.

Ease of use is further enhanced by the integrated 9.2-inch touchscreen, which provides intuitive, step-by-step guidance throughout the milling process. This enables even new users to operate the machine with confidence.

Clear

- ✓ Extra-large door for excellent visibility and easy access
- ✓ Round milling chamber for easy cleaning

Intuitive

- ✓ 9.2-inch display with an intuitive step-by-step user interface

Precise

- ✓ Zero point clamping system and optimized CAM strategies deliver consistently precise results with minimal rework

Simply. Spotless.

The smart design of the milling chamber, combined with the integrated ionizer and new jet stream extraction system, significantly reduces cleaning effort.

This saves valuable time during daily maintenance and makes PrograMill Dry the ideal choice for PMMA processing.



PrograMill® Dry Base Smart storage. Clean integration.

Like all Ivoclar milling machines, PrograMill® Dry is available with a dedicated base unit. PrograMill® Dry Base features a low-noise suction unit, as well as an integrated drawer with a specially designed insert to keep materials, tools and accessories neatly organized and within easy reach.

With all the essentials.
At the forefront.

PrograMill Dry is the ideal choice for efficient, flexible and cost-effective milling. Its low weight and well-thought-out design make it easy to transport, install and maintain – wherever you are.



Expertise.
For professionals.

Join the CAM Academy: Whatever your goals, our trainers will guide you through customized CAM strategies and setup options. This enables you to fabricate exactly what you need, the way you want – and to master even complex cases with confidence.

Expand your expertise in both theory and practice with hands-on training, live demonstrations and real-world case studies using the PrograMill® portfolio and PrograMill CAM software.



