



Simply connect and get started!

YOUNG PROXEO TWIST

wireless handpiece with foot pedal

The Young Proxeo Twist Cordless enables rotary polishing without any restrictions – in more ways than one. **Complete cordless freedom of movement** ensures maximum comfort during use. The handpiece has a cordless foot control and is compatible with two different cup systems: Young polishing cups with screw connection or hygienic disposable contra-angles.

Recommended rotation speeds
1.250 - 2.500 U/min



Versatile
easy to move from room to room

Cordless and light-weight

Modular
Choose from various cup systems



Well balanced
Reduces hand fatigue

Reliable
Li-ion battery lasts all day

Your Advantages:

Ideal rotation speeds protect the tooth from overheating and prevent abrasion. They also ensure a comfortable treatment and reduce splashing of the polishing paste.

The **wireless foot control** allows maximum freedom of movement and quick position changes around the patient. Easy to operate thanks to WiFi connection. The speed can be adjusted continuously from 0 to 2,500 rpm.

Ergonomic treatments are made possible by **draught-free working**. No cable – no draught on the device, resulting in less hand fatigue and a significantly reduced risk of developing carpal tunnel syndrome. The cordless handpiece also enables a back-friendly working posture. **Also suitable for left-handers**.

| Item Number | Description |
|-------------|--|
| YI-PL-64SET | Young Proxeo Twist Cordless Set, consisting of handpiece drive, handpiece holder, charger incl. adaptor, foot control and Screw-Type attachment YI-PL-64 |
| YI-PL-44SET | Young Proxeo Twist Cordless Set, consisting of handpiece drive, handpiece holder, charger incl. adaptor, foot control and attachment YI-PL-44 for disposable prophy angles (Doriot system) |
| YI-PL-64 | Young Proxeo Twist Cordless Attachment for Screw-Type prophy cups |
| YI-PL-44 | Young Proxeo Twist Cordless Attachment for disposable prophy angles, (Doriot system) |

The following products are included in the kits (YI-PL-64SET & YI-PL-44SET)



Discover all the advantages now



YI-GB-PTCH/1125