

Sicherheitsdatenblatt Repin (white)

ABSCHNITT 1: BEZEICHNUNG DES STOFFS BEZIEHUNGSWEISE DES GEMISCHS UND DES UNTERNEHMENS

1.1. Produktidentifikator

Produktname : Repin (white)

1.2. Relevante identifizierte Verwendungen des Stoffs oder Gemischs und Verwendungen, von denen abgeraten wird

Relevante identifizierte Verwendungen

Hauptverwendungskategorie : Professionelle Anwendung
 Funktions- oder Verwendungskategorie : Material ist für den Einsatz im Dentalbereich.

Verwendungen, von denen abgeraten wird

Keine weiteren Informationen verfügbar

1.3. Einzelheiten zum Lieferanten, der das Sicherheitsdatenblatt bereitstellt

Hersteller

SpofaDental a.s.
 Markova 238
 CZ-506 46 Jicin Czech Republic
 T +420 493 583 204

Ansprechpartner : safety@kavokerr.com - tel. +41 91 610 06 00 (08.00-17.00)

1.4. Notrufnummer

Notrufnummer : CHEMTREC® Emergency Call Center. Emergency Telephone Number (for USA only) 001-800-424-9300 International and Maritime Telephone Number +1 (703) 527-3887

Land	Organisation/Firma	Anschrift	Notrufnummer
Belgien	Centre Anti-Poisons/Antigifocentrum c/o Hôpital Central de la Base - Reine Astrid	Rue Bruyn 1 1120 Bruxelles/Brussel	+32 70 245 245
Deutschland	Clinical Toxicology and Berlin Poison Information Centre www.giftnotruf.de	Institute of Toxicology, Oranienburger Str 285, Berlin	+49 30 192 40 +49 30 3068 6711
Luxemburg	Centre Anti-Poisons/Antigifocentrum c/o Hôpital Central de la Base - Reine Astrid	Rue Bruyn 1 1120 Bruxelles/Brussel	+352 8002 5500

ABSCHNITT 2: MÖGLICHE GEFAHREN

2.1. Einstufung des Stoffs oder Gemischs

Einstufung gemäß Verordnung (EG) Nr. 1272/2008 [CLP]

Aquatic Acute 1 H400
 Aquatic Chronic 1 H410

Volltext der Gefahrenklassen und Gefahrenhinweise: siehe Kapitel 16

2.2. Kennzeichnungselemente

Kennzeichnung gemäß Verordnung (EG) Nr. 1272/2008 [CLP]

Gefahrenpiktogramme (CLP) :



GHS09

Signalwort (CLP) : Achtung
 Gefahrenhinweise (CLP) : H410 - Sehr giftig für Wasserorganismen mit langfristiger Wirkung
 Sicherheitshinweise (CLP) : P273 - Freisetzung in die Umwelt vermeiden

Zusätzliche Sätze

P391 - Verschüttete Mengen aufnehmen
P501 - Inhalt/Behälter einer Sammelstelle für gefährliche Abfälle oder Sondermüll, gemäß den lokalen, regionalen, nationalen und/oder internationalen Vorschriften zuführen

: Das Produkt wird als Medizinprodukt betrachtet und unterliegt daher nicht der Kennzeichnung (EU-Verordnung 1272/2008, Artikel 1, Absatz 5d).

2.3. Sonstige Gefahren

Weitere Gefahren ohne Einfluss auf die Einstufung : Unter normalen Umständen kein(e).

Dieser Stoff/Gemisch erfüllt nicht die PBT-Kriterien der REACH-Verordnung, Annex XIII.

Dieser Stoff/Gemisch erfüllt nicht die vPvB-Kriterien der REACH-Verordnung, Annex XIII.

ABSCHNITT 3: ZUSAMMENSETZUNG/ANGABEN ZU BESTANDTEILEN

3.1. Stoffe

Nicht anwendbar

3.2. Gemische

Name	Produktidentifikator	%	Einstufung gemäß Verordnung (EG) Nr. 1272/2008 [CLP]
Zinkoxid	(CAS-Nr.) 1314-13-2 (EG-Nr.) 215-222-5 (EG Index-Nr.) 030-013-00-7 (REACH-Nr) 01-2119463881-32	50 - 75	Aquatic Acute 1, H400 Aquatic Chronic 1, H410

Wortlaut der H-Sätze: siehe unter Abschnitt 16

ABSCHNITT 4: ERSTE-HILFE-MAßNAHMEN

4.1. Beschreibung der Erste-Hilfe-Maßnahmen

Erste-Hilfe-Maßnahmen allgemein : Keine besonderen/spezifischen Maßnahmen erforderlich.

Erste-Hilfe-Maßnahmen nach Einatmen : Die betroffene Person an die frische Luft bringen und in einer Position ruhigstellen, die das Atmen erleichtert. Bei Unwohlsein ärztlichen Rat einholen/ärztliche Hilfe hinzuziehen.

Erste-Hilfe-Maßnahmen nach Hautkontakt : Behutsam mit viel Wasser und Seife waschen. Bei Hautreizung: Ärztlichen Rat einholen/ärztliche Hilfe hinzuziehen.

Erste-Hilfe-Maßnahmen nach Augenkontakt : BEI KONTAKT MIT DEN AUGEN: Einige Minuten lang behutsam mit Wasser spülen. Eventuell vorhandene Kontaktlinsen nach Möglichkeit entfernen. Weiter spülen. Bei Unwohlsein ärztlichen Rat einholen/ärztliche Hilfe hinzuziehen.

Erste-Hilfe-Maßnahmen nach Verschlucken : Nach Verschlucken, Mund mit Wasser ausspülen (nur wenn die betroffene Person bei Bewusstsein ist). Bei Unwohlsein GIFTINFORMATIONSZENTRUM oder Arzt anrufen.

4.2. Wichtigste akute und verzögert auftretende Symptome und Wirkungen

Symptome/Schäden : Bei üblichen Gebrauchsbedingungen keine nennenswerte Gefährdung zu erwarten.

4.3. Hinweise auf ärztliche Soforthilfe oder Spezialbehandlung

Keine spezifischen Maßnahmen festgestellt.

ABSCHNITT 5: MAßNAHMEN ZUR BRANDBEKÄMPFUNG

5.1. Löschmittel

Geeignete Löschmittel : Für Umgebungsbrände geeignete Löschmittel verwenden. Schaum, Kohlendioxid (CO₂) und Trockenlöschpulver.

Ungünstige Löschmittel : Keinen starken Wasserstrahl benutzen.

5.2. Besondere vom Stoff oder Gemisch ausgehende Gefahren

Brandgefahr : Nicht brennbar.

Gefährliche Zerfallsprodukte im Brandfall : Kohlenstoffoxide (CO, CO₂). Zinkoxid.

5.3. Hinweise für die Brandbekämpfung

Löschanweisungen : Zur Kühlung exponierter Behälter Wassersprühstrahl oder -nebel benutzen. Eindringen von Löschwasser in die Umwelt vermeiden (verhindern).

Schutz bei der Brandbekämpfung : Brandabschnitt nicht ohne ausreichende Schutzausrüstung, einschließlich Atemschutz betreten.

ABSCHNITT 6: MAßNAHMEN BEI UNBEABSICHTIGTER FREISETZUNG

6.1. Personenbezogene Vorsichtsmaßnahmen, Schutzausrüstungen und in Notfällen anzuwendende Verfahren

Allgemeine Maßnahmen : Berührung mit den Augen und der Haut vermeiden.

Nicht für Notfälle geschultes Personal

Schutzrüstung : Siehe Abschnitt 8.

Einsatzkräfte

Keine weiteren Informationen verfügbar

6.2. Umweltschutzmaßnahmen

Die Einleitung in Flüsse oder Kanalisation ist verboten.

6.3. Methoden und Material für Rückhaltung und Reinigung

Zur Rückhaltung : Abfälle in geeigneten und gekennzeichneten Behältern sammeln und unter Beachtung der örtlichen Gesetze entsorgen.

Reinigungsverfahren : Das Produkt mechanisch aufnehmen. Mit viel Wasser abwaschen.

6.4. Verweis auf andere Abschnitte

Weitere Angaben zur Entsorgung siehe Abschnitt 13.

ABSCHNITT 7: HANDHABUNG UND LAGERUNG**7.1. Schutzmaßnahmen zur sicheren Handhabung**

Schutzmaßnahmen zur sicheren Handhabung : Berührung mit den Augen und der Haut vermeiden. Vorgeschriebene persönliche Schutzrüstung verwenden.

Hygienemaßnahmen : Vor dem Essen, Trinken oder Rauchen und beim Verlassen des Arbeitsplatzes die Hände und andere exponierte Körperstellen mit milder Seife und Wasser waschen. Bei Verwendung dieses Produkts nicht essen, trinken oder rauchen.

7.2. Bedingungen zur sicheren Lagerung unter Berücksichtigung von Unverträglichkeiten

Lagerbedingungen : An einem trockenen Ort aufbewahren. In einem geschlossenen Behälter aufbewahren. In der Originalverpackung aufbewahren.

Lager : An einem gut belüfteten Ort aufbewahren.

7.3. Spezifische Endanwendungen

Zusätzliche Informationen beim Lieferanten erfragen.

ABSCHNITT 8: BEGRENZUNG UND ÜBERWACHUNG DER EXPOSITION/PERSÖNLICHE SCHUTZAUSRÜSTUNGEN**8.1. Zu überwachende Parameter**

Keine weiteren Informationen verfügbar

8.2. Begrenzung und Überwachung der Exposition

Geeignete technische Steuerungseinrichtungen : Für eine gute Belüftung des Arbeitsplatzes sorgen. Augen-Notduschen sollten in unmittelbarer Nähe einer möglichen Exposition verfügbar sein.

Persönliche Schutzrüstung : Handschuhe.

Handschutz : Geeignete Schutzhandschuhe tragen. Nitrilkautschukhandschuhe. Materialdicke: 0,09mm. Durchbruchzeit: >480 min. STANDARD EN 374.

Augenschutz : Spritzschutzbrille tragen, wenn Augenkontakt durch Verspritzen möglich ist. STANDARD EN 166.

Haut- und Körperschutz : Bei der Arbeit geeignete Schutzkleidung tragen

Atemschutz : Atemschutzmaske nicht erforderlich.



Sonstige Angaben

: Während der Verwendung nicht essen, trinken oder rauchen. Vor dem Essen, Trinken oder Rauchen und beim Verlassen des Arbeitsplatzes die Hände und andere exponierte Körperstellen mit milder Seife und Wasser waschen. Die persönliche Schutzrüstung muss in Übereinstimmung mit den geltenden CEN-Normen und nach Absprache mit dem Lieferanten der persönlichen Schutzrüstung gewählt werden.

ABSCHNITT 9: PHYSIKALISCHE UND CHEMISCHE EIGENSCHAFTEN**9.1. Angaben zu den grundlegenden physikalischen und chemischen Eigenschaften**

Aggregatzustand : Flüssigkeit

Aussehen : Paste.

Farbe : Weiß.

Geruch : Nach Minze.

Geruchsschwelle	: Keine Daten verfügbar
pH-Wert	: Keine Daten verfügbar
Verdunstungsgrad (Butylacetat=1)	: Keine Daten verfügbar
Schmelzpunkt	: Keine Daten verfügbar
Gefrierpunkt	: Keine Daten verfügbar
Siedepunkt	: Keine Daten verfügbar
Flammpunkt	: Keine Daten verfügbar
Selbstentzündungstemperatur	: Keine Daten verfügbar
Zersetzungstemperatur	: Keine Daten verfügbar
Entzündlichkeit (fest, gasförmig)	: Keine Daten verfügbar
Dampfdruck	: Keine Daten verfügbar
Relative Dampfdichte bei 20 °C	: Keine Daten verfügbar
Relative Dichte	: Keine Daten verfügbar
Löslichkeit	: Keine Daten verfügbar
Log Pow	: Keine Daten verfügbar
Viskosität, kinematisch	: Keine Daten verfügbar
Viskosität, dynamisch	: Keine Daten verfügbar
Explosive Eigenschaften	: Das Produkt ist nicht explosionsgefährlich.
Brandfördernde Eigenschaften	: Nicht brennbar.
Explosionsgrenzen	: Keine Daten verfügbar

9.2. Sonstige Angaben

Zusätzliche Hinweise : Nach unserer Kenntnis, keine

ABSCHNITT 10: STABILITÄT UND REAKTIVITÄT**10.1. Reaktivität**

Das Produkt ist nicht reaktiv unter normalen Gebrauchs-, Lagerungs- und Transportbedingungen.

10.2. Chemische Stabilität

Stabil unter normalen Bedingungen.

10.3. Möglichkeit gefährlicher Reaktionen

Keine Polymerisation.

10.4. Zu vermeidende Bedingungen

Nach unserer Kenntnis, keine.

10.5. Unverträgliche Materialien

Nach unserer Kenntnis, keine.

10.6. Gefährliche Zersetzungsprodukte

Keine Zersetzung bei normaler Verwendung und Lagerung.

ABSCHNITT 11: TOXIKOLOGISCHE ANGABEN**11.1. Angaben zu toxikologischen Wirkungen**

Akute Toxizität : Nichteingestuft

Zinkoxid (1314-13-2)	
LD50 oral Ratte	> 5000 mg/kg
LD50 Dermal Ratte	> 2000 mg/kg
LC50 Inhalation Ratte (mg/l)	0,4 mg/l/4h

Ätz-/Reizwirkung auf die Haut	: Nichteingestuft
Schwere Augenschädigung/-reizung	: Nichteingestuft
Sensibilisierung der Atemwege/Haut	: Nichteingestuft
Keimzell-Mutagenität	: Nichteingestuft
Karzinogenität	: Nichteingestuft
Reproduktionstoxizität	: Nichteingestuft
Spezifische Zielorgan-Toxizität bei einmaliger Exposition	: Nichteingestuft
Spezifische Zielorgan-Toxizität bei wiederholter Exposition	: Nichteingestuft
Aspirationsgefahr	: Nichteingestuft

Mögliche schädliche Wirkungen auf den Menschen und mögliche Symptome : Bei üblichen Verwendungsbedingungen wurden keine gesundheitsgefährdenden Wirkungen festgestellt.

ABSCHNITT 12: UMWELTBEZOGENE ANGABEN

12.1. Toxizität

Ökologie - Wasser : Sehr giftig für Wasserorganismen mit langfristiger Wirkung.

Zinkoxid (1314-13-2)	
LC50 Fische 1	1,1 mg/l (96 Stunden -Regenboreforelle)
EC50 Daphnia 1	24,6 mg/l (48 Stunden - Daphnia magna)

12.2. Persistenz und Abbaubarkeit

Repin (white)	
Persistenz und Abbaubarkeit	Kann längerfristig schädliche Wirkungen auf die Umwelt haben.

12.3. Bioakkumulationspotenzial

Repin (white)	
Bioakkumulationspotenzial	Keine Angaben.
Zinkoxid (1314-13-2)	
Log Pow	< 0

12.4. Mobilität im Boden

Repin (white)	
Ökologie - Boden	Keine Daten verfügbar.

12.5. Ergebnisse der PBT- und vPvB-Beurteilung

Repin (white)	
Dieser Stoff/Gemisch erfüllt nicht die PBT-Kriterien der REACH-Verordnung, Annex XIII.	
Dieser Stoff/Gemisch erfüllt nicht die vPvB-Kriterien der REACH-Verordnung, Annex XIII.	

12.6. Andere schädliche Wirkungen

Andere schädliche Wirkungen : Nach unserer Kenntnis, keine.
Zusätzliche Hinweise : Keine weiteren Auswirkungen bekannt

ABSCHNITT 13: HINWEISE ZUR ENTSORGUNG

13.1. Verfahren der Abfallbehandlung

Örtliche Vorschriften (Abfall) : Als gefährlichen Abfall entsorgen.
Verfahren der Abfallbehandlung : Inhalt/Behälter gemäß den Sortieranweisungen des zugelassenen Einsammlers entsorgen
Empfehlungen für die Produkt-/Verpackung-Abfallentsorgung : Inhalt/Behälter einer Sammelstelle für gefährliche Abfälle oder Sonderabfälle zuführen. Freisetzung in die Umwelt vermeiden.
EAK-Code : 18 01 06* - Chemikalien, die aus gefährlichen Stoffen bestehen oder solche enthalten

ABSCHNITT 14: ANGABEN ZUM TRANSPORT

Entsprechend den Anforderungen von ADR / RID / IMDG / IATA / ADN

ADR	IMDG	IATA	RID
14.1. UN-Nummer			
3082	3082	3082	3082
14.2. Ordnungsgemäße UN-Versandbezeichnung			
UMWELTGEFÄHRDENDER STOFF, FLÜSSIG, N.A.G.	ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S.	Environmentally hazardous substance, liquid, n.o.s.	UMWELTGEFÄHRDENDER STOFF, FLÜSSIG, N.A.G.
Eintragung in das Beförderungspapier			
UN 3082 UMWELTGEFÄHRDENDER STOFF, FLÜSSIG, N.A.G. (Zinkoxid(1314-13-2)), 9, III, (E)	UN 3082 ENVIRONMENTALLY HAZARDOUS SUBSTANCE, LIQUID, N.O.S., 9, III, MARINE POLLUTANT		
14.3. Transportgefahrenklassen			
9	9	9	9

ADR	IMDG	IATA	RID
14.4. Verpackungsgruppe			
III	III	III	III
14.5. Umweltgefahren			
Umweltgefährlich : Ja	Umweltgefährlich : Ja Meeresschadstoff : Ja	Umweltgefährlich : Ja	Umweltgefährlich : Ja
Keine zusätzlichen Informationen verfügbar			

14.6. Besondere Vorsichtsmaßnahmen für den Verwender**- Landtransport**

Klassifizierungscode (ADR)	: M6
Sonderbestimmung (ADR)	: 274, 335, 601, 375
Begrenzte Mengen (ADR)	: 5L
Freigestellte Mengen (ADR)	: E1
Verpackungsanweisungen (ADR)	: P001, IBC03, LP01, R001
Sondervorschriften für die Verpackung (ADR)	: PP1
Sondervorschriften für die Zusammenpackung (ADR)	: MP19
Anweisungen für Tankfahrzeuge und Schüttgutcontainer (ADR)	: T4
Besondere Bestimmungen für Tankfahrzeuge und Schüttgutcontainer (ADR)	: TP1, TP29
Tankcodierung (ADR)	: LGBV
Tanktransportfahrzeug	: AT
Beförderungskategorie (ADR)	: 3
Besondere Beförderungsbestimmungen - Pakete (ADR)	: V12
Besondere Bestimmungen für die Beförderung - Be-, Entladen und Handhabung (ADR)	: CV13
Gefahr-Nr. (Kemlerzahl)	: 90
Orangefarbene Tafeln	:



Tunnelbeschränkungscode (ADR)	: E
-------------------------------	-----

- Seeschifftransport

Sonderbestimmung (IMDG)	: 274, 335, 969
Begrenzte Mengen (IMDG)	: 5L
Freigestellte Mengen (IMDG)	: E1
Verpackungsanweisungen (IMDG)	: P001, LP01
Sondervorschriften für die Verpackung (IMDG)	: PP1
IBC-Verpackungsanweisungen (IMDG)	: IBC03
Tankanweisungen (IMDG)	: T4
Besondere Bestimmungen für Tanks (IMDG)	: TP2, TP29
EmS-Nr. (Brand)	: F-A
EmS-Nr. (Unbeabsichtigte Freisetzung)	: S-F
Ladungskategorie (IMDG)	: A

- Lufttransport

PCA freigestellte Mengen (IATA)	: E1
PCA begrenzte Mengen (IATA)	: Y964
PCA begrenzte max. Nettomenge (IATA)	: 30kgG
PCA Verpackungsvorschriften (IATA)	: 964
Max. PCA Nettomenge (IATA)	: 450L
CAO Verpackungsvorschriften (IATA)	: 964
Max. CAO Nettomenge (IATA)	: 450L
Sonderbestimmung (IATA)	: A97, A158, A197
ERG-Code (IATA)	: 9L

Bahntransport

Keine Daten verfügbar

14.7. Massengutbeförderung gemäß Anhang II des MARPOL-Übereinkommens 73/78 und gemäß IBC-Code

Nicht anwendbar

ABSCHNITT 15: RECHTSVORSCHRIFTEN**15.1. Vorschriften zu Sicherheit, Gesundheits- und Umweltschutz/spezifische Rechtsvorschriften für den Stoff oder das Gemisch****EU-Verordnungen**

Enthält keinen Stoff, der den Beschränkungen von Anhang XVII der REACH-Verordnung unterliegt

Enthält keinen REACH-Kandidatenstoff

Enthält keinen in REACH-Anhang XIV gelisteten Stoff

Nationale Vorschriften

Verordnung 2015/830 /EC (CLP), 1907/2006/EC (REACH), 1272/2008/EC, 790/2009/EC. Landtransport ADR/RID und GGVS/GGVE (grenzüberschreitend/Inland), Seeschifftransport IMDG/GGVSee, Lufttransport ICAO-TI und IATA-DGR, Grenzwerte am Arbeitsplatz, MAK-Werte. EG-abfallkatalog.

Deutschland

VwVwS, Verweis auf Anhang : Wassergefährdungsklasse (WGK) 2, Wassergefährdend (Einstufung nach VwVwS, Anhang 4)

Störfall-Verordnung - 12. BImSchV : Unterliegt nicht der 12. BImSchV (Bundes-Immissionsschutzverordnung) (Störfall-Verordnung)

15.2. Stoffsicherheitsbeurteilung

Für die Stoffe oder Mischung wurde keine Stoffsicherheitsbeurteilung durch den Lieferanten durchgeführt

ABSCHNITT 16: SONSTIGE ANGABEN

Änderungshinweise:

1.3	Kontakt Informationen	Aktualisiert
Datenquellen	: Verordnung 2015/830 /EC (CLP), 1907/2006/EC (REACH), 1272/2008/EC, 790/2009/EC. Landtransport ADR/RID und GGVS/GGVE (grenzüberschreitend/Inland), Seeschifftransport IMDG/GGVSee, Lufttransport ICAO-TI und IATA-DGR, Grenzwerte am Arbeitsplatz, MAK-Werte. EG-abfallkatalog. Einstufung gemäß Verordnung (EG) Nr. 1272/2008 [CLP].	
Ausgabedatum	: 29.06.2016	
Überarbeitungsdatum	: 19.02.2020	
Ersetzt	: 03.01.2017	
Datum der totalrevision	: 29.06.2016	
Version	: 1.2	

Vollständiger Wortlaut der H- und EUH-Sätze:

Aquatic Acute 1	Akut gewässergefährdend, Kategorie 1
Aquatic Chronic 1	Chronisch gewässergefährdend, Kategorie 1
H400	Sehr giftig für Wasserorganismen
H410	Sehr giftig für Wasserorganismen mit langfristiger Wirkung

Die Angaben dieses Datenblatts werden in Anbetracht der gegenwärtigen Kenntnisse und Erfahrungen als korrekt angesehen, es kann jedoch keine Vollständigkeitsgarantie hinsichtlich der Informationen gewährleistet werden. Deswegen liegt es im Interesse des Verbrauchers, Sicherheit darüber zu erhalten, dass die Angaben in Bezug auf den vorgesehenen Anwendungsbereich ausreichen.

Repin (brown)

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878
Issue date: 29/06/2016 Revision date: 09/10/2022 Supersedes version of: 04/01/2020 Version: 2.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture
Product name : Repin (brown)

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses

Main use category : Professional use
Use of the substance/mixture : Preparation intended for dental medical use
Function or use category : Dental materials.

Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

SpofaDental a.s
Markova 238
CZ-506 01 Jičín - Czech Republic
T +420 493 583 204 / Information Phone Number: 1-800-841-1428 (Customer Service)
safety@envistaco.com

1.4. Emergency telephone number

Emergency number : CHEMTREC® Emergency Call Center. Emergency Telephone Number (for USA only) 001-800-424-9300
International and Maritime Telephone Number +1 (703) 527-3887

Country	Official advisory body	Address	Emergency number	Comment
Gibraltar	GHA Call Centre Zone 2, Level3, St Bernard's Hospital	Harbour Views Road	+350 200 79700 +350 200 72266	
Ireland	National Poisons Information Centre Beaumont Hospital	PO Box 1297 Beaumont Road 9 Dublin	+353 1 809 2566 (Healthcare professionals-24/7) +353 1 809 2166 (public, 8am - 10pm, 7/7)	
Malta	Medicines & Poisons Info Office	Mater Dei Hospital MSD 2090 Msida	+356 2545 6508	
United Kingdom	National Poisons Information Service (Newcastle Unit)	Claremont Place Newcastle-upon-Tyne, Newcastle	+44 191 2606182 +44 191 2606180	Hours of operation: 24hrs

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Eye Irrit. 2 H319
Skin Sens. 1 H317

Full text of hazard classes and H-statements : see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :



GHS07

Safety Data Sheet

Repin (brown)

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

Signal word (CLP)	: Warning
Contains	: Clove bud oil, rosin; colophony
Hazard statements (CLP)	: H317 - May cause an allergic skin reaction. H319 - Causes serious eye irritation.
Precautionary statements (CLP)	: P261 - Avoid breathing spray, vapours, fume, mist. P280 - Wear eye protection, protective gloves. P305+P351+P338+P310 - IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor. P337+P313 - If eye irritation persists: Get medical advice/attention. P333+P313 - If skin irritation or rash occurs: Get medical advice/attention. P362+P364 - Take off contaminated clothing and wash it before reuse.
Extra phrases	: The product is seen as a medical device (Directive 90/385/EEC, 93/42/EEC, 98/79/EC) and therefore not subject to labelling (EU-regulation 1272/2008, article 1, paragraph 5d). A safety data sheet is not required for this product under Article 31 of REACH. This Product Safety Information Sheet has been created on a voluntary basis.

2.3. Other hazards

Other hazards which do not result in classification : None under normal conditions.

This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII

This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

The mixture does not contain substance(s) included in the list established in accordance with Article 59(1) of REACH for having endocrine disrupting properties, or is not identified as having endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at a concentration equal to or greater than 0,1 %

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	Conc. (% w/w)	Classification according to Regulation (EC) No. 1272/2008 [CLP]
rosin; colophony	(CAS-No.) 8050-09-7 (EC-No.) 232-475-7 (EC Index-No.) 650-015-00-7 (REACH-no) 01-2119480418-32	50 – 75	Skin Sens. 1, H317
Clove bud oil	(CAS-No.) 8000-34-8 (EC-No.) 616-772-2	< 20	Eye Irrit. 2, H319 Skin Sens. 1, H317

Full text of H- and EUH-statements: see section 16

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general	: In all cases of doubt, or when symptoms persist, seek medical attention. Never give anything by mouth to an unconscious person. If you feel unwell, seek medical advice (show the label where possible).
First-aid measures after inhalation	: Remove person to fresh air and keep comfortable for breathing. Get medical advice/attention if you feel unwell.
First-aid measures after skin contact	: Gently wash with plenty of soap and water. If skin irritation occurs: Get medical advice/attention.
First-aid measures after eye contact	: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.
First-aid measures after ingestion	: If swallowed, rinse mouth with water (only if the person is conscious). Call a POISON CENTER/doctor if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after skin contact	: May cause an allergic skin reaction.
Symptoms/effects after eye contact	: Causes serious eye irritation.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically. In all cases of doubt, or when symptoms persist, seek medical attention.

Safety Data Sheet

Repin (brown)

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Use extinguishing media appropriate for surrounding fire. Foam, carbon dioxide (CO₂) and powder.
Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

Fire hazard : Non flammable.
Explosion hazard : Product is not explosive.
Hazardous decomposition products in case of fire : Carbon dioxide. Carbon monoxide. Toxic fumes may be released.

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Exercise caution when fighting any chemical fire. Prevent fire fighting water from entering the environment.
Protection during firefighting : Do not enter fire area without proper personal protective equipment, including respiratory protection (EN137).

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Avoid contact with skin and eyes. Wear chemically protective gloves, lab coat or apron to prevent prolonged or repeated skin contact. Ensure adequate ventilation, especially in confined areas.

6.1.1. For non-emergency personnel

Protective equipment : See Section 8.
Emergency procedures : Evacuate unnecessary personnel.

6.1.2. For emergency responders

Protective equipment : Equip cleanup and emergency crew with proper protection.
Emergency procedures : Ventilate area.

6.2. Environmental precautions

Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

For containment : Collect all waste in suitable and labelled containers and dispose according to local legislation.
Methods for cleaning up : Collect spillage. Mechanically recover the product. Sweep or shovel into suitable containers. This material and its container must be disposed of in a safe way, and as per local legislation.

6.4. Reference to other sections

For further information refer to section 13. See Section 8. Exposure controls and personal protection.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Avoid contact with skin and eyes. Use personal protective equipment as required.
Hygiene measures : Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Do not eat, drink or smoke when using this product.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Keep container tightly closed. Store in original container. Store in a dry place. Store in a well-ventilated place.

7.3. Specific end use(s)

No additional data.

Safety Data Sheet

Repin (brown)

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

8.1.1 National occupational exposure and biological limit values

rosin; colophony (8050-09-7)	
United Kingdom - Occupational Exposure Limits	
Local name	Rosin-based solder flux fume
WEL TWA (OEL TWA) [1]	0.05 mg/m ³
WEL STEL (OEL STEL)	0.15 mg/m ³
Remark	Sen (Capable of causing occupational asthma)
Regulatory reference	EH40/2005 (Fourth edition, 2020). HSE

8.1.2. Recommended monitoring procedures

No additional information available

8.1.3. Air contaminants formed

No additional information available

8.1.4. DNEL and PNEC

No additional information available

8.1.5. Control banding

No additional information available

8.2. Exposure controls

8.2.1. Appropriate engineering controls

Appropriate engineering controls:

Ensure good ventilation of the work station.

8.2.2. Personal protection equipment

Personal protective equipment:

Gloves. Safety glasses.



8.2.2.1. Eye and face protection

Eye protection:

Use splash goggles when eye contact due to splashing is possible. STANDARD EN 166.

8.2.2.2. Skin protection

Skin and body protection:

Wear suitable protective clothing

Hand protection:

Wear suitable gloves. Nitrile rubber gloves. Layer thickness : 0,11mm. Breakthrough time : >480 min. STANDARD EN 374.

8.2.2.3. Respiratory protection

Respiratory protection:

Not required

8.2.2.4. Thermal hazards

No additional information available

8.2.3. Environmental exposure controls

Other information:

Do not eat, drink or smoke during use. Wash hands and other exposed areas with mild soap and water before eating, drinking or smoking and when leaving work. Personal protective equipment should be chosen according to the CEN standards and in discussion with the supplier of the protective equipment.

Safety Data Sheet

Repin (brown)

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Liquid
Colour	: Yellowish.
Appearance	: Viscous.
Odour	: scent of cloves.
Odour threshold	: Not determined.
Melting point	: Not determined.
Freezing point	: Not determined.
Boiling point	: Not determined.
Flammability	: Not available
Explosive properties	: Product is not explosive.
Oxidising properties	: Non flammable.
Explosive limits	: Not determined.
Lower explosive limit (LEL)	: Not available
Upper explosive limit (UEL)	: Not available
Flash point	: Not determined.
Auto-ignition temperature	: Not determined.
Decomposition temperature	: Not determined.
pH	: Not determined.
Viscosity, kinematic	: Not determined.
Viscosity, dynamic	: Not determined.
Solubility	: Not determined.
Partition coefficient n-octanol/water (Log Kow)	: Not available
Partition coefficient n-octanol/water (Log Pow)	: Not determined.
Vapour pressure	: Not determined.
Vapour pressure at 50 °C	: Not available
Density	: Not available
Relative density	: Not determined.
Relative vapour density at 20 °C	: Not determined.
Particle size	: Not applicable
Particle size distribution	: Not applicable
Particle shape	: Not applicable
Particle aspect ratio	: Not applicable
Particle aggregation state	: Not applicable
Particle agglomeration state	: Not applicable
Particle specific surface area	: Not applicable
Particle dustiness	: Not applicable

9.2. Other information

9.2.1. Information with regard to physical hazard classes

No additional information available

9.2.2. Other safety characteristics

Relative evaporation rate (butylacetate=1) : Not determined.

SECTION 10: Stability and reactivity

10.1. Reactivity

The product is non-reactive under normal conditions of use, storage and transport.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No polymerization.

10.4. Conditions to avoid

None to our knowledge.

10.5. Incompatible materials

None to our knowledge.

10.6. Hazardous decomposition products

No decomposition if stored and used normally.

Safety Data Sheet

Repin (brown)

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute toxicity (oral)	: Not classified
Acute toxicity (dermal)	: Not classified
Acute toxicity (inhalation)	: Not classified
Additional information	: Based on available data, the classification criteria are not met
Skin corrosion/irritation	: Not classified pH: Not determined.
Additional information	: Based on available data, the classification criteria are not met
Serious eye damage/irritation	: Causes serious eye irritation. pH: Not determined.
Respiratory or skin sensitisation	: May cause an allergic skin reaction.
Germ cell mutagenicity	: Not classified
Additional information	: Based on available data, the classification criteria are not met
Carcinogenicity	: Not classified
Additional information	: Based on available data, the classification criteria are not met
Reproductive toxicity	: Not classified
Additional information	: Based on available data, the classification criteria are not met
STOT-single exposure	: Not classified
Additional information	: Based on available data, the classification criteria are not met
STOT-repeated exposure	: Not classified
Additional information	: Based on available data, the classification criteria are not met
Aspiration hazard	: Not classified
Additional information	: Based on available data, the classification criteria are not met

Repin (brown)	
Viscosity, kinematic	Not determined.

11.2. Information on other hazards

11.2.1. Endocrine disrupting properties

Adverse health effects caused by endocrine disrupting properties	: The mixture does not contain substance(s) included in the list established in accordance with Article 59(1) of REACH for having endocrine disrupting properties, or is not identified as having endocrine disrupting properties in accordance with the criteria set out in Commission Delegated Regulation (EU) 2017/2100 or Commission Regulation (EU) 2018/605 at a concentration equal to or greater than 0,1 %
--	--

11.2.2 Other information

Potential adverse human health effects and symptoms	: For further information see section 4
---	---

SECTION 12: Ecological information

12.1. Toxicity

Hazardous to the aquatic environment, short-term (acute)	: Not classified
Hazardous to the aquatic environment, long-term (chronic)	: Not classified

rosin; colophony (8050-09-7)	
LC50 - Fish [1]	10000 mg/l (96 hours - Brachydanio rerio, zebra-fish)
EC50 - Crustacea [1]	4.5 mg/l (48 hours - Daphnia magna)
ErC50 algae	400 mg/l (IC50, 72 hours -Selenastrum)

12.2. Persistence and degradability

Repin (brown)	
Persistence and degradability	No data available.

Safety Data Sheet

Repin (brown)

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

rosin; colophony (8050-09-7)	
Persistence and degradability	Not readily biodegradable.
Biodegradation	70 – 80 % (OECD 301C method)

12.3. Bioaccumulative potential

Repin (brown)	
Partition coefficient n-octanol/water (Log Pow)	Not determined.
Bioaccumulative potential	No data.

12.4. Mobility in soil

Repin (brown)	
Ecology - soil	No data available.

12.5. Results of PBT and vPvB assessment

Repin (brown)	
This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII	
This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII	

12.6. Endocrine disrupting properties

No additional information available

12.7. Other adverse effects

Other adverse effects : None to our knowledge.
Additional information : No other effects known

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods : Must not be disposed together with household garbage. Do not allow product to reach the sewage system.
Product/Packaging disposal recommendations : Dispose in a safe manner in accordance with local/national regulations. Empty container completely. After cleaning, recycle or dispose of at an authorised site.
European List of Waste (LoW) code : 18 01 06* - chemicals consisting of or containing dangerous substances

SECTION 14: Transport information

In accordance with ADR / IMDG / IATA / ADN / RID /

ADR	IMDG	IATA	ADN	RID
14.1. UN number or ID number				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
14.2. UN proper shipping name				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
14.3. Transport hazard class(es)				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
14.4. Packing group				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
14.5. Environmental hazards				
Not applicable	Not applicable	Not applicable	Not applicable	Not applicable
No supplementary information available				

Safety Data Sheet

Repin (brown)

according to Regulation (EC) No. 1907/2006 (REACH) with its amendment Regulation (EU) 2020/878

14.6. Special precautions for user

Overland transport

Not applicable

Transport by sea

Not applicable

Air transport

Not applicable

Inland waterway transport

Not applicable

Rail transport

Not applicable

14.7. Maritime transport in bulk according to IMO instruments

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains no substance on the REACH candidate list

Contains no REACH Annex XIV substances

Contains no substance subject to Regulation (EU) No 649/2012 of the European Parliament and of the Council of 4 July 2012 concerning the export and import of hazardous chemicals.

Contains no substance subject to Regulation (EU) No 2019/1021 of the European Parliament and of the Council of 20 June 2019 on persistent organic pollutants

15.1.2. National regulations

REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006

15.2. Chemical safety assessment

No chemical safety assessment has been carried out for the substance or the mixture by the supplier

SECTION 16: Other information

Indication of changes:

Section	Changed item	Change	Comments
2.2	Precautionary statements (CLP)	Modified	

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Full text of H- and EUH-statements:

Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
H317	May cause an allergic skin reaction.
H319	Causes serious eye irritation.
Skin Sens. 1	Skin sensitisation, Category 1

The information in this safety data sheet is based on information from the manufacturer/supplier, present european and national legislation, and presupposes that the product is used within the specified area of application.