



SAFETY DATA SHEET

page 1/5
SHOFU DISK ZR Lucent
SHOFU DISK ZR Lucent Supra
Printing date: October 25, 2021

SECTION 1. Identification of the substance/mixture and of the company/undertaking

- 1.1 Product identifier
Trade Name:
SHOFU DISK ZR Lucent
SHOFU DISK ZR Lucent Supra
- 1.2 Relevant identified uses of the substance or mixture and uses advised against
Relevant identified uses: Dental material
Uses advised against: No further data
- 1.3 Details of the supplier of the safety data sheet
Company/Undertaking identification
Manufacturer's Name: SHOFU DENTAL GmbH
Address: An der Pönt 70, 40885 Ratingen, Germany
Phone: +49 (0) 2102-8664-0
Fax: +49 (0) 2102-8664-64
E-Mail: info@shofu.de
Section in Charge: Quality Management & Regulatory Affairs
- 1.4 Emergency Telephone Number
+49-2102-8664-53 (SHOFU DENTAL GmbH) 24 hours / 7 days

SECTION 2. Hazards identification

- 2.1 Classification of the substance or mixture
CLASSIFICATION (EC 1272/2008)
Not applicable.
- 2.2 Label elements
LABEL IN ACCORDANCE WITH (EC) NO.1272/2008
Void
HAZARD-DETERMINING COMPONENTS OF LABELLING
Void
SIGNAL WORD
Void
HAZARD STATEMENTS
Void
PRECAUTIONARY STATEMENTS
Void
- 2.3 Other hazards
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.

SECTION 3. Composition/information on ingredients

- 3.1 Chemical characterization: Mixtures
3.2 Description: Mixture of substances listed below with nonhazardous additions.

(Contd. on page2)



SAFETY DATA SHEET

page 2/5
SHOFU DISK ZR Lucent
SHOFU DISK ZR Lucent Supra
Printing date: October 25, 2021

(Contd. of page1)

- 3.3 Dangerous components: Void
Ceramic materials
- 3.4 Additional information: For the wording of the listed risk phrases refer to section 2.

SECTION 4. First-aid measures

- 4.1 Description of first aid measures
(Grinding dust)
Eye contact: Rinse cautiously with water for several minutes. If eye irritation persists, get medical advice/attention.
Skin contact: Rinse with water.
Ingestion: Rinse mouth and seek medical advice if necessary.
Inhalation: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
If experiencing respiratory symptoms: Call a POISON CENTER/doctor.
- 4.2 Most important symptoms and effects, both acute and delayed
No further relevant information available.
- 4.3 Indication of any immediate medical attention and special treatment needed
No further relevant information available.

SECTION 5. Fire-fighting measures

- 5.1 Extinguishing Media:
The product is not flammable.
In case of fire, use fire extinguishing methods suitable to surrounding conditions.
- 5.2 Special hazards arising from the substance or mixture:
No further relevant information available.
- 5.3 Advice for firefighters:
No further relevant information available.

SECTION 6. Accidental release measures

- 6.1 Personal precautions, protective equipment and emergency procedures:
Do not breathe dust.
- 6.2 Environmental Precautions:
No special measures required.
- 6.3 Methods and material for containment and cleaning Up:
Pick up mechanically.
Dispose of absorbed material in accordance with local/national/international regulation.
- 6.4 Reference to other section:
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

SECTION 7. Handling and storage

- 7.1 Precautions for safe handling:
Use in accordance with instruction. In case of grinding and polishing, wear dust mask and safety glasses, and using dust respirator and local exhaust.

(Contd. on page3)



SAFETY DATA SHEET

page 3/5
SHOFU DISK ZR Lucent
SHOFU DISK ZR Lucent Supra
Printing date: October 25, 2021

(Contd. of page2)

7.2 Conditions for safe storage, including any incompatibilities:

Store in a cool, dry and dark place.

7.3 Specific end use(s):

No further relevant information available.

SECTION 8. Exposure controls/personal protection

8.1 Control parameters:

Not further data; see item 7.

8.2 Exposure controls:

Respiratory Protection:	Dust respirator and/or local exhaust.
Hand Protection:	Not required.
Eye Protection:	Safety goggles (May be required.)
Skin and body Protection:	Not required.

SECTION 9. Physical and chemical properties

9.1 Information on basic physical and chemical properties:

Appearance/Odor/Colour:	Odorless solid.
Odour threshold	Not determined.
pH	Not determined.
Melting point/freezing point	Not determined.
Boiling Point:	Not determined.
Flash point:	Not determined.
Evaporation rate	Not determined.
Flammability (solid, gas)	Not applicable.
Upper/lower flammability or explosive limits	Not determined.
Vapour pressure	Not determined.
Vapour density	Not determined.
Relative Density:	Not determined.
Solubility: water solubility	Insoluble
Partition coefficient: n-octanol/water	Not determined.
Auto-ignition temperature	Not determined.
Decomposition temperature	Not determined.
Viscosity	Not determined.
Explosive properties	Not applicable.
Oxidising properties	Not applicable.

9.2 Other information:

No further relevant information available.

SECTION 10. Stability and reactivity

10.1 Reactivity:

No further relevant information available.

10.2 Chemical stability:

Stable under normal temperatures and pressures.

10.3 Possibility of hazardous reactions:

No dangerous reactions known.

(Contd. on page4)



SAFETY DATA SHEET

page 4/5
SHOFU DISK ZR Lucent
SHOFU DISK ZR Lucent Supra
Printing date: October 25, 2021

(Contd. of page3)

10.4 Condition to Avoid:

No further relevant information available.

10.5 Incompatible materials:

No further relevant information available.

10.6 Hazardous Decomposition Products:

None under normal conditions of storage and use.

SECTION 11. Toxicological information

11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008:

Acute toxicity: Based on available data, the classification criteria are not met.

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

Eye damage/irritation: Based on available data, the classification criteria are not met.

Sensitization to the respiratory tract:

Based on available data, the classification criteria are not met.

Skin sensitization: Based on available data, the classification criteria are not met.

Germ cell mutagenicity/Genotoxicity:

Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met.

Reproductive toxicity: Based on available data, the classification criteria are not met.

Effects on or via lactation: Lack of data.

Specific target organ toxicity (single exposure):

Based on available data, the classification criteria are not met.

Specific target organ toxicity (repeated exposure):

Based on available data, the classification criteria are not met.

Aspiration hazard: Based on available data, the classification criteria are not met.

11.2 Information on other hazards:

No further relevant information available.

SECTION 12. Ecological information

12.1 Toxicity:

No further relevant information available.

12.2 Persistence and degradability:

No further relevant information available.

12.3 Bioaccumulative potential:

No further relevant information available.

12.4 Mobility in soil:

No further relevant information available.

12.5 Results of PBT and vPvB assessment:

Not applicable.

12.6 Endocrine disrupting properties:

This product does not contain any known or suspected endocrine disruptors.

12.7 Other adverse effects:

No further relevant information available.

(Contd. on page5)

Version Number 2

Revision date: October 25, 2021



SAFETY DATA SHEET

page 5/5
SHOFU DISK ZR Lucent
SHOFU DISK ZR Lucent Supra
Printing date: October 25, 2021

(Contd. of page4)

SECTION 13. Disposal considerations

13.1 Waste treatment methods:

Dispose of contents/container to in accordance with local/regional/national/international regulations.

SECTION 14. Transport information

- | | |
|---|--|
| 14.1 UN number or ID number: | Void |
| 14.2 UN proper shipping name: | Void |
| 14.3 Transport hazard class(es): | Void |
| 14.4 Packing group: | Void |
| 14.5 Environmental hazards: | No further relevant information available. |
| 14.6 Special precautions for user: | Not applicable. |
| 14.7 Maritime transport in bulk according to IMO instruments: | Not applicable. |

SECTION 15. Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture:

- EU REGULATIONS: See Section 2
- Other regulations, limitations and prohibitive regulations:

The product is a medical device according to the EC-directive 93/42/EEC.

15.2 Chemical safety assessment:

A Chemical Safety Assessment has not been carried out.

SECTION 16. Other information

This product is intended for use by dental professionals. (instrument/material)

Abbreviations and acronyms:

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative