

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 03/01/2006

SDS No : 32-001.7

Date Revised : 09/23/2016

Revision No : 6

Ultradent® EDTA 18% Solution**SECTION 1 : Identification of the substance/preparation and of the company/undertaking****1.1. Product identifier**

Product code : ET/66009
 Product name : Ultradent® EDTA 18% Solution
 Product description : Root canal chelating conditioner and cleanser

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses : Professional Dental Root Canal Cleansing & Conditioning Agent

1.3. Details of the supplier of the safety data sheet**Manufacturer**

Ultradent Products, Inc.
 505 W. 10200 S.
 South Jordan, UT 84095

EC Responsible Person

Ultradent Products GmbH
 Am Westhover Berg 30
 51149 Cologne Germany
 Email: infoDE@ultradent.com
Emergency Phone : +49(0)2203-35-92-0

1.4. Emergency telephone number

CHEMTREC (NORTH AMERICA) : (800) 424 - 9300
 (INTERNATIONAL) : +1(703) 527 - 3887

SECTION 2: Hazards identification**2.1. Classification of the substance or mixture****Classification according to Regulation(EC) No 1272/2008 [CLP]**

: The mixture has been assessed and/or tested for its physical, health and environmental hazards and the following classification applies.

Health : Acute Toxicity (Oral), Category 3
 Eye Irritation, Category 2A

2.2. Label elements**Classification according to Regulation(EC) No 1272/2008 [CLP]**

Hazard pictogram(s) :



Exclamation
mark

Signal Word : WARNING

Hazard statement(s) : H319: Causes serious eye irritation.
 H302: Harmful if swallowed.

Precautionary statement(s)

Prevention : P280: Wear protective gloves/protective clothing/eye protection/face protection.

Response : P305: IF IN EYES:

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010
Dental Use



Date Issued : 03/01/2006
SDS No : 32-001.7
Date Revised : 09/23/2016
Revision No : 6

Ultradent® EDTA 18% Solution

P351: Rinse cautiously with water for several minutes.
P337+P313: If eye irritation persists: Get medical advice/attention.
P302+P352: IF ON SKIN: Wash with plenty of soap and water.
P332+P313: If skin irritation occurs: Get medical advice/attention.
P301: IF SWALLOWED:
P331: Do NOT induce vomiting.
P301+P312: IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.

Storage : P403+P235: Store in a well-ventilated place. Keep cool.
P405: Store locked up.

Disposal : P501: Dispose of in compliance with governmental regulation. (EC1975L0442-20/11/2003)

2.3. Other hazards

Conclusion The substance does not fulfill the PBT and vPvB criteria for screening assessment; there are no indications of P or B properties.

SECTION 3: Composition / information on ingredients

3.1. Substances

Not Applicable

3.2. Mixtures

Chemical Name	CAS	EINECS No.	Wt. %	Classification according to Regulation(EC) No 1272/2008 [CLP]
EDTA Tetrasodium Salt	64-02-8	200-573-9	< 20	Eye Irr., Cat. 2A; H319

For full text of H-statements: see SECTION 16.

SECTION 4: First aid measures

4.1. Description of first aid measures

Following eyes : Immediately flush eyes with plenty of water. Get medical attention, if irritation persists.

Following skin : Wash with soap and water. Get medical attention if irritation develops or persists.

Following ingestion : Give large volumes of water. Do not induce vomiting. Consult a physician if you feel unwell.

Following inhalation : No specific treatment is necessary since this material is not likely to be hazardous by inhalation.

4.2. Most important symptoms and effects, both acute and delayed

Eyes : Irritating to eyes.

Skin : May cause slight irritation.

Ingestion : Irritating to mouth, throat, and stomach.

Inhalation : None expected for this product.

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010
Dental Use



Date Issued : 03/01/2006
SDS No : 32-001.7
Date Revised : 09/23/2016
Revision No : 6

Ultradent® EDTA 18% Solution

4.3. Indication of any immediate medical attention and special treatment needed

Notes to physician : Not Applicable

SECTION 5: Fire fighting measures

5.1. Extinguishing media

Extinguishing media : Foam, dry chemical, carbon dioxide (CO₂).

5.2. Special hazards arising from the substance or mixture

Hazardous combustion products : None known.

5.3. Advice for firefighters

Fire fighting procedures : General: Evacuate all personnel; use protective equipment for fire-fighting. Use self-contained breathing apparatus when the product is involved in fire.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General procedures : Refer to Section 8 for Personal Protective Equipment.

6.2. Environmental precautions

Water spill : Do not allow to enter sewers or drains that may lead to waterways.

6.3. Methods and material for containment and cleaning up

Small Spill : Clean up spills immediately, observing precautions in Protective Equipment section.

Large Spill : Absorb with inert, damp non-combustible material, then flush area with water.

6.4. Reference to other sections

Reference to other sections : Not Applicable

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Handling : Use suitable protective equipment.

Storage : See product labeling.

7.2. Conditions for safe storage, including any incompatibilities

Shelf life : See product labeling

7.3. Specific end use(s)

Specific end use(s) : Professional Dental Root Canal Cleansing & Conditioning Agent

SECTION 8: Exposure controls / personal protection

8.1. Control parameters

Control parameters : None

8.2. Exposure controls

Eye/face protection : Wear eye protection

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010
Dental Use



Date Issued : 03/01/2006
SDS No : 32-001.7
Date Revised : 09/23/2016
Revision No : 6

Ultradent® EDTA 18% Solution

Skin protection : Wear suitable protective clothing and gloves.
Respiratory protection : Good general ventilation should be sufficient to control airborne levels.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Low viscosity liquid
Appearance : Clear
Colour : Colorless
Odour : Odorless or no characteristic odor
pH : 10.9 - 11.1
Specific gravity : 1.1 to 1.2

9.2. Other information

Percent volatile : Not Volatile

SECTION 10: Stability and reactivity

10.1. Reactivity

Reactivity : Stable

10.2. Chemical stability

Chemical stability : Stable when stored and handled under recommended conditions.

10.3. Possibility of hazardous reactions

Hazardous polymerization : None

10.4. Conditions to avoid

Conditions to avoid : Strong acids, strong bases, and oxidizing agents.

10.5. Incompatible materials

Incompatible materials : Not Determined

10.6. Hazardous decomposition products

Hazardous decomposition products : Carbon dioxide, carbon monoxide, nitrogen oxides, sodium hydroxide, and ammonia

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute

Dermal LD₅₀ : Not able to calculate ATE mix(dermal) due to lack of data on the components.

Oral LD₅₀ : > 150 mg/kg Mouse

Inhalation LC₅₀ : Not able to calculate ATE mix (inhalation) due to lack of data on the components.

Other information : Device is biocompatible when used as directed by dental professionals per ISO 10993-1.

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 03/01/2006

SDS No : 32-001.7

Date Revised : 09/23/2016

Revision No : 6

Ultradent® EDTA 18% Solution

SECTION 12: Ecological information

12.1. Toxicity

Aquatic toxicity (acute) : Do not allow to enter sewers or drains that may lead to waterways.

12.2. Persistence and degradability

Persistence and degradability : Not classified

12.3. Bioaccumulative potential

Bioaccumulative potential : Not classified

12.4. Mobility in soil

Mobility in soil : Not Available

12.5. Results of PBT and vPvB assessment

Results of PBT and vPvB assessment : Conclusion The substance does not fulfill the PBT and vPvB criteria for screening assessment; there are no indications of P or B properties.

12.6. Other adverse effects

Environmental data : Do not flush to sewer.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal method : Dispose of in compliance with governmental regulation.(EC 1975L0442-20/11/2003)

SECTION 14: Transport information

14.1. UN number

UN number : NAP

14.2. UN proper shipping name

UN proper shipping name : Not Dangerous Goods

14.3. Transport hazard class(es)

Hazard classification : NAP

14.4. Packing group

Packing group : NAP

14.5. Environmental hazards

Marine pollutant #1 : NAP

14.6. Special precautions for user

ADR - road : NAP

RID - rail : NAP

IMDG - sea : NAP

IATA - air : NAP

14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

Transport in bulk : NAP

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010
Dental Use



Date Issued : 03/01/2006
SDS No : 32-001.7
Date Revised : 09/23/2016
Revision No : 6

Ultradent® EDTA 18% Solution

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

General comments : Please refer to Medical Devices Directive 93/42/EEC.

15.2. Chemical safety assessment

Chemical safety assessment : See Section 11

Additional information : None of the ingredients are listed as a carcinogen by IARC, NTP, or CA Prop 65.
This product is not considered dangerous for transportation.

SECTION 16: Other information

Relevant H-statements (number and full text) : Eye Irr., Cat. 2A: Eye Irritation, Category 2A
H319: Causes serious eye irritation.

Prepared by : Anu Kattoju

Revision summary : This SDS replaces the 04/10/2015 SDS. Revised: **Section 1:** MSDS No. **Section 2:** EMERGENCY OVERVIEW - IMMEDIATE CONCERNS 2.3. Other hazards. **Section 4:** INGESTION, INHALATION. **Section 5:** EXPLOSION HAZARDS, EXTINGUISHING MEDIA, FIRE EXPLOSION, FIRE FIGHTING EQUIPMENT, SENSITIVE TO STATIC DISCHARGE, SENSITIVITY TO IMPACT. **Section 6:** WATER SPILL. **Section 8:** . **Section 9:** pH, PHYSICAL STATE, SPECIFIC GRAVITY. **Section 11:** ACUTE (INHALATION LC₅₀ (rat), ORAL LD₅₀ (rat), ORAL LD₅₀ (rat), ORAL LD₅₀ (rat) ACUTE, GENERAL COMMENTS. **Section 12:** CHEMICAL FATE INFORMATION, DISTRIBUTION, ENVIRONMENTAL DATA, BIOACCUMULATION/ACCUMULATION, AQUATIC TOXICITY (ACUTE), AQUATIC TOXICITY (ACUTE) (48-HOUR EC₅₀, 96-HOUR EC₅₀), SECTION 12: Ecological information. **Section 14:** ROAD AND RAIL (ADR/RID), VESSEL (IMO/IMDG) - MARINE POLLUTANT #1 ROAD AND RAIL (UK only) (CDG) - SPECIAL PROVISIONS SECTION 14: Transport information, UN number, Packing group, ADR - road, IMDG - sea, IATA - air. **Section 15:** COMMENTS, GENERAL COMMENTS, RoHS. **Section 16:** GENERAL STATEMENTS.

General statements : NAP = Not Applicable
NE = Not Established
TLV = Threshold Limit Value
PEL = Permissible exposure limits
MAK = Maximum Workplace Concentration
STEL = Short-term exposure limit
STEV = Short-term exposure value
TWA = Time Weighted Average
PPE = Personal Protective Equipment

Manufacturer disclaimer : FOR DENTAL USE ONLY: Use as directed. The information and recommendations are taken from sources (raw material SDS(s) and manufacturer's knowledge) believed to be accurate; however, the manufacturer, makes no warranty with respect to the accuracy of the information or the suitability of the recommendation and assumes no liability to any user thereof.

SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 03/01/2006

SDS No : 32-001.7

Date Revised : 09/23/2016

Revision No : 6

Ultradent® EDTA 18% Solution

Each user should review these recommendations in the specific context of the intended use and determine whether they are appropriate.