



**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 28.09.2017

Version number 3

Revision: 28.09.2017

**SECTION 1: Identification of the substance/mixture and of the company/undertaking**

· **1.1 Product identifier**

· Trade name: **Universal Adhesive**

· **1.2 Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.

· **Application of the substance / the mixture**  
Tray adhesive for use with silicone impression materials

· **1.3 Details of the supplier of the safety data sheet**

· **Manufacturer/Supplier:**

Kulzer GmbH

Leipziger Straße 2, 63450 Hanau (Germany)

Tel.: +49 (0)800 4372522

· **Informing department:** E-Mail: [msds@kulzer-dental.com](mailto:msds@kulzer-dental.com)

· **1.4 Emergency telephone number:** Emergency CONTACT (24-Hour-Number): +49 (0)6132-84463

**SECTION 2: Hazards identification**

· **2.1 Classification of the substance or mixture**

· **Classification according to Regulation (EC) No 1272/2008**

Flam. Liq. 2 H225 Highly flammable liquid and vapour.

Skin Irrit. 2 H315 Causes skin irritation.

Eye Dam. 1 H318 Causes serious eye damage.

· **2.2 Label elements**

· **Labelling according to Regulation (EC) No 1272/2008**

The product is classified and labelled according to the CLP regulation.

· **Hazard pictograms**



GHS02 GHS05

· **Signal word** Danger

· **Hazard-determining components of labelling:**

titanium tetrabutanolate

· **Hazard statements**

H225 Highly flammable liquid and vapour.

H315 Causes skin irritation.

H318 Causes serious eye damage.

· **Precautionary statements**

P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

P241 Use explosion-proof electrical/ventilating/lighting/equipment.

P303+P361+P353 IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER/doctor.

P321 Specific treatment (see on this label).

· **2.3 Other hazards -**

· **Results of PBT and vPvB assessment**

· **PBT:** Not applicable.

(Contd. on page 2)

GB

**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 28.09.2017

Version number 3

Revision: 28.09.2017

**Trade name: Universal Adhesive**· **vPvB:** Not applicable.

(Contd. of page 1)

**SECTION 3: Composition/information on ingredients**· **3.2 Chemical characterisation: Mixtures**· **Description:** -· **Dangerous components:**

|                                                                    |                                                                                                                |        |
|--------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------|--------|
| CAS: 141-78-6<br>EINECS: 205-500-4                                 | ethyl acetate<br>Flam. Liq. 2, H225; Eye Irrit. 2, H319; STOT SE 3, H336                                       | 10-25% |
| CAS: 1330-20-7<br>EINECS: 215-535-7<br>Reg.nr.: 2119488216-32-xxxx | xylene, mixed isomers, pure<br>Flam. Liq. 3, H226; Acute Tox. 4, H312; Acute Tox. 4, H332; Skin Irrit. 2, H315 | 5-10%  |
| CAS: 5593-70-4<br>EINECS: 227-006-8                                | titanium tetrabutanolate<br>Flam. Liq. 3, H226; Eye Dam. 1, H318; Skin Irrit. 2, H315                          | 0-5%   |
| CAS: 100-41-4<br>EINECS: 202-849-4                                 | ethylbenzene<br>Flam. Liq. 2, H225; STOT RE 2, H373; Asp. Tox. 1, H304; Acute Tox. 4, H332                     | 0-5%   |
| CAS: 556-67-2<br>EINECS: 209-136-7                                 | octamethylcyclotetrasiloxane<br>Flam. Liq. 3, H226; Repr. 2, H361f; Aquatic Chronic 4, H413                    | < 1%   |

· **Additional information** For the wording of the listed hazard phrases refer to section 16.**SECTION 4: First aid measures**· **4.1 Description of first aid measures**· **After inhalation** Supply fresh air; consult doctor in case of symptoms.· **After skin contact**

Instantly rinse with water.

If skin irritation continues, consult a doctor.

· **After eye contact**

Rinse opened eye for several minutes under running water. Then consult doctor.

· **After swallowing**

Rinse out mouth and then drink plenty of water.

In case of persistent symptoms consult doctor.

· **4.2 Most important symptoms and effects, both acute and delayed**

No further relevant information available.

· **4.3 Indication of any immediate medical attention and special treatment needed**

No further relevant information available.

**SECTION 5: Firefighting measures**· **5.1 Extinguishing media**· **Suitable extinguishing agents**CO<sub>2</sub>, extinguishing powder or water jet. Fight larger fire with alcohol-resistant foam.· **For safety reasons unsuitable extinguishing agents** Water with a full water jet.· **5.2 Special hazards arising from the substance or mixture**

Can form explosive gas-air mixtures.

Formation of toxic gases is possible during heating or in case of fire.

· **5.3 Advice for firefighters**· **Protective equipment:** No special measures required.

(Contd. on page 3)

GB



**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 28.09.2017

Version number 3

Revision: 28.09.2017

**Trade name: Universal Adhesive**

· **Additional information -**

(Contd. of page 2)

**SECTION 6: Accidental release measures**

- **6.1 Personal precautions, protective equipment and emergency procedures**  
Wear protective equipment. Keep unprotected persons away.
- **6.2 Environmental precautions:** Prevent material from reaching sewage system, holes and cellars.
- **6.3 Methods and material for containment and cleaning up:**  
Absorb with liquid-binding material (diatomite, universal binders, for small amounts tissues).  
Do not flush with water or aqueous cleansing agents
- **6.4 Reference to other sections**  
No dangerous materials are released.  
See Section 7 for information on safe handling  
See Section 8 for information on personal protection equipment.  
See Section 13 for information on disposal.  
-

**SECTION 7: Handling and storage**

- **7.1 Precautions for safe handling**  
Keep containers tightly sealed.  
Ensure good ventilation/exhaustion at the workplace.  
Prevent formation of aerosols.
  - **Information about protection against explosions and fires:**  
Keep ignition sources away - Do not smoke.  
Protect against electrostatic charges.
- **7.2 Conditions for safe storage, including any incompatibilities**
  - **Storage**
    - **Requirements to be met by storerooms and containers:** Store in cool location.
    - **Information about storage in one common storage facility:** Not required.
    - **Further information about storage conditions:**  
Store in cool, dry conditions in well sealed containers.
- **7.3 Specific end use(s)** No further relevant information available.

**SECTION 8: Exposure controls/personal protection**

- **Additional information about design of technical systems:** No further data; see item 7.
- **8.1 Control parameters**

| · <b>Components with critical values that require monitoring at the workplace:</b> |                                                                                                                  |
|------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|
| <b>141-78-6 ethyl acetate</b>                                                      |                                                                                                                  |
| WEL                                                                                | Short-term value: 400 ppm<br>Long-term value: 200 ppm                                                            |
| <b>1330-20-7 xylene, mixed isomers, pure</b>                                       |                                                                                                                  |
| WEL                                                                                | Short-term value: 441 mg/m <sup>3</sup> , 100 ppm<br>Long-term value: 220 mg/m <sup>3</sup> , 50 ppm<br>Sk; BMGV |
| <b>100-41-4 ethylbenzene</b>                                                       |                                                                                                                  |
| WEL                                                                                | Short-term value: 552 mg/m <sup>3</sup> , 125 ppm<br>Long-term value: 441 mg/m <sup>3</sup> , 100 ppm<br>Sk      |

(Contd. on page 4)



**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 28.09.2017

Version number 3

Revision: 28.09.2017

**Trade name: Universal Adhesive**

(Contd. of page 3)

|                                                    |                                                                                                          |
|----------------------------------------------------|----------------------------------------------------------------------------------------------------------|
| · <b>Ingredients with biological limit values:</b> |                                                                                                          |
| <b>1330-20-7 xylene, mixed isomers, pure</b>       |                                                                                                          |
| BMGV                                               | 650 mmol/mol creatinine<br>Medium: urine<br>Sampling time: post shift<br>Parameter: methyl hippuric acid |

· **Additional information:** The lists that were valid during the compilation were used as basis.

· **8.2 Exposure controls**

· **Personal protective equipment**

· **General protective and hygienic measures**

Avoid contact with the eyes.  
Keep away from foodstuffs, beverages and food.  
Instantly remove any soiled and impregnated garments.  
Wash hands during breaks and at the end of the work.  
Avoid contact with the eyes and skin.

· **Breathing equipment:**

Not necessary with efficient local exhaust. If exposition to vapours is possible, use breathing protective mask (filter A).

· **Protection of hands:**

Check protective gloves prior to each use for their proper condition.  
recommended

· **Material of gloves**

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· **Penetration time of glove material**

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· **For the permanent contact of a maximum of 15 minutes gloves made of the following materials are suitable:**

Butyl rubber, BR  
Nitrile rubber, NBR

· **Eye protection:** not absolutely necessary

· **Body protection:** Light weight protective clothing

**SECTION 9: Physical and chemical properties**

· **9.1 Information on basic physical and chemical properties**

· **General Information**

· **Appearance:**

· **Form:** Fluid  
· **Colour:** Blue  
· **Smell:** Ester-like  
· **Odour threshold:** Not determined.

· **pH-value:** Not determined.

· **Change in condition**

· **Melting point/Melting range:** Not determined  
· **Boiling point/Boiling range:** 77 °C

· **Flash point:** 3 °C

(Contd. on page 5)



**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 28.09.2017

Version number 3

Revision: 28.09.2017

**Trade name: Universal Adhesive**

(Contd. of page 4)

|                                                   |                                                                                            |
|---------------------------------------------------|--------------------------------------------------------------------------------------------|
| · <b>Inflammability (solid, gaseous)</b>          | Not applicable.                                                                            |
| · <b>Ignition temperature:</b>                    | 460 °C                                                                                     |
| · <b>Decomposition temperature:</b>               | Not determined.                                                                            |
| · <b>Self-inflammability:</b>                     | Product is not selfigniting.                                                               |
| · <b>Danger of explosion:</b>                     | Product is not explosive. However, formation of explosive air/vapour mixtures is possible. |
| · <b>Critical values for explosion:</b>           |                                                                                            |
| · <b>Lower:</b>                                   | 2.1 Vol %                                                                                  |
| · <b>Upper:</b>                                   | 11.5 Vol %                                                                                 |
| · <b>Steam pressure at 20 °C:</b>                 | 102 hPa                                                                                    |
| · <b>Density at 20 °C</b>                         | 1.20169 g/cm <sup>3</sup>                                                                  |
| · <b>Relative density</b>                         | Not determined.                                                                            |
| · <b>Vapour density</b>                           | Not determined.                                                                            |
| · <b>Evaporation rate</b>                         | Not determined.                                                                            |
| · <b>Solubility in / Miscibility with</b>         |                                                                                            |
| · <b>Water:</b>                                   | Not miscible or difficult to mix                                                           |
| · <b>Partition coefficient (n-octanol/water):</b> | Not determined.                                                                            |
| · <b>Viscosity:</b>                               |                                                                                            |
| · <b>dynamic:</b>                                 | Not determined.                                                                            |
| · <b>kinematic:</b>                               | Not determined.                                                                            |
| · <b>9.2 Other information</b>                    | No further relevant information available.                                                 |

### SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
  - **Conditions to be avoided:** No decomposition if used and stored according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known
- **10.4 Conditions to avoid** No further relevant information available.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** None
- **Additional information:** -

### SECTION 11: Toxicological information

- **11.1 Information on toxicological effects**
  - **Acute toxicity** Based on available data, the classification criteria are not met.

| · <b>LD/LC50 values that are relevant for classification:</b> |         |                                                      |
|---------------------------------------------------------------|---------|------------------------------------------------------|
| <b>141-78-6 ethyl acetate</b>                                 |         |                                                      |
| Oral                                                          | LD50    | 5600 mg/kg (rat)                                     |
| Dermal                                                        | LD50    | 18000 mg/kg (rabbit)                                 |
| Sensitisation                                                 | Sensib. | nicht sensib. (guinea pig) (OECD-Prüfrichtlinie 406) |
| <b>1330-20-7 xylene, mixed isomers, pure</b>                  |         |                                                      |
| Oral                                                          | LD50    | 8700 mg/kg (rat)                                     |
| Dermal                                                        | LD50    | 2000 mg/kg (rabbit)                                  |

(Contd. on page 6)

**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 28.09.2017

Version number 3

Revision: 28.09.2017

**Trade name: Universal Adhesive**

(Contd. of page 5)

|                              |          |                      |
|------------------------------|----------|----------------------|
| Inhalative                   | LC50/4 h | 6350 mg/l (rat)      |
| <b>100-41-4 ethylbenzene</b> |          |                      |
| Oral                         | LD50     | 3500 mg/kg (rat)     |
| Dermal                       | LD50     | 17800 mg/kg (rabbit) |

- **Primary irritant effect:**
  - **Skin corrosion/irritation**  
Causes skin irritation.
  - **Serious eye damage/irritation**  
Causes serious eye damage.
- **Respiratory or skin sensitisation**  
Based on available data, the classification criteria are not met.
- **CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)**
  - **Germ cell mutagenicity** Based on available data, the classification criteria are not met.
  - **Carcinogenicity** Based on available data, the classification criteria are not met.
  - **Reproductive toxicity** Based on available data, the classification criteria are not met.
- **STOT-single exposure** Based on available data, the classification criteria are not met.
- **STOT-repeated exposure** Based on available data, the classification criteria are not met.
- **Aspiration hazard** Based on available data, the classification criteria are not met.

**SECTION 12: Ecological information**· **12.1 Toxicity**· **Aquatic toxicity:****141-78-6 ethyl acetate**

|          |                                |
|----------|--------------------------------|
| EC50/48h | 3300 mg/l (algae)              |
|          | 717 mg/l (daphnia) (DIN 38412) |
| LC50/96h | 230 mg/l (fish)                |

- **12.2 Persistence and degradability** No further relevant information available.
- **12.3 Bioaccumulative potential** No further relevant information available.
- **12.4 Mobility in soil** No further relevant information available.
- **12.5 Results of PBT and vPvB assessment**
  - **PBT:** Not applicable.
  - **vPvB:** Not applicable.
- **12.6 Other adverse effects** No further relevant information available.

**SECTION 13: Disposal considerations**· **13.1 Waste treatment methods**· **Recommendation**

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

· **Uncleaned packagings:**

· **Recommendation:** Disposal must be made according to official regulations.

**SECTION 14: Transport information**· **14.1 UN-Number**· **ADR, IMDG, IATA**

UN1866

(Contd. on page 7)

GB



**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 28.09.2017

Version number 3

Revision: 28.09.2017

**Trade name: Universal Adhesive**

(Contd. of page 6)

· **14.2 UN proper shipping name**

- ADR
- IMDG, IATA

1866 RESIN SOLUTION, special provision 640D  
RESIN SOLUTION

· **14.3 Transport hazard class(es)**

- ADR



- Class
- Label

3 (F1) Flammable liquids.  
3

- IMDG, IATA



- Class
- Label

3 Flammable liquids.  
3

· **14.4 Packing group**

- ADR, IMDG, IATA

II

· **14.5 Environmental hazards:**

- Marine pollutant:

No

· **14.6 Special precautions for user**

- Kemler Number:
- EMS Number:
- Stowage Category

Warning: Flammable liquids.  
33  
F-E, S-E  
B

· **14.7 Transport in bulk according to Annex II  
of Marpol and the IBC Code**

Not applicable.

- Transport/Additional information:

-

- ADR

- Limited quantities (LQ)
- Excepted quantities (EQ)

5L  
Code: E2  
Maximum net quantity per inner  
packaging: 30 ml  
Maximum net quantity per outer  
packaging: 500 ml

- Transport category
- Tunnel restriction code

2  
D/E

- IMDG

- Limited quantities (LQ)
- Excepted quantities (EQ)

5L  
Code: E2  
Maximum net quantity per inner  
packaging: 30 ml  
Maximum net quantity per outer  
packaging: 500 ml

(Contd. on page 8)

GB



**Safety data sheet  
according to 1907/2006/EC, Article 31**

Printing date 28.09.2017

Version number 3

Revision: 28.09.2017

**Trade name: Universal Adhesive**

(Contd. of page 7)

· **UN "Model Regulation":**

UN 1866 RESIN SOLUTION, SPECIAL  
PROVISION 640D, 3, II

### SECTION 15: Regulatory information

· **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**

· **Directive 2012/18/EU**

- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **Seveso category P5c** FLAMMABLE LIQUIDS
- **Qualifying quantity (tonnes) for the application of lower-tier requirements** 5,000 t
- **Qualifying quantity (tonnes) for the application of upper-tier requirements** 50,000 t
- **REGULATION (EC) No 1907/2006 ANNEX XVII** Conditions of restriction: 3

· **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### SECTION 16: Other information

These data are based on our present knowledge. However, they shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Relevant phrases**

- H225 Highly flammable liquid and vapour.
- H226 Flammable liquid and vapour.
- H304 May be fatal if swallowed and enters airways.
- H312 Harmful in contact with skin.
- H315 Causes skin irritation.
- H318 Causes serious eye damage.
- H319 Causes serious eye irritation.
- H332 Harmful if inhaled.
- H336 May cause drowsiness or dizziness.
- H361f Suspected of damaging fertility.
- H373 May cause damage to organs through prolonged or repeated exposure.
- H413 May cause long lasting harmful effects to aquatic life.

· **Abbreviations and acronyms:**

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

vPvB: very Persistent and very Bioaccumulative

Flam. Liq. 2: Flammable liquids – Category 2

Flam. Liq. 3: Flammable liquids – Category 3

Acute Tox. 4: Acute toxicity – Category 4

Skin Irrit. 2: Skin corrosion/irritation – Category 2

Eye Dam. 1: Serious eye damage/eye irritation – Category 1

Eye Irrit. 2: Serious eye damage/eye irritation – Category 2

Repr. 2: Reproductive toxicity – Category 2

STOT SE 3: Specific target organ toxicity (single exposure) – Category 3

STOT RE 2: Specific target organ toxicity (repeated exposure) – Category 2

Asp. Tox. 1: Aspiration hazard – Category 1

Aquatic Chronic 4: Hazardous to the aquatic environment - long-term aquatic hazard – Category 4