



SAFETY DATA SHEET

page 1/5
VINTAGE Halo OPAQUE
PORCELAIN POWDER
Printing date: June 23. 2016

SECTION 1. Identification of the substance/mixture and of the company/undertaking

- 1.1 Product identifier
Trade Name:
VINTAGE Halo OPAQUE PORCELAIN "POWDER"
VINTAGE Halo RED SHIFT OPAQUE POCELAIN "POWDER"
VINTAGE Halo OPAQUE MODIFIER
- 1.2 Relevant identified uses of the substance or mixture and uses advised against
Relevant identified uses: Dental material
Uses advised against: No further data
- 1.3 Details of the supplier of the safety data sheet
Company/Undertaking identification
Manufacturer's Name: SHOFU DENTAL GmbH
Address: Am Brüll 17, 40878 Ratingen, Germany
Phone: +49 (0) 2102-8664-0
Fax: +49 (0) 2102-8664-64
E-Mail: info@shofu.de
Section in Charge: Quality Management & Regulatory Affairs
- 1.4 Emergency Telephone Number
+49-2102-8664-53 (SHOFU DENTAL GmbH) 24 hours / 7 days

SECTION 2. Hazards identification

- 2.1 Classification of the substance or mixture
CLASSIFICATION (EC 1272/2008)
Not applicable.
- 2.2 Label elements
LABEL IN ACCORDANCE WITH (EC) NO.1272/2008
Void
HAZARD-DETERMINING COMPONENTS OF LABELLING
SIGNAL WORD
Void
HAZARD STATEMENTS
Void
PRECAUTIONARY STATEMENTS
Void
- 2.3 Other hazards
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.

SECTION 3. Composition/information on ingredients

- 3.1 Chemical characterization: Mixtures
- 3.2 Description: Mixture of substances listed below with nonhazardous additions.
- 3.3 Dangerous components: Void
Ceramic materials
- 3.4 Additional information: For the wording of the listed risk phrases refer to section 2.

(Contd. on page2)



SAFETY DATA SHEET

page 2/5
VINTAGE Halo OPAQUE
PORCELAIN POWDER
Printing date: June 23, 2016

(Contd. of page1)

SECTION 4. First-aid measures

- 4.1 Description of first aid measures
- Eye contact: Rinse cautiously with water for several minutes. If eye irritation persists, get medical advice/attention.
- Skin contact: Rinse with water.
- Ingestion: Rinse mouth and seek medical advice if necessary.
- Inhalation: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
- 4.2 Most important symptoms and effects, both acute and delayed
- No further relevant information available.
- 4.3 Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

SECTION 5. Fire-fighting measures

- 5.1 Extinguishing Media:
- The product is not flammable.
- 5.2 Special hazards arising from the substance or mixture:
- No further relevant information available.
- 5.3 Advice for firefighters:
- No further relevant information available.

SECTION 6. Accidental release measures

- 6.1 Personal precautions, protective equipment and emergency procedures:
- Don't breathe dust.
- 6.2 Environmental Precautions:
- No special measures required.
- 6.3 Methods and material for containment and cleaning up:
- Sweep and wipe up.
- 6.4 Reference to other section:
- See Section 7 for information on safe handling.
- See Section 8 for information on personal protection equipment.
- See Section 13 for disposal information.

SECTION 7. Handling and storage

- 7.1 Precautions for safe handling:
- Use in accordance with instruction. In case of grinding and polishing, wear dust mask and safety glasses, and using dust respirator and local exhaust.
- 7.2 Conditions for safe storage, including any incompatibilities:
- Store in a cool and dry conditions with lid tightly closed.
- 7.3 Specific end use(s):
- No further relevant information available.

SECTION 8. Exposure controls/personal protection

- 8.1 Control parameters:
- Ingredients with limit values that require monitoring at the workplace:

(Contd. on page3)



SAFETY DATA SHEET

page 3/5
VINTAGE Halo OPAQUE
PORCELAIN POWDER
Printing date: June 23. 2016

(Contd. of page2)

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

8.2 Exposure controls:

Respiratory Protection:

Dust mask

Skin Protection: Hand Protection

Not required.

· Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye Protection: Safety goggles

SECTION 9. Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance/Odor/Colour:	Odorless powder
Odour threshold	Not determined.
pH	Not determined.
Melting point/freezing point	Not determined.
Boiling Point:	Not determined.
Flash point:	Not determined.
Evaporation rate	Not determined.
Flammability (solid, gas)	Not applicable.
Upper/lower flammability or explosive limits	Not determined.
Vapour pressure	Not determined.
Vapour density	Not determined.
Relative Density:	≈ 2.2 (water=1)
Solubility: water solubility	Insoluble
Partition coefficient: n-octanol/water	Not determined.
Auto-ignition temperature	Not determined.
Decomposition temperature	Not determined.
Viscosity	Not determined.
Explosive properties	Not applicable.
Oxidising properties	Not applicable.

9.2 Other information

No further relevant information available.

SECTION 10. Stability and reactivity

10.1 Reactivity:

No further relevant information available.

10.2 Chemical stability:

(Contd. on page4)



SAFETY DATA SHEET

page 4/5
VINTAGE Halo OPAQUE
PORCELAIN POWDER
Printing date: June 23, 2016

(Contd. of page3)

Stable under normal temperatures and pressures.

10.3 Possibility of hazardous reactions:

No dangerous reactions known.

10.4 Condition to Avoid:

No further relevant information available.

10.5 Incompatible materials:

No further relevant information available.

10.6 Hazardous Decomposition Products:

None under normal conditions of storage and use.

SECTION 11. Toxicological information

11.1 Information on toxicological effects:

Acute toxicity: Based on available data, the classification criteria are not met.

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

Eye damage/irritation: Based on available data, the classification criteria are not met.

Sensitization to the respiratory tract:

Based on available data, the classification criteria are not met.

Skin sensitization: Based on available data, the classification criteria are not met.

Germ cell mutagenicity/Genotoxicity:

Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met.

Reproductive toxicity: Based on available data, the classification criteria are not met.

Effects on or via lactation: Lack of data.

Specific target organ toxicity (single exposure):

Based on available data, the classification criteria are not met.

Specific target organ toxicity (repeated exposure):

Based on available data, the classification criteria are not met.

Aspiration hazard:

Based on available data, the classification criteria are not met.

SECTION 12. Ecological information

12.1 Toxicity:

No further relevant information available.

12.2 Persistence and degradability:

No further relevant information available.

12.3 Bioaccumulative potential:

No further relevant information available.

12.4 Mobility in soil:

No further relevant information available.

12.5 Results of PBT and vPvB assessment:

Not applicable.

12.6 Other adverse effects:

No further relevant information available.

(Contd. on page5)



SAFETY DATA SHEET

page 5/5
VINTAGE Halo OPAQUE
PORCELAIN POWDER
Printing date: June 23, 2016

(Contd. of page4)

SECTION 13. Disposal considerations

13.1 Waste treatment methods:

Dispose of contents/container to in accordance with local/regional/national/international regulations.

SECTION 14. Transport information

- | | |
|---|--|
| 14.1 UN number: | Void |
| 14.2 UN proper shipping name: | Void |
| 14.3 Transport hazard class(es): | Void |
| 14.4 Packing group: | Void |
| 14.5 Environmental hazards: | No further relevant information available. |
| 14.6 Special precautions for user: | Not applicable. |
| 14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code: | Not applicable. |

SECTION 15. Regulatory information

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture:
- EU REGULATIONS: See Section 2
 - Other regulations, limitations and prohibitive regulations:
The product is a medical device according to the EC-directive 93/42/EEC.
- 15.2 Chemical safety assessment:
A Chemical Safety Assessment has not been carried out.

SECTION 16. Other information

This product is intended for use by dental professionals. (instrument/material)

Abbreviations and acronyms:

PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative



SAFETY DATA SHEET

page 1/5
VINTAGE Halo
PORCELAIN POWDER
Printing date: June 23. 2016

SECTION 1. Identification of the substance/mixture and of the company/undertaking

- 1.1 Product identifier
Trade Name:
VINTAGE Halo PORCELAIN "POWDER"
**(Body/Opal/Color Powder(Effect)/Translucent/Add on/Cervical/Margin/
Opaque Dentin/Red Shift/Value plus)**
- 1.2 Relevant identified uses of the substance or mixture and uses advised against
Relevant identified uses: Dental material
Uses advised against: No further data
- 1.3 Details of the supplier of the safety data sheet
Company/Undertaking identification
Manufacturer's Name: SHOFU DENTAL GmbH
Address: Am Brüll 17, 40878 Ratingen, Germany
Phone: +49 (0) 2102-8664-0
Fax: +49 (0) 2102-8664-64
E-Mail: info@shofu.de
Section in Charge: Quality Management & Regulatory Affairs
- 1.4 Emergency Telephone Number
+49-2102-8664-53 (SHOFU DENTAL GmbH) 24 hours / 7 days

SECTION 2. Hazards identification

- 2.1 Classification of the substance or mixture
CLASSIFICATION (EC 1272/2008)
Not applicable.
- 2.2 Label elements
LABEL IN ACCORDANCE WITH (EC) NO.1272/2008
Void
HAZARD-DETERMINING COMPONENTS OF LABELLING
SIGNAL WORD
Void
HAZARD STATEMENTS
Void
PRECAUTIONARY STATEMENTS
Void
- 2.3 Other hazards
Results of PBT and vPvB assessment
PBT: Not applicable.
vPvB: Not applicable.

SECTION 3. Composition/information on ingredients

- 3.1 Chemical characterization: Mixtures
- 3.2 Description: Mixture of substances listed below with nonhazardous additions.
- 3.3 Dangerous components: Void
Ceramic materials
- 3.4 Additional information: For the wording of the listed risk phrases refer to section 2.

(Contd. on page2)



SAFETY DATA SHEET

page 2/5
VINTAGE Halo
PORCELAIN POWDER
Printing date: June 23, 2016

(Contd. of page1)

SECTION 4. First-aid measures

- 4.1 Description of first aid measures
 - Eye contact: Rinse cautiously with water for several minutes. If eye irritation persists, get medical advice/attention.
 - Skin contact: Rinse with water.
 - Ingestion: Rinse mouth and seek medical advice if necessary.
 - Inhalation: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
- 4.2 Most important symptoms and effects, both acute and delayed
 - No further relevant information available.
- 4.3 Indication of any immediate medical attention and special treatment needed
 - No further relevant information available.

SECTION 5. Fire-fighting measures

- 5.1 Extinguishing Media:
 - The product is not flammable.
- 5.2 Special hazards arising from the substance or mixture:
 - No further relevant information available.
- 5.3 Advice for firefighters:
 - No further relevant information available.

SECTION 6. Accidental release measures

- 6.1 Personal precautions, protective equipment and emergency procedures:
 - Don't breathe dust.
- 6.2 Environmental Precautions:
 - No special measures required.
- 6.3 Methods and material for containment and cleaning up:
 - Sweep and wipe up.
- 6.4 Reference to other section:
 - See Section 7 for information on safe handling.
 - See Section 8 for information on personal protection equipment.
 - See Section 13 for disposal information.

SECTION 7. Handling and storage

- 7.1 Precautions for safe handling:
 - Use in accordance with instruction. In case of grinding and polishing, wear dust mask and safety glasses, and using dust respirator and local exhaust.
- 7.2 Conditions for safe storage, including any incompatibilities:
 - Store in a cool and dry conditions with lid tightly closed.
- 7.3 Specific end use(s):
 - No further relevant information available.

SECTION 8. Exposure controls/personal protection

- 8.1 Control parameters:
 - Ingredients with limit values that require monitoring at the workplace:

(Contd. on page3)



SAFETY DATA SHEET

page 3/5
VINTAGE Halo
PORCELAIN POWDER
Printing date: June 23. 2016

(Contd. of page2)

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

8.2 Exposure controls:

Respiratory Protection:

Dust mask

Skin Protection: Hand Protection

Not required.

· Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

Eye Protection: Safety goggles

SECTION 9. Physical and chemical properties

9.1 Information on basic physical and chemical properties

Appearance/Odor/Colour:	Odorless powder
Odour threshold	Not determined.
pH	Not determined.
Melting point/freezing point	Not determined.
Boiling Point:	Not determined.
Flash point:	Not determined.
Evaporation rate	Not determined.
Flammability (solid, gas)	Not applicable.
Upper/lower flammability or explosive limits	Not determined.
Vapour pressure	Not determined.
Vapour density	Not determined.
Relative Density:	≈ 2.2 (water=1)
Solubility: water solubility	Insoluble
Partition coefficient: n-octanol/water	Not determined.
Auto-ignition temperature	Not determined.
Decomposition temperature	Not determined.
Viscosity	Not determined.
Explosive properties	Not applicable.
Oxidising properties	Not applicable.

9.2 Other information

No further relevant information available.

SECTION 10. Stability and reactivity

10.1 Reactivity:

No further relevant information available.

10.2 Chemical stability:

(Contd. on page4)



SAFETY DATA SHEET

page 4/5
VINTAGE Halo
PORCELAIN POWDER
Printing date: June 23, 2016

(Contd. of page3)

Stable under normal temperatures and pressures.

10.3 Possibility of hazardous reactions:

No dangerous reactions known.

10.4 Condition to Avoid:

No further relevant information available.

10.5 Incompatible materials:

No further relevant information available.

10.6 Hazardous Decomposition Products:

None under normal conditions of storage and use.

SECTION 11. Toxicological information

11.1 Information on toxicological effects:

Acute toxicity: Based on available data, the classification criteria are not met.

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

Eye damage/irritation: Based on available data, the classification criteria are not met.

Sensitization to the respiratory tract:

Based on available data, the classification criteria are not met.

Skin sensitization: Based on available data, the classification criteria are not met.

Germ cell mutagenicity/Genotoxicity:

Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met.

Reproductive toxicity: Based on available data, the classification criteria are not met.

Effects on or via lactation: Lack of data.

Specific target organ toxicity (single exposure):

Based on available data, the classification criteria are not met.

Specific target organ toxicity (repeated exposure):

Based on available data, the classification criteria are not met.

Aspiration hazard:

Based on available data, the classification criteria are not met.

SECTION 12. Ecological information

12.1 Toxicity:

No further relevant information available.

12.2 Persistence and degradability:

No further relevant information available.

12.3 Bioaccumulative potential:

No further relevant information available.

12.4 Mobility in soil:

No further relevant information available.

12.5 Results of PBT and vPvB assessment:

Not applicable.

12.6 Other adverse effects:

No further relevant information available.

(Contd. on page5)



SAFETY DATA SHEET

page 5/5
VINTAGE Halo
PORCELAIN POWDER
Printing date: June 23, 2016

(Contd. of page4)

SECTION 13. Disposal considerations

13.1 Waste treatment methods:

Dispose of contents/container to in accordance with local/regional/national/international regulations.

SECTION 14. Transport information

- | | |
|---|--|
| 14.1 UN number: | Void |
| 14.2 UN proper shipping name: | Void |
| 14.3 Transport hazard class(es): | Void |
| 14.4 Packing group: | Void |
| 14.5 Environmental hazards: | No further relevant information available. |
| 14.6 Special precautions for user: | Not applicable. |
| 14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code: | Not applicable. |

SECTION 15. Regulatory information

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture:
- EU REGULATIONS: See Section 2
 - Other regulations, limitations and prohibitive regulations:
The product is a medical device according to the EC-directive 93/42/EEC.
- 15.2 Chemical safety assessment:
A Chemical Safety Assessment has not been carried out.

SECTION 16. Other information

This product is intended for use by dental professionals. (instrument/material)

Abbreviations and acronyms:

PBT: Persistent, Bioaccumulative and Toxic
vPvB: very Persistent and very Bioaccumulative