




SOOO STABLE
IS **VISALYS**[®] CORE

CORE BUILD-UP & POST CEMENTATION



KETTENBACHDENTAL

Simply intelligent

VISALYS® CORE – EASY TREATMENT

Dual-curing core build-up & root post cementation material

Visalys® Core – especially for difficult and freehand core build-ups.



EXCELLENT STABILITY

- with very good flow properties at the same time



UNIQUE ACTIVE-CONNECT- TECHNOLOGY

- for optimized adhesive bonding
- for use with common adhesives

DENTIN-LIKE GRINDABILITY



FREE OF BISPHENOL A




Distinguished by the
experts “The Dental
Advisor”

VISALYS® CORE – STABLE RESULTS




Adhesive bonding between Visalys® Core and dentin with different adhesives



Visalys® Core achieves the clinically relevant adhesive bond of approx. 10 MPa.

Average removal depth when grinding various core build-up materials



Visalys® Core can be precisely processed similar to dentin.

VISALYS® CORE – THE RANGE

Visalys® Core

Intro pack	White Dentin	5 mL automix syringe, 10 brown mixing tips, 5 intraoral tips, 5 endo tips
Normal pack	White Dentin	2 x 5 mL automix syringe, 20 brown mixing tips, 10 intraoral tips, 10 endo tips
Normal pack	White Dentin	25 mL cartridge, 20 yellow mixing tips, 20 intraoral tips



Kettenbach GmbH & Co. KG
 Im Heerfeld 7
 35713 Eschenburg, Germany
 Phone +49 (0)2774 7050
 info@kettenbach.com
 www.kettenbach-dental.com

022740_4020


KETTENBACHDENTAL
 Simply intelligent