

vista-tec

update!

protective
shields
since 1989

practical and comfortable

elegant design

CE marked PPE protection

glasses and loupes
compatible



polydentia
swiss manufacture 

polydentia.ch




protective
face shields



update! vista-tec line

protective
face shields



REV.01 10.2022

Vista-Tec

the original face shields
Since 1989

Work protected with the stylish and ultra-light Vista-Tec face shield.
Vista-Tec is CE-certified Class II PPE protective face shield, comfortable to wear in long working days. Its easy-to-clean, durable and replaceable shields are made of premium optical quality polycarbonate for relaxed and clear vision.

- ultra-light
- comfortable to wear, adjustable fit
- premium optical quality
- disinfectable, durable and replaceable shields
- class II PPE CE EN166:2001 certified

Vista-Tec kit

REF 5621N 1 frame, 5 shields, 1 headband

Vista-Tec shield

REF 5626N 5 shields



for more
information

Clipon & L-Clipon

protective shields for
thin-frame spectacles and loupes

Vista-Tec Clipon is a practical clip system designed to easily mount a protective shield on prescription glasses.
The clips provide support for both Clipon shields and the larger L-shields, ideal for Galilean loupes. This system is ideal for spectacles and loupes with temple arms of a maximum height of 1cm.

- ideal protection for thin-frame spectacles and loupes
- easy to fit
- durable, disinfectable, and replaceable shields
- premium optical quality

Vista-Tec Clipon kit

REF 5622 1 pair Clipon clips, 5 pcs. shields

Refill

REF 5647 5 pcs. shields
REF 5646 10 pcs. shields
REF 5623 10 pcs. Clipon clips

Vista-Tec L-Clipon kit

REF 5630 1 pair Clipon clips, 5 pcs. L- shields

Refill

REF 5632 10 pcs. L- shields
REF 5623 10 pcs. Clipon clips

Lupon & L-Lupon

protective shields
for loupes

Vista-Tec LupOn protects your lenses and loupes thanks to its wide and durable polycarbonate shields. Available in 2 sizes to best fit prismatic and galilean loupes, the protective shield easily mounts on the frame through a system of clips and oRings. Vista-Tec LupOn is ideal for spectacles and loupes with temple arms over 1cm high

- ideal protection for thick-frame spectacles and loupes
- easy to fit
- premium optical quality
- disinfectable, durable and replaceable shields
- autoclavable clips and oRings

Vista-Tec LupOn kit

REF 5633
1 pair LupOn clips, 1 M-Shield
1 L-Shield, 3 pairs oRings S, M, L

Refill

Vista-Tec LupOn M-Shield
REF 5635 5 M-Shields

Vista-Tec LupOn L-Shield
REF 5636 5 L-Shields

Vista-Tec LupOn oRings
REF 5634 3 pairs oRings S, M, L