

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**els duobond base**

Print date: 16.06.2015

Product code: D300103

Page 1 of 7

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

els duobond base

**1.2. Relevant identified uses of the substance or mixture and uses advised against****1.3. Details of the supplier of the safety data sheet**

|                         |                            |                             |
|-------------------------|----------------------------|-----------------------------|
| Company name:           | Saremco Dental AG          |                             |
| Street:                 | Gewerbestrasse 4           |                             |
| Place:                  | CH-9445 Rebstein           |                             |
| Telephone:              | +41 71 775 80 90           | Telefax: +41 71 775 80 99   |
| e-mail:                 | info@saremco.ch            |                             |
| Contact person:         | Dr. Christoph Evers        | Telephone: +41 71 775 80 60 |
| e-mail:                 | Christoph.evers@saremco.ch |                             |
| Internet:               | www.saremco.ch             |                             |
| Responsible Department: | Vigilance                  |                             |

**1.4. Emergency telephone number:**CH: 145  
Ausland: +41 44 251 51 51**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture**Indications of danger: Xi - Irritant  
R phrases:  
Irritating to eyes, respiratory system and skin.  
May cause sensitisation by skin contact.**GHS classification**Hazard categories:  
Skin corrosion/irritation: Skin Corr. 1B  
Serious eye damage/eye irritation: Eye Dam. 1  
Hazard Statements:  
Causes severe skin burns and eye damage.**2.2. Label elements****Hazardous components which must be listed on the label**

10-Methacryloyl-oxydecyl-dihydrogenphosphat

Signal word: Danger

Pictograms: GHS05

**Hazard statements**

H314 Causes severe skin burns and eye damage.

**Precautionary statements**

|                |                                                                                                      |
|----------------|------------------------------------------------------------------------------------------------------|
| P260           | Do not breathe dust/fume/gas/mist/vapours/spray.                                                     |
| P264           | Wash ... thoroughly after handling.                                                                  |
| P280           | Wear protective gloves/protective clothing/eye protection/face protection.                           |
| P301+P330+P331 | IF SWALLOWED: rinse mouth. Do NOT induce vomiting.                                                   |
| P303+P361+P353 | IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower. |

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## els duobond base

Print date: 16.06.2015

Product code: D300103

Page 2 of 7

|                |                                                                                                                                  |
|----------------|----------------------------------------------------------------------------------------------------------------------------------|
| P363           | Wash contaminated clothing before reuse.                                                                                         |
| P304+P340      | IF INHALED: Remove person to fresh air and keep comfortable for breathing.                                                       |
| P310           | Immediately call a POISON CENTER/doctor.                                                                                         |
| P321           | Specific treatment (see ... on this label).                                                                                      |
| P305+P351+P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P310           | Immediately call a POISON CENTER/doctor.                                                                                         |
| P405           | Store locked up.                                                                                                                 |
| P501           | Dispose of contents/container to ....                                                                                            |

### 2.3. Other hazards

No information available.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures

#### Hazardous components

| EC No      | Chemical name                                                                     | Quantity    |
|------------|-----------------------------------------------------------------------------------|-------------|
| CAS No     | Classification                                                                    |             |
| Index No   | GHS classification                                                                |             |
| REACH No   |                                                                                   |             |
|            | ethox. BisGMA                                                                     | 45 - < 50 % |
| 41637-38-1 | Xi - Irritant R36/37/38-43                                                        |             |
|            |                                                                                   |             |
|            | 10-Methacryloyl-oxydecyl-dihydrogenphosphat                                       | 40 - < 45 % |
| 85590-00-7 |                                                                                   |             |
|            | Skin Corr. 1B; H314                                                               |             |
|            | Campherchinon                                                                     | 1 - < 5 %   |
| 10373-78-1 | Xi - Irritant R36/37/38                                                           |             |
|            |                                                                                   |             |
| 204-881-4  | 2,6-Di-tert.-butyl-p-cresol                                                       | < 1 %       |
| 128-37-0   | Xn - Harmful, N - Dangerous for the environment R22-36/38-51-53                   |             |
|            | Acute Tox. 4, Skin Irrit. 2, Eye Irrit. 2, Aquatic Chronic 2; H302 H315 H319 H411 |             |

Full text of R and H phrases: see Section 16.

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### After inhalation

Provide fresh air. When in doubt or if symptoms are observed, get medical advice.

#### After contact with skin

After contact with skin, wash immediately with polyethylene glycol, followed by plenty of water. Take off immediately all contaminated clothing and wash it before reuse. Medical treatment necessary.

#### After contact with eyes

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then consult an ophthalmologist immediately.

#### After ingestion

Rinse mouth immediately and drink plenty of water.

### 4.2. Most important symptoms and effects, both acute and delayed

No information available.

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

**els duobond base**

Print date: 16.06.2015

Product code: D300103

Page 3 of 7

**4.3. Indication of any immediate medical attention and special treatment needed**

Treat symptomatically.

**SECTION 5: Firefighting measures****5.1. Extinguishing media****Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

**5.2. Special hazards arising from the substance or mixture**

Non-flammable.

**5.3. Advice for firefighters**

Wear a self-contained breathing apparatus and chemical protective clothing. Full protection suit.

**Additional information**

Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures**

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

**6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

**6.3. Methods and material for containment and cleaning up**

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

**6.4. Reference to other sections**

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

**Advice on protection against fire and explosion**

No special fire protection measures are necessary.

**7.2. Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

Keep container tightly closed.

**Advice on storage compatibility**

No special measures are necessary.

**SECTION 8: Exposure controls/personal protection****8.1. Control parameters**

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## els duobond base

Print date: 16.06.2015

Product code: D300103

Page 4 of 7

### Exposure limits (EH40)

| CAS No   | Substance                  | ppm | mg/m <sup>3</sup> | fibres/ml | Category      | Origin |
|----------|----------------------------|-----|-------------------|-----------|---------------|--------|
| 128-37-0 | 2,6-Di-tert-butyl-p-cresol | -   | 10                |           | TWA (8 h)     | WEL    |
|          |                            | -   | -                 |           | STEL (15 min) | WEL    |

### 8.2. Exposure controls



#### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

#### Protective and hygiene measures

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

#### Eye/face protection

Suitable eye protection: goggles.

#### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

#### Skin protection

Wear suitable protective clothing.

#### Respiratory protection

In case of inadequate ventilation wear respiratory protection.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

|                 |           |
|-----------------|-----------|
| Physical state: | liquid    |
| Colour:         | yellow    |
| Odour:          | odourless |

#### Test method

|           |                |
|-----------|----------------|
| pH-Value: | not determined |
|-----------|----------------|

#### Changes in the physical state

|                                          |                |
|------------------------------------------|----------------|
| Melting point:                           | not determined |
| Initial boiling point and boiling range: | not determined |
| Flash point:                             | > 150 °C       |

#### Flammability

|                         |                |
|-------------------------|----------------|
| Solid:                  | not applicable |
| Gas:                    | not applicable |
| Lower explosion limits: | not determined |
| Upper explosion limits: | not determined |

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## els duobond base

Print date: 16.06.2015

Product code: D300103

Page 5 of 7

### Auto-ignition temperature

Solid:

not applicable

Gas:

not applicable

Decomposition temperature:

not determined

### Oxidizing properties

Not oxidizing.

Vapour pressure:

not determined

Density (at 20 °C):

1.3 g/cm<sup>3</sup>

Water solubility:

insoluble

### Solubility in other solvents

not determined

Partition coefficient:

not determined

Vapour density:

not determined

Evaporation rate:

not determined

Solvent content:

0.00 %, Water: 0.03 %

### 9.2. Other information

Solid content:

8.97 %

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

### 10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

### 10.3. Possibility of hazardous reactions

No known hazardous reactions.

### 10.4. Conditions to avoid

none

### 10.5. Incompatible materials

No information available.

### 10.6. Hazardous decomposition products

No known hazardous decomposition products.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

#### Acute toxicity

| CAS No     | Chemical name                               |        |              |         |        |
|------------|---------------------------------------------|--------|--------------|---------|--------|
|            | Exposure routes                             | Method | Dose         | Species | Source |
| 85590-00-7 | 10-Methacryloyl-oxydecyl-dihydrogenphosphat |        |              |         |        |
|            | oral                                        | LD50   | > 2000 mg/kg | Rat     |        |
| 128-37-0   | 2,6-Di-tert.-butyl-p-cresol                 |        |              |         |        |
|            | oral                                        | LD50   | 890.0 mg/kg  | Rat     |        |

#### Additional information on tests

This mixture is classified as hazardous according to 1999/45/EC.

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**els duobond base**

Print date: 16.06.2015

Product code: D300103

Page 6 of 7

**SECTION 12: Ecological information****12.1. Toxicity**

Aquatic toxicity: The classification criteria for this hazard class are not met by definition.

**12.2. Persistence and degradability**

The product has not been tested.

**12.3. Bioaccumulative potential**

The product has not been tested.

**12.4. Mobility in soil**

The product has not been tested.

**12.5. Results of PBT and vPvB assessment**

The product has not been tested.

**12.6. Other adverse effects**

No information available.

**Further information**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil.

**SECTION 13: Disposal considerations****13.1. Waste treatment methods****Advice on disposal**

Do not allow to enter into surface water or drains. Do not allow to enter into soil/subsoil. Dispose of waste according to applicable legislation.

**Contaminated packaging**

Non-contaminated packages may be recycled. Handle contaminated packages in the same way as the substance itself.

**SECTION 14: Transport information****Land transport (ADR/RID)****Other applicable information (land transport)**

No dangerous good in sense of these transport regulations.

**Marine transport (IMDG)****Other applicable information (marine transport)**

No dangerous good in sense of these transport regulations.

**Air transport (ICAO)****Other applicable information (air transport)**

No dangerous good in sense of these transport regulations.

**14.6. Special precautions for user**

No information available.

**14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code**

not applicable

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National regulatory information**

Employment restrictions:

Observe employment restrictions for young people.

Water contaminating class (D):

3 - highly water contaminating

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**els duobond base**

Print date: 16.06.2015

Product code: D300103

Page 7 of 7

**15.2. Chemical safety assessment**

Chemical safety assessments for substances in this mixture were not carried out.

**SECTION 16: Other information****Abbreviations and acronyms**

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road )  
IMDG: International Maritime Code for Dangerous Goods  
IATA: International Air Transport Association  
GHS: Globally Harmonized System of Classification and Labelling of Chemicals  
EINECS: European Inventory of Existing Commercial Chemical Substances  
ELINCS: European List of Notified Chemical Substances  
CAS: Chemical Abstracts Service  
LC50: Lethal concentration, 50%  
LD50: Lethal dose, 50%

**Full text of R phrases referred to under Sections 2 and 3**

22 Harmful if swallowed.  
36/37/38 Irritating to eyes, respiratory system and skin.  
36/38 Irritating to eyes and skin.  
43 May cause sensitisation by skin contact.  
51 Toxic to aquatic organisms.  
53 May cause long-term adverse effects in the aquatic environment.

**Full text of H statements referred to under Sections 2 and 3**

H302 Harmful if swallowed.  
H314 Causes severe skin burns and eye damage.  
H315 Causes skin irritation.  
H319 Causes serious eye irritation.  
H411 Toxic to aquatic life with long lasting effects.

**Further Information**

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

**els duobond Catalyst**

Print date: 16.06.2015

Product code: D300105

Page 1 of 8

**SECTION 1: Identification of the substance/mixture and of the company/undertaking****1.1. Product identifier**

els duobond Catalyst

**1.2. Relevant identified uses of the substance or mixture and uses advised against****1.3. Details of the supplier of the safety data sheet**

|                         |                            |                             |
|-------------------------|----------------------------|-----------------------------|
| Company name:           | Saremco Dental AG          |                             |
| Street:                 | Gewerbestrasse 4           |                             |
| Place:                  | CH-9445 Rebstein           |                             |
| Telephone:              | +41 71 775 80 90           | Telefax: +41 71 775 80 99   |
| e-mail:                 | info@saremco.ch            |                             |
| Contact person:         | Dr. Christoph Evers        | Telephone: +41 71 775 80 60 |
| e-mail:                 | Christoph.evers@saremco.ch |                             |
| Internet:               | www.saremco.ch             |                             |
| Responsible Department: | Vigilance                  |                             |

**1.4. Emergency telephone** CH: 145  
**number:** Ausland: +41 44 251 51 51

**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture**

Indications of danger: Xi - Irritant  
R phrases:  
Irritating to eyes.  
Repeated exposure may cause skin dryness or cracking.  
Vapours may cause drowsiness and dizziness.

**GHS classification**

Hazard categories:  
Flammable liquid: Flam. Liq. 3  
Serious eye damage/eye irritation: Eye Irrit. 2  
Specific target organ toxicity - single exposure: STOT SE 3  
Hazard Statements:  
Flammable liquid and vapour.  
Causes serious eye irritation.  
May cause drowsiness or dizziness.

**2.2. Label elements****Hazardous components which must be listed on the label**

acetone; propan-2-one; propanone

Signal word: Warning  
Pictograms: GHS02-GHS07

**Hazard statements**

|      |                                    |
|------|------------------------------------|
| H226 | Flammable liquid and vapour.       |
| H319 | Causes serious eye irritation.     |
| H336 | May cause drowsiness or dizziness. |

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## els duobond Catalyst

Print date: 16.06.2015

Product code: D300105

Page 2 of 8

### Precautionary statements

|                |                                                                                                                                  |
|----------------|----------------------------------------------------------------------------------------------------------------------------------|
| P210           | Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.                                   |
| P233           | Keep container tightly closed.                                                                                                   |
| P240           | Ground/bond container and receiving equipment.                                                                                   |
| P241           | Use explosion-proof electrical/ventilating/lighting equipment.                                                                   |
| P242           | Use only non-sparking tools.                                                                                                     |
| P243           | Take precautionary measures against static discharge.                                                                            |
| P261           | Avoid breathing dust/fume/gas/mist/vapours/spray.                                                                                |
| P264           | Wash ... thoroughly after handling.                                                                                              |
| P271           | Use only outdoors or in a well-ventilated area.                                                                                  |
| P280           | Wear protective gloves/protective clothing/eye protection/face protection.                                                       |
| P303+P361+P353 | IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower.                             |
| P304+P340      | IF INHALED: Remove person to fresh air and keep comfortable for breathing.                                                       |
| P312           | Call a POISON CENTER/doctor if you feel unwell.                                                                                  |
| P305+P351+P338 | IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. |
| P337+P313      | If eye irritation persists: Get medical advice/attention.                                                                        |
| P370+P378      | In case of fire: Use ... to extinguish.                                                                                          |
| P235           | Keep cool.                                                                                                                       |
| P403+P233      | Store in a well-ventilated place. Keep container tightly closed.                                                                 |
| P405           | Store locked up.                                                                                                                 |
| P501           | Dispose of contents/container to ....                                                                                            |

### 2.3. Other hazards

No information available.

## SECTION 3: Composition/information on ingredients

### 3.2. Mixtures

#### Hazardous components

| EC No        | Chemical name                                         | Quantity    |
|--------------|-------------------------------------------------------|-------------|
| CAS No       | Classification                                        |             |
| Index No     | GHS classification                                    |             |
| REACH No     |                                                       |             |
| 200-662-2    | acetone; propan-2-one; propanone                      |             |
| 67-64-1      | F - Highly flammable, Xi - Irritant R11-36-66-67      | 65 - < 70 % |
| 606-001-00-8 | Flam. Liq. 2, Eye Irrit. 2, STOT SE 3; H225 H319 H336 |             |

Full text of R and H phrases: see Section 16.

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### General information

First aider: Pay attention to self-protection! Move victim out of danger zone.

#### After inhalation

Provide fresh air. In case of accident or unwellness, seek medical advice immediately (show directions for use or safety data sheet if possible).

#### After contact with skin

Wash with plenty of water. Take off contaminated clothing and wash it before reuse.

#### After contact with eyes

After contact with the eyes, rinse with water with the eyelids open for a sufficient length of time, then

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**els duobond Catalyst**

Print date: 16.06.2015

Product code: D300105

Page 3 of 8

consult an ophthalmologist immediately.

**After ingestion**

Rinse mouth immediately and drink plenty of water.

**4.2. Most important symptoms and effects, both acute and delayed**

No information available.

**4.3. Indication of any immediate medical attention and special treatment needed**

Treat symptomatically.

**SECTION 5: Firefighting measures****5.1. Extinguishing media****Suitable extinguishing media**

Co-ordinate fire-fighting measures to the fire surroundings.

**5.2. Special hazards arising from the substance or mixture**

Non-flammable. Vapours can form explosive mixtures with air.

**5.3. Advice for firefighters**

In case of fire: Wear self-contained breathing apparatus.

**Additional information**

Use water spray jet to protect personnel and to cool endangered containers. Suppress gases/vapours/mists with water spray jet. Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

**SECTION 6: Accidental release measures****6.1. Personal precautions, protective equipment and emergency procedures**

Provide adequate ventilation. Do not breathe gas/fumes/vapour/spray. Avoid contact with skin, eyes and clothes. Use personal protection equipment.

**6.2. Environmental precautions**

Do not allow to enter into surface water or drains.

**6.3. Methods and material for containment and cleaning up**

Absorb with liquid-binding material (e.g. sand, diatomaceous earth, acid- or universal binding agents). Treat the recovered material as prescribed in the section on waste disposal.

**6.4. Reference to other sections**

Safe handling: see section 7

Personal protection equipment: see section 8

Disposal: see section 13

**SECTION 7: Handling and storage****7.1. Precautions for safe handling****Advice on safe handling**

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

**Advice on protection against fire and explosion**

No special fire protection measures are necessary.

**7.2. Conditions for safe storage, including any incompatibilities****Requirements for storage rooms and vessels**

Keep container tightly closed.

**Advice on storage compatibility**

No special measures are necessary.

**SECTION 8: Exposure controls/personal protection**

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## els duobond Catalyst

Print date: 16.06.2015

Product code: D300105

Page 4 of 8

### 8.1. Control parameters

#### Exposure limits (EH40)

| CAS No  | Substance | ppm  | mg/m <sup>3</sup> | fibres/ml | Category      | Origin |
|---------|-----------|------|-------------------|-----------|---------------|--------|
| 67-64-1 | Acetone   | 500  | 1210              |           | TWA (8 h)     | WEL    |
|         |           | 1500 | 3620              |           | STEL (15 min) | WEL    |

### 8.2. Exposure controls



#### Appropriate engineering controls

If handled uncovered, arrangements with local exhaust ventilation have to be used. Do not breathe gas/fumes/vapour/spray.

#### Protective and hygiene measures

Remove contaminated, saturated clothing immediately. Draw up and observe skin protection programme. Wash hands and face before breaks and after work and take a shower if necessary. When using do not eat or drink.

#### Eye/face protection

Suitable eye protection: goggles.

#### Hand protection

When handling with chemical substances, protective gloves must be worn with the CE-label including the four control digits. The quality of the protective gloves resistant to chemicals must be chosen as a function of the specific working place concentration and quantity of hazardous substances. For special purposes, it is recommended to check the resistance to chemicals of the protective gloves mentioned above together with the supplier of these gloves.

#### Skin protection

Wear suitable protective clothing.

#### Respiratory protection

In case of inadequate ventilation wear respiratory protection.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

Physical state: liquid  
 Colour: yellow  
 Odour: like: Acetone

#### Test method

pH-Value: not determined

#### Changes in the physical state

Melting point: not determined

Initial boiling point and boiling range: 56 °C

Flash point: < -20 °C

#### Flammability

Solid: not applicable

Gas: not applicable

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## els duobond Catalyst

Print date: 16.06.2015

Product code: D300105

Page 5 of 8

Lower explosion limits: 2,5 vol. %  
 Upper explosion limits: 15 vol. %  
 Ignition temperature: 400 °C

### Auto-ignition temperature

Solid: not applicable  
 Gas: not applicable  
 Decomposition temperature: not determined

### Oxidizing properties

Not oxidizing.

Vapour pressure: 246 hPa  
 (at 20 °C)

Vapour pressure: 814 hPa  
 (at 50 °C)

Density: not determined

Water solubility: easily soluble

### Solubility in other solvents

not determined

Partition coefficient: not determined

Vapour density: not determined

Evaporation rate: not determined

Solvent content: 21.77 %, Water: 31.90 %

### 9.2. Other information

Solid content: 3.11 %

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

No hazardous reaction when handled and stored according to provisions.

### 10.2. Chemical stability

The product is stable under storage at normal ambient temperatures.

### 10.3. Possibility of hazardous reactions

No known hazardous reactions.

### 10.4. Conditions to avoid

none

### 10.5. Incompatible materials

No information available.

### 10.6. Hazardous decomposition products

No known hazardous decomposition products.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**els duobond Catalyst**

Print date: 16.06.2015

Product code: D300105

Page 6 of 8

**Acute toxicity**

| CAS No  | Chemical name                    |        |             |           |        |
|---------|----------------------------------|--------|-------------|-----------|--------|
|         | Exposure routes                  | Method | Dose        | Species   | Source |
| 67-64-1 | acetone; propan-2-one; propanone |        |             |           |        |
|         | oral                             | LD50   | 5800 mg/kg  | Ratte     | RTECS  |
|         | dermal                           | LD50   | 20000 mg/kg | Kaninchen | IUCLID |
|         | inhalative (4 h) vapour          | LC50   | 76 mg/l     | Ratte     |        |

**Additional information on tests**

This mixture is classified as hazardous according to 1999/45/EC.

**SECTION 12: Ecological information****12.1. Toxicity**

Aquatic toxicity: The classification criteria for this hazard class are not met by definition.

| CAS No  | Chemical name                    |        |           |           |                      |        |
|---------|----------------------------------|--------|-----------|-----------|----------------------|--------|
|         | Aquatic toxicity                 | Method | Dose      | [h]   [d] | Species              | Source |
| 67-64-1 | acetone; propan-2-one; propanone |        |           |           |                      |        |
|         | Acute fish toxicity              | LC50   | 5540 mg/l | 96 h      | Onchorhynchus mykiss |        |
|         | Acute crustacea toxicity         | EC50   | 6100 mg/l | 48 h      | Daphnia magna        |        |

**12.2. Persistence and degradability**

The product has not been tested.

**12.3. Bioaccumulative potential**

The product has not been tested.

**Partition coefficient n-octanol/water**

| CAS No  | Chemical name                    | Log Pow |
|---------|----------------------------------|---------|
| 67-64-1 | acetone; propan-2-one; propanone | -0,24   |

**12.4. Mobility in soil**

The product has not been tested.

**12.5. Results of PBT and vPvB assessment**

The product has not been tested.

**12.6. Other adverse effects**

No information available.

**Further information**

No special environmental measures are necessary.

**SECTION 13: Disposal considerations****13.1. Waste treatment methods****Advice on disposal**

Do not allow to enter into surface water or drains. Dispose of waste according to applicable legislation.

**Contaminated packaging**

Wash with plenty of water. Completely emptied packages can be recycled.

**SECTION 14: Transport information****Land transport (ADR/RID)**

# Safety Data Sheet

according to Regulation (EC) No 1907/2006

## els duobond Catalyst

Print date: 16.06.2015

Product code: D300105

Page 7 of 8

### Other applicable information (land transport)

No dangerous good in sense of these transport regulations.

### Marine transport (IMDG)

#### Other applicable information (marine transport)

No dangerous good in sense of these transport regulations.

### Air transport (ICAO)

#### Other applicable information (air transport)

No dangerous good in sense of these transport regulations.

### 14.6. Special precautions for user

No information available.

### 14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code

not applicable

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### EU regulatory information

1999/13/EC (VOC): 69 %

#### National regulatory information

Employment restrictions: Observe employment restrictions for young people.

Water contaminating class (D): 1 - slightly water contaminating

### 15.2. Chemical safety assessment

Chemical safety assessments for substances in this mixture were not carried out.

## SECTION 16: Other information

### Abbreviations and acronyms

ADR: Accord européen sur le transport des marchandises dangereuses par Route  
(European Agreement concerning the International Carriage of Dangerous Goods by Road)

IMDG: International Maritime Code for Dangerous Goods

IATA: International Air Transport Association

GHS: Globally Harmonized System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

ELINCS: European List of Notified Chemical Substances

CAS: Chemical Abstracts Service

LC50: Lethal concentration, 50%

LD50: Lethal dose, 50%

### Full text of R phrases referred to under Sections 2 and 3

- 11 Highly flammable.
- 36 Irritating to eyes.
- 66 Repeated exposure may cause skin dryness or cracking.
- 67 Vapours may cause drowsiness and dizziness.

### Full text of H statements referred to under Sections 2 and 3

- H225 Highly flammable liquid and vapour.
- H226 Flammable liquid and vapour.
- H319 Causes serious eye irritation.
- H336 May cause drowsiness or dizziness.

### Further Information

The above information describes exclusively the safety requirements of the product and is based on our present-day knowledge. The information is intended to give you advice about the safe handling of

**Safety Data Sheet**

according to Regulation (EC) No 1907/2006

**els duobond Catalyst**

Print date: 16.06.2015

Product code: D300105

Page 8 of 8

the product named in this safety data sheet, for storage, processing, transport and disposal. The information cannot be transferred to other products. In the case of mixing the product with other products or in the case of processing, the information on this safety data sheet is not necessarily valid for the new made-up material.

---

*(The data for the hazardous ingredients were taken respectively from the last version of the sub-contractor's safety data sheet.)*