

neodisher Septo Active

Version: 2 / GB

Replaces Version: 1 /
GB

Date revised: 20.05.2016

Print date: 28.06.17

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

neodisher Septo Active

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/preparation

disinfectants, Washing and cleaning products (including solvent based products)

1.3. Details of the supplier of the safety data sheet

Address:

Chemische Fabrik Dr. Weigert GmbH & Co. KG
Mühlenhagen 85
D-20539 Hamburg
Telephone no. +49 40 789 60 0
Fax no. +49 40 789 60 120
www.drweigert.com

E-mail address of person responsible for this SDS:

sida@drweigert.de

1.4. Emergency telephone number

GBK/ Infotrac: (USA domestic) 1 800 535 5053 or international +1 352 323 3500

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification (Regulation (EC) No. 1272/2008)

Classification (Regulation (EC) No. 1272/2008)

Acute Tox. 4	H302
Eye Dam. 1	H318

2.2. Label elements

Labelling according to regulation (EC) No 1272/2008

Hazard pictograms



Signal word

Danger

Hazard statements

H302	Harmful if swallowed.
H318	Causes serious eye damage.

Precautionary statements

P280	Wear protective gloves/protective clothing/eye protection/face protection.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor.

neodisher Septo Active

Version: 2 / GB

Replaces Version: 1 /
GB

Date revised: 20.05.2016

Print date: 28.06.17

Dispose only when container is empty and closed. For disposal of product residues, refer to Safety Data Sheet.

Hazardous component(s) to be indicated on label (Regulation (EC) No. 1272/2008)
contains sodium carbonate peroxyhydrate

2.3. Other hazards

No special hazards have to be mentioned.

SECTION 3: Composition/information on ingredients

3.2. Mixtures

Hazardous ingredients

sodium carbonate peroxyhydrate

CAS No.	15630-89-4			
EINECS no.	239-707-6			
Registration no.	01-2119457268-30			
Concentration	>= 25	<	50	%
Classification (Regulation (EC) No. 1272/2008)				
	Ox. Sol. 3		H272	
	Acute Tox. 4		H302	
	Eye Dam. 1		H318	

Concentration limits (Regulation (EC) No. 1272/2008)

Eye Dam. 1	H318	>= 25 %
Eye Irrit. 2	H319	>= 7,5 < 25 %

sodium carbonate

CAS No.	497-19-8			
EINECS no.	207-838-8			
Registration no.	01-2119485498-19			
Concentration	>= 1	<	10	%
Classification (Regulation (EC) No. 1272/2008)				
	Eye Irrit. 2		H319	

citric acid, anhydrous

CAS No.	77-92-9			
EINECS no.	201-069-1			
Registration no.	01-2119457026-42			
Concentration	>= 10	<	25	%
Classification (Regulation (EC) No. 1272/2008)				
	Eye Irrit. 2		H319	

fatty alcohols, ethoxylated, propoxylated

CAS No.	68439-51-0			
Concentration	>= 1	<	5	%
Classification (Regulation (EC) No. 1272/2008)				
	Aquatic Chronic 3		H412	

SECTION 4: First aid measures

4.1. Description of first aid measures

General information

Remove contaminated clothing immediately and dispose of safely. In any case show the physician the Safety Data Sheet.

neodisher Septo Active

Version: 2 / GB

Replaces Version: 1 /
GB

Date revised: 20.05.2016

Print date: 28.06.17

After inhalation

Ensure supply of fresh air. When dust is intensively inhaled, seek medical help immediately.

After skin contact

Wash off immediately with soap and water. Consult a doctor if skin irritation persists.

After eye contact

Separate eyelids, wash the eyes thoroughly with water (15 min.). Summon a doctor immediately.

After ingestion

If swallowed, seek medical advice immediately and show this container or label. Rinse mouth thoroughly with water. Let plenty of water be drunk in small gulps. Do not induce vomiting.

Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

4.3. Indication of any immediate medical attention and special treatment needed

Hints for the physician / hazards

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Product itself is non-combustible; adapt fire extinguishing measures to surrounding areas.

Non suitable extinguishing media

Full water jet

5.2. Special hazards arising from the substance or mixture

In case of combustion evolution of dangerous gases possible.

5.3. Advice for firefighters

Special protective equipment for fire-fighting

Do not inhale explosion and/or combustion gases. In case of combustion use a suitable breathing apparatus.

Other information

Collect contaminated fire-fighting water separately, must not be discharged into the drains. Fire residues and contaminated fire-fighting water must be disposed of in accordance with the local regulations.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Avoid contact with skin, eyes and clothing.

6.2. Environmental precautions

Do not discharge into the drains/surface waters/groundwater. Knock down dust with water spray jet.

6.3. Methods and material for containment and cleaning up

Pick up mechanically. Dispose of absorbed material in accordance with the regulations.

6.4. Reference to other sections

Refer to protective measures listed in Sections 7 and 8.

SECTION 7: Handling and storage

neodisher Septo Active

Version: 2 / GB

Replaces Version: 1 /
GB

Date revised: 20.05.2016

Print date: 28.06.17

7.1. Precautions for safe handling

Advice on safe handling

Avoid the formation and deposition of dust. Keep container tightly closed.

Advice on protection against fire and explosion

Keep away from combustible material.

7.2. Conditions for safe storage, including any incompatibilities

Recommended storage temperature

Value $\geq$ 0 $<$ 25 °C

Requirements for storage rooms and vessels

Keep in original packaging, tightly closed.

Hints on storage assembly

Do not store with combustible materials.

Storage class according to TRGS 510

Storage class according to TRGS 510 13 Non- combustible solids

Further information on storage conditions

Protect from heat and direct sunlight.

7.3. Specific end use(s)

no data

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Other information

Contains no substances with occupational exposure limit values.

8.2. Exposure controls

General protective and hygiene measures

Do not inhale dust/fumes/aerosols. Avoid contact with skin and eyes. Do not eat, drink or smoke during work time. Wash hands before breaks and after work. Clean skin thoroughly after work; apply skin cream.

Respiratory protection

Use breathing apparatus in dust-laden atmosphere. Dust mask

Hand protection

Chemical resistant gloves (EN 374)

Use Permanent hand contact
Appropriate Material neoprene

Material thickness $\geq$ 0,65 mm
Breakthrough time $>$ 480 min

Appropriate Material butyl
Material thickness $\geq$ 0,7 mm
Breakthrough time $>$ 480 min

Appropriate Material nitrile
Material thickness $\geq$ 0,4 mm
Breakthrough time $>$ 480 min

Use Short-term hand contact
Appropriate Material nitrile
Material thickness $\geq$ 0,11 mm

Eye protection

neodisher Septo Active

Version: 2 / GB

Replaces Version: 1 /
GB

Date revised: 20.05.2016

Print date: 28.06.17

Safety glasses with side protection shield (EN 166)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Form	solid		
Colour	white		
Odour	characteristic		
Odour threshold			
Remarks	not determined		
pH value			
Value	7,6	to	7,9
Concentration/H ₂ O	2	%	
Temperature	20	°C	
Melting point			
Remarks	not determined		
Freezing point			
Remarks	not determined		
Initial boiling point and boiling range			
Remarks	not determined		
Flash point			
Remarks	Not applicable		
Evaporation rate (ether = 1) :			
Remarks	not determined		
Flammability (solid, gas)			
evaluation	not determined		
Upper/lower flammability or explosive limits			
Remarks	not determined		
Vapour pressure			
Remarks	not determined		
Vapour density			
Remarks	not determined		
Density			
Remarks	not determined		
Solubility in water			
Remarks	not determined		
Solubility(ies)			
Remarks	not determined		
Partition coefficient: n-octanol/water			
Remarks	not determined		
Ignition temperature			
Remarks	not determined		
Decomposition temperature			
Remarks	not determined		
Viscosity			
Remarks	Not applicable		

neodisher Septo Active

Version: 2 / GB

Replaces Version: 1 /
GB

Date revised: 20.05.2016

Print date: 28.06.17

Explosive properties

evaluation not determined

Oxidising properties

evaluation oxidizing

9.2. Other information

Bulk density

Value appr. 0,80 kg/l

Burning test

Remarks Not applicable

Other information

None known

SECTION 10: Stability and reactivity

10.1. Reactivity

No hazardous reactions when stored and handled according to prescribed instructions.

10.2. Chemical stability

No hazardous reactions known.

10.3. Possibility of hazardous reactions

No hazardous reactions known.

10.4. Conditions to avoid

Protect from heat and direct sunlight.

Decomposition temperature

Remarks not determined

10.5. Incompatible materials

Acids, Reducing agents

10.6. Hazardous decomposition products

Oxygen, Irritant gases/vapours

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute oral toxicity

Species	rat	
LD50	appr. 2000	mg/kg
Method	calculated value (Regulation (EC) No. 1272/2008)	

Acute oral toxicity (Components)

sodium carbonate peroxyhydrate

Species	rat	
LD50	1034	mg/kg
Method	Value taken from the literature	

Acute dermal toxicity

Remarks Based on available data, the classification criteria are not met.

Acute dermal toxicity (Components)

sodium carbonate peroxyhydrate

Species	rabbit	
LD50	> 2000	mg/kg

neodisher Septo Active

Version: 2 / GB

Replaces Version: 1 /
GB

Date revised: 20.05.2016

Print date: 28.06.17

Method OECD 402

Acute inhalational toxicity

Remarks Based on available data, the classification criteria are not met.

Skin corrosion/irritation

Remarks Based on available data, the classification criteria are not met.

Serious eye damage/irritation

evaluation irritant - risk of serious damage to eyes

Serious eye damage/irritation (Components)

sodium carbonate peroxyhydrate

Species rabbit eye
evaluation irritant - risk of serious damage to eyes
Method OECD 405

Sensitization

Remarks Based on available data, the classification criteria are not met.

Sensitization (Components)

sodium carbonate peroxyhydrate

Route of exposure dermal
Species guinea pig
evaluation non-sensitizing
Method OECD 406

Subacute, subchronic, chronic toxicity

Remarks Based on available data, the classification criteria are not met.

Mutagenicity

Remarks Based on available data, the classification criteria are not met.

Reproductive toxicity

Remarks Based on available data, the classification criteria are not met.

Carcinogenicity

Remarks Based on available data, the classification criteria are not met.

Specific Target Organ Toxicity (STOT)

Remarks not determined

Aspiration hazard

No special hazards have to be mentioned.

Experience in practice

Inhalation of dusts may irritate the respiratory tract.

Other information

There is no data available on the product apart from the information given in this subsection.

SECTION 12: Ecological information ***

12.1. Toxicity

General information

not determined

Fish toxicity (Components)

sodium carbonate peroxyhydrate

Species Fathead minnow (*Pimephales promelas*)
LC50 70,7 mg/l
Duration of exposure 96 h

neodisher Septo Active

Version: 2 / GB

Replaces Version: 1 /
GB

Date revised: 20.05.2016

Print date: 28.06.17

Daphnia toxicity

Reference substance sodium carbonate peroxyhydrate

Daphnia toxicity (Components)

sodium carbonate peroxyhydrate

Species	Daphnia pulex	
EC50	4,9	mg/l
Duration of exposure	48	h

sodium carbonate peroxyhydrate

Species	Daphnia pulex	
NOEC	2	mg/l
Duration of exposure	48	h

Bacteria toxicity (Components)

sodium carbonate peroxyhydrate

Species	activated sludge	
EC50	466	mg/l
Duration of exposure	30	min

12.2. Persistence and degradability

General information

not determined

12.3. Bioaccumulative potential

General information

not determined

Partition coefficient: n-octanol/water

Remarks not determined

12.4. Mobility in soil

General information

not determined

12.5. Results of PBT and vPvB assessment

General information

not determined

Evaluation of persistence and bioaccumulation potential

The product contains no PBT or vPvB substances.

12.6. Other adverse effects

General information

not determined

General information / ecology

Do not allow to enter soil, waterways or waste water canal.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal recommendations for the product

Allocation of a waste code number, according to the European Waste Catalogue (EWC), should be carried out in agreement with the regional waste disposal company.

Disposal recommendations for packaging

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal

neodisher Septo Active

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 20.05.2016

Print date: 28.06.17

company.

SECTION 14: Transport information

	Land transport ADR/RID	Marine transport IMDG/GGVSee	Air transport ICAO/IATA
14.1. UN number	Non-dangerous goods	The product does not constitute a hazardous substance in sea transport.	The product does not constitute a hazardous substance in air transport.

Information for all modes of transport

14.6. Special precautions for user

See Sections 6 to 8

Other information

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

Ingredients (Regulation (EC) No 648/2004)

30 % and more:

oxygen-based bleaching agents

less than 5 %:

non-ionic surfactants, phosphates

Water Hazard Class (Germany)

Water Hazard Class (Germany) WGK 1

Remarks Classification according to Annex 4 VwVwS

VOC

VOC (EU) 0 %

Other information

The product does not contain substances of very high concern (SVHC).

15.2. Chemical safety assessment

For this preparation a chemical safety assessment has not been carried out.

SECTION 16: Other information

Hazard statements listed in Chapter 3

H272	May intensify fire; oxidiser.
H302	Harmful if swallowed.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H412	Harmful to aquatic life with long lasting effects.

CLP categories listed in Chapter 3

Acute Tox. 4	Acute toxicity, Category 4
Aquatic Chronic 3	Hazardous to the aquatic environment, chronic, Category 3
Eye Dam. 1	Serious eye damage, Category 1
Eye Irrit. 2	Eye irritation, Category 2
Ox. Sol. 3	Oxidising solid, Category 3

neodisher Septo Active

Version: 2 / GB

Replaces Version: 1 /
GB

Date revised: 20.05.2016

Print date: 28.06.17

Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: ***
This information is based on our present state of knowledge. However, it should not constitute a
guarantee for any specific product properties and shall not establish a legally valid relationship.