

neodisher Septo Fin

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 22.12.2017

Print date: 22.12.17

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

neodisher Septo Fin

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/preparation

disinfectants

1.3. Details of the supplier of the safety data sheet

Address:

Chemische Fabrik Dr. Weigert GmbH & Co. KG
Mühlenhagen 85
D-20539 Hamburg
Telephone no. +49 40 789 60 0
Fax no. +49 40 789 60 120
www.drweigert.com

E-mail address of person responsible for this SDS:

sida@drweigert.de

1.4. Emergency telephone number

GBK/ Infotrac: (USA domestic) +1 800 535 5053 or international +1 352 323 3500

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification (Regulation (EC) No. 1272/2008)

Classification (Regulation (EC) No. 1272/2008)

Acute Tox. 4	H302
Acute Tox. 4	H332
Skin Sens. 1	H317
Resp. Sens. 1	H334
Skin Corr. 1B	H314
Eye Dam. 1	H318
STOT SE 3	H335
Aquatic Chronic 3	H412

Route of exposure: oral

Route of exposure: inhalative

2.2. Label elements

Labelling according to regulation (EC) No 1272/2008

Hazard pictograms



Signal word

Danger

Hazard statements

H302+H332
H314

Harmful if swallowed or if inhaled.
Causes severe skin burns and eye damage.

neodisher Septo Fin

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 22.12.2017

Print date: 22.12.17

H317	May cause an allergic skin reaction.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H412	Harmful to aquatic life with long lasting effects.
EUH071	Corrosive to the respiratory tract.

Precautionary statements

P260	Do not breathe dust/fume/gas/mist/vapours/spray.
P273	Avoid release to the environment.
P280	Wear protective gloves/protective clothing/eye protection/face protection.
P303+P361+P353	IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
P304+P340	IF INHALED: Remove person to fresh air and keep comfortable for breathing.
P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310	Immediately call a POISON CENTER or doctor. Dispose only when container is empty and closed. For disposal of product residues, refer to Safety Data Sheet.

Hazardous component(s) to be indicated on label (Regulation (EC) No. 1272/2008)

contains glutaral

2.3. Other hazards

No special hazards have to be mentioned.

SECTION 3: Composition/information on ingredients

3.2. Mixtures

Hazardous ingredients

alkylether carboxylic acid

CAS No.	53563-70-5			
Concentration	>= 1	< 10	%	
Classification (Regulation (EC) No. 1272/2008)	Eye Dam. 1	H318		
	Skin Irrit. 2	H315		

glutaral

CAS No.	111-30-8			
EINECS no.	203-856-5			
Registration no.	01-2119455549-26			
Concentration	>= 10	< 25	%	
Classification (Regulation (EC) No. 1272/2008)	Skin Sens. 1A	H317		
	Aquatic Acute 1	H400		
	Resp. Sens. 1	H334		
	Skin Corr. 1B	H314		
	Acute Tox. 3	H301		Route of exposure: oral
	Acute Tox. 2	H330		Route of exposure: inhalative
	Aquatic Chronic 2	H411		
	STOT SE 3	H335		

Concentration limits (Regulation (EC) No. 1272/2008)

STOT SE 3 H335 >= 0,5 < 5

fatty alcohols, ethoxylated

CAS No.	68439-50-9			
Concentration	< 1	%		
Classification (Regulation (EC) No. 1272/2008)	Acute Tox. 4	H302		

neodisher Septo Fin

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 22.12.2017

Print date: 22.12.17

Eye Dam. 1 H318
Aquatic Chronic 3 H412

Concentration limits (Regulation (EC) No. 1272/2008)

Eye Dam. 1 H318 > 10 %
Eye Irrit. 2 H319 > 1 <= 10 %

1,2,3-benzotriazole

CAS No. 95-14-7

EINECS no. 202-394-1

Concentration < 1 %

Classification (Regulation (EC) No. 1272/2008)

Acute Tox. 4 H302 Route of exposure: oral
Eye Irrit. 2 H319
Aquatic Chronic 2 H411

Other information

Complete text of hazard statements in chapter 16

SECTION 4: First aid measures

4.1. Description of first aid measures

General information

Remove contaminated, soaked clothing immediately and dispose of safely. Adhere to personal protective measures when giving first aid. Clean body thoroughly (bath, shower). In any case show the physician the Safety Data Sheet.

After inhalation

Ensure supply of fresh air. Remove affected person from danger area. Seek medical advice immediately.

After skin contact

Wash off immediately with soap and water. Seek medical advice immediately.

After eye contact

Separate eyelids, wash the eyes thoroughly with water (15 min.). Take medical treatment.

After ingestion

Call in a physician immediately and show him the Safety Data Sheet. Rinse mouth thoroughly with water. Let plenty of water be drunk in small gulps. Do not induce vomiting.

Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

4.3. Indication of any immediate medical attention and special treatment needed

Hints for the physician / hazards

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media

Dry powder, Foam, Water spray jet

Non suitable extinguishing media

Full water jet

neodisher Septo Fin

Version: 3 / GB

Replaces Version: 2 / GB

Date revised: 22.12.2017

Print date: 22.12.17

5.2. Special hazards arising from the substance or mixture

In case of combustion evolution of dangerous gases possible.

5.3. Advice for firefighters

Special protective equipment for fire-fighting

Do not inhale explosion and/or combustion gases. In case of combustion use a suitable breathing apparatus. Wear full protective suit.

Other information

Collect contaminated fire-fighting water separately, must not be discharged into the drains. Fire residues and contaminated fire-fighting water must be disposed of in accordance with the local regulations.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Use breathing apparatus if exposed to vapours/dust/aerosol. Avoid contact with skin, eyes and clothing. Refer to protective measures listed in Sections 7 and 8.

6.2. Environmental precautions

Prevent spread over a wide area (e.g. by containment or oil barriers). Do not discharge into the drains/surface waters/groundwater. Do not discharge into the subsoil/soil. Retain and dispose of contaminated wash water. In case of gas escape or of entry into waterways, soil or drains, inform the responsible authorities.

6.3. Methods and material for containment and cleaning up

Pick up with absorbent material. Clean contaminated floors and objects thoroughly with water and detergents, observing environmental regulations. Containers in which spilt substance has been collected must be adequately labelled. Dispose of absorbed material in accordance with the regulations.

6.4. Reference to other sections

Refer to protective measures listed in Sections 7 and 8.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Advice on safe handling

Avoid formation of aerosols. Perform filling operations only at stations with exhaust ventilation facilities. Provide suitable exhaust ventilation at the processing machines. If workplace limits are exceeded, a respiratory protection approved for this particular job must be worn. Keep container tightly closed.

7.2. Conditions for safe storage, including any incompatibilities

Recommended storage temperature

Value > 0 < 30 °C

Requirements for storage rooms and vessels

Keep in original packaging, tightly closed. Storage rooms must be properly ventilated. Containers which are opened must be carefully resealed and kept upright to prevent leakage. Provide solvent-resistant and impermeable floor.

Hints on storage assembly

Do not store together with foodstuffs.

Storage class according to TRGS 510

Storage class according to TRGS 510 8B Non-combustible corrosive hazardous substances

Further information on storage conditions

Keep under lock and key or accessible only to specialists or people who are authorized.

neodisher Septo Fin

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 22.12.2017

Print date: 22.12.17

7.3. Specific end use(s)

no data

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Exposure limit values

glutaral

List	EH40			
Type	WEL			
Value	0.2	mg/m ³	0.05	ppm(V)
Short term exposure limit	0.2	mg/m ³	0.05	ppm(V)
Status: 2011				

Other information

There are not known any further control parameters.

8.2. Exposure controls

General protective and hygiene measures

Hold emergency shower available. Hold eye wash fountain available. Do not inhale gases/vapours/aerosols. Avoid contact with skin and eyes. Do not eat, drink or smoke during work time. Storage of foodstuffs in work rooms is forbidden. Wash hands before breaks and after work. Clean skin thoroughly after work; apply skin cream.

Respiratory protection

If workplace limits are exceeded, a respiratory protection approved for this particular job must be worn.

Hand protection

Chemical resistant gloves (EN 374)

Use	Permanent hand contact
Appropriate Material	neoprene
Material thickness	>= 0,65 mm
Breakthrough time	> 480
Appropriate Material	nitrile
Material thickness	>= 0,4 mm
Breakthrough time	> 480 min
Appropriate Material	butyl
Material thickness	>= 0,7 mm
Breakthrough time	> 480 min
Use	Short-term hand contact
Appropriate Material	nitrile
Material thickness	>= 0,11 mm

Eye protection

Safety glasses with side protection shield (EN 166)

Body protection

Clothing as usual in the chemical industry. Protective shoes

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Form	liquid
Colour	light yellow
Odour	characteristic
Odour threshold	
Remarks	not determined

neodisher Septo Fin

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 22.12.2017

Print date: 22.12.17

pH value

Value	appr.	4,1	
Temperature		20	°C

Melting point

Remarks	not determined
---------	----------------

Freezing point

Remarks	not determined
---------	----------------

Initial boiling point and boiling range

Remarks	not determined
---------	----------------

Flash point

Remarks	Not applicable
---------	----------------

Evaporation rate (ether = 1) :

Remarks	not determined
---------	----------------

Flammability (solid, gas)

evaluation	not determined
------------	----------------

Upper/lower flammability or explosive limits

Remarks	not determined
---------	----------------

Vapour pressure

Remarks	not determined
---------	----------------

Vapour density

Remarks	not determined
---------	----------------

Density

Value	1,05		g/cm ³
Temperature	20	°C	

Solubility in water

Remarks	miscible in all proportions
---------	-----------------------------

Solubility(ies)

Remarks	not determined
---------	----------------

Partition coefficient: n-octanol/water

Remarks	not determined
---------	----------------

Ignition temperature

Remarks	not determined
---------	----------------

Decomposition temperature

Remarks	not determined
---------	----------------

Viscosity

dynamic

Value	<	50		mPa.s
Temperature		20	°C	

Explosive properties

evaluation	not determined
------------	----------------

Oxidising properties

evaluation	None known
------------	------------

9.2. Other information

Other information

None known

neodisher Septo Fin

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 22.12.2017

Print date: 22.12.17

SECTION 10: Stability and reactivity

10.1. Reactivity

No hazardous reactions when stored and handled according to prescribed instructions.

10.2. Chemical stability

No hazardous reactions known.

10.3. Possibility of hazardous reactions

No hazardous reactions known.

10.4. Conditions to avoid

No hazardous reactions known.

Decomposition temperature

Remarks not determined

10.5. Incompatible materials

None known

10.6. Hazardous decomposition products

Toxic gases/vapours, Irritant gases/vapours

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute oral toxicity

Species	rat			
ATE	300	to	2000	mg/kg
Method	calculated value (Regulation (EC) No. 1272/2008)			

Acute oral toxicity (Components)

glutaral

Species	rat			
	77			mg/kg
Method	OECD 401			

alkylether carboxylic acid

Reference substance	alkylether carboxylic acid			
Species	rat			
LD50	>	2000		mg/kg

Acute dermal toxicity

Remarks Based on available data, the classification criteria are not met.

Acute dermal toxicity (Components)

glutaral

Species	rabbit			
	>	2000		mg/kg

Acute inhalational toxicity

Species	rat			
ATE	1	to	5	mg/l
Administration/Form	Dust/Mist			
Method	calculated value (Regulation (EC) No. 1272/2008)			

Acute inhalative toxicity (Components)

glutaral

Species	rat			
LC50	0,28	to	0,48	mg/l

neodisher Septo Fin

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 22.12.2017

Print date: 22.12.17

Duration of exposure 4 h
Administration/Form Dust/Mist
Method OECD 403

Skin corrosion/irritation

evaluation corrosive

Serious eye damage/irritation

evaluation irritant - risk of serious damage to eyes

Sensitization

Route of exposure dermal
evaluation sensitizing
Route of exposure inhalative
evaluation sensitizing

Subacute, subchronic, chronic toxicity

Remarks Based on available data, the classification criteria are not met.

Mutagenicity

Remarks Based on available data, the classification criteria are not met.

Reproductive toxicity

Remarks Based on available data, the classification criteria are not met.

Carcinogenicity

Remarks Based on available data, the classification criteria are not met.

Specific Target Organ Toxicity (STOT)

Remarks Inhalation may lead to irritation of the respiratory tract.

Aspiration hazard

No special hazards have to be mentioned.

Experience in practice

Inhalation may lead to irritation of the respiratory tract.

Other information

There is no data available on the product apart from the information given in this subsection.

SECTION 12: Ecological information

12.1. Toxicity

General information

not determined

Fish toxicity (Components)

glutaral

Species Fathead minnow (*Pimephales promelas*)
LC50 5,4 mg/l
Duration of exposure 96 h

alkylether carboxylic acid

Reference substance alkylether carboxylic acid
LC50 > 100 mg/l
Duration of exposure 96 h

Daphnia toxicity (Components)

glutaral

Species Daphnia magna
EC50 5 mg/l
Duration of exposure 48 h

neodisher Septo Fin

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 22.12.2017

Print date: 22.12.17

Algae toxicity (Components)

glutaral

Species	Selenastrum capricornutum	
LC50	0,81	mg/l
Duration of exposure	120	h

Bacteria toxicity (Components)

glutaral

Species	activated sludge	
EC20	appr. 15	mg/l
Duration of exposure	30	min
Method	OECD 209	

12.2. Persistence and degradability

General information

not determined

12.3. Bioaccumulative potential

General information

not determined

Partition coefficient: n-octanol/water

Remarks not determined

12.4. Mobility in soil

General information

not determined

12.5. Results of PBT and vPvB assessment

Evaluation of persistence and bioaccumulation potential

The product contains no PBT or vPvB substances.

12.6. Other adverse effects

General information

not determined

General information / ecology

Do not allow to enter soil, waterways or waste water canal. Avoid release into the atmosphere.

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Disposal recommendations for the product

Allocation of a waste code number, according to the European Waste Catalogue (EWC), should be carried out in agreement with the regional waste disposal company.

Disposal recommendations for packaging

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

SECTION 14: Transport information

	Land transport ADR/RID	Marine transport IMDG/GGVSee	Air transport ICAO/IATA
14.1. UN number	1903	1903	1903

neodisher Septo Fin

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 22.12.2017

Print date: 22.12.17

14.2. UN proper shipping name	DISINFECTANT, LIQUID, CORROSIVE, N.O.S. (glutaral)	DISINFECTANT, LIQUID, CORROSIVE, N.O.S. (glutaral)	DISINFECTANT, LIQUID, CORROSIVE, N.O.S. (glutaral)
14.3. Transport hazard class(es)	8	8	8
Label			
14.4. Packing group	III	III	III
Limited Quantity	5 l		
Transport category	3		
14.5. Environmental hazards		no	
Tunnel restriction code	E		
IMDG-Code segregation group		0 Not applicable	

Information for all modes of transport

14.6. Special precautions for user

See Sections 6 to 8

Other information

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

Water Hazard Class (Germany)

Water Hazard Class (Germany) WGK 3

Remarks Classification according to AwSV

VOC

VOC (EU) 0 %

Other information

The product does not contain substances of very high concern (SVHC).

15.2. Chemical safety assessment

For this preparation a chemical safety assessment has not been carried out.

SECTION 16: Other information

Hazard statements listed in Chapter 3

H301 Toxic if swallowed.
H302 Harmful if swallowed.
H314 Causes severe skin burns and eye damage.
H315 Causes skin irritation.
H317 May cause an allergic skin reaction.

neodisher Septo Fin

Version: 3 / GB

Replaces Version: 2 /
GB

Date revised: 22.12.2017

Print date: 22.12.17

H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H330	Fatal if inhaled.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.
H400	Very toxic to aquatic life.
H411	Toxic to aquatic life with long lasting effects.
H412	Harmful to aquatic life with long lasting effects.

CLP categories listed in Chapter 3

Acute Tox. 2	Acute toxicity, Category 2
Acute Tox. 3	Acute toxicity, Category 3
Acute Tox. 4	Acute toxicity, Category 4
Aquatic Acute 1	Hazardous to the aquatic environment, acute, Category 1
Aquatic Chronic 2	Hazardous to the aquatic environment, chronic, Category 2
Aquatic Chronic 3	Hazardous to the aquatic environment, chronic, Category 3
Eye Dam. 1	Serious eye damage, Category 1
Eye Irrit. 2	Eye irritation, Category 2
Resp. Sens. 1	Respiratory sensitization, Category 1
Skin Corr. 1B	Skin corrosion, Category 1B
Skin Irrit. 2	Skin irritation, Category 2
Skin Sens. 1A	Skin sensitization, Category 1A
STOT SE 3	Specific target organ toxicity - single exposure, Category 3

Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: ***
This information is based on our present state of knowledge. However, it should not constitute a
guarantee for any specific product properties and shall not establish a legally valid relationship.