

Safety Data Sheet

according to regulation (EU) No. 1907/2006

Printing date: 04/29/2026

pico-bimsclean

Page 1 of 8

1. Identification of the substance/preparation and of the company/undertaking

1.1. Information on the product

Trade name: pico-bimsclean

1.2. Relevant identified uses of the substance/preparation and uses that are inadvisable

No further relevant information available

1.3. Details about the supplier of the safety data sheet

Manufacturer/Supplier:

Company name: picodent Dental Produktions- und Vertriebs-GmbH
Street: Lüdenscheider Str. 24-26
City: D-51688 Wipperfürth
Telephone: +49 2267 6580-0
E-Mail: picodent@picodent.de
Internet: www.picodent.de
Department for information: picodent Dental Produktions- und Vertriebs-GmbH
Telephone-No. +49 2267 6580-0

1.4. Emergency contact number:

Medical Emergency information in case of poisoning:
Poison Information Center Mainz

-24 h – Phone: +49 6131 19240
advisory service in German or English language

2. Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008:

The product is not classified according to the CLP regulation.

2.2. Label elements

Labelling according to Regulation (EC) No 1272/2008:

Void

Hazard pictograms:

Void

Signal word:

Void

Hazard statements:

Void

Additional information:

Safety data sheet available on request.

2.3. Other hazards:

Results of PBT and vPvB assessment

PBT: Not applicable.
vPvB: Not applicable.

Safety Data Sheet

according to regulation (EU) No. 1907/2006

Printing date: 04/29/2026

pico-bimsclean

Page 2 of 8

3. Composition/information on ingredients *

3.1. Chemical characterization: Mixtures
Description: Mixture of substances listed below with nonhazardous additions.

Dangerous components:
CAS: 57-55-6 Propylene glycol 10-25%
EINECS: 200-338-0 substance with a Community workplace exposure limit
RTECS: TY 2000000
Reg.nr.: 01-2119456809-23-xxxx

Additional information: For the wording of the listed risk phrases refer to section 16.

4. First aid measures *

4.1. Description of first aid measures
After inhalation: Supply fresh air; consult doctor in case of complaints.
After skin contact: Generally the product does not irritate the skin.
After eye contact: Rinse opened eye for several minutes under running water.
After swallowing: If symptoms persist consult doctor.

4.2. Most important symptoms and effects, both acute and delayed: No further relevant information available.

4.3. Indication of any immediate medical attention and special treatment needed: No further relevant information available.

5. Fire-fighting measures

5.1. Extinguishing media
Suitable extinguishing agents: Use fire extinguishing methods suitable to surrounding conditions.

5.2. Special hazards arising from the substance or mixture: No further relevant information available.

5.3. Advice for firefighters
Protective equipment: No special measures required.

6. Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures: Not required.

6.2. Environmental precautions: Dilute with plenty of water. Do not allow to enter sewers/ surface or ground water.

6.3. Methods and material for containment and cleaning up: Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).

Safety Data Sheet

according to regulation (EU) No. 1907/2006

Printing date: 04/29/2026

pico-bimsclean

Page 3 of 8

6.4. Reference to other sections:

No dangerous substances are released.
See Section 7 for information on safe handling.
See Section 8 for information on personal protection equipment.
See Section 13 for disposal information.

7. Handling and storage

7.1. Handling

Precautions for safe handling:

No special precautions are necessary if used correctly.

Information about fire - and explosion protection:

No special measures required.

Conditions for safe storage, including any incompatibilities

7.2. Storage

Requirements to be met by storerooms and receptacles:

Store only in the original receptacle.

Information about storage in one common storage facility:

Not required.

Further information about storage conditions:

Store in upright position.

7.3. Specific end use(s):

No further relevant information available.

8. Exposure controls/personal protection *

8.1. Control parameters

Ingredients with limit values that require monitoring at the workplace:

57-55-6 Propylene glycol

WEL Long-term value: 474* 10** mg/m³, 150* ppm*total vapour and particulates **particulates

Additional information:

The lists valid during the making were used as basis.

8.2. Exposure controls

Appropriate engineering controls:

No further data; see item 7.

Individual protection measures, such as personal protective equipment

General protective and hygienic measures:

Wash hands before breaks and at the end of work.

Respiratory protection:

Not required.

Protection of hands:

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.
Due to missing tests no recommendation to the glove material can be given for the product/ the preparation/ the chemical mixture.
Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Safety Data Sheet

according to regulation (EU) No. 1907/2006

Printing date: 04/29/2026

pico-bimsclean

Page 4 of 8

Material of gloves.

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material cannot be calculated in advance and has therefore to be checked prior to the application.

Penetration time of glove material:

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

For the permanent contact in work areas without heightened risk of injury (e.g. Laboratory) gloves made of the following material are suitable:

Rubber gloves

For the permanent contact gloves made of the following materials are suitable:

Neoprene gloves

As protection from splashes gloves made of the following materials are suitable:

Nitrile rubber, NBR

Not suitable are gloves made of the following materials:

Butyl rubber, BR
Nitrile rubber, NBR
Natural rubber, NR

Eye/face protection:

Goggles recommended during refilling.

9. Physical and chemical properties *

9.1. Information on basic physical and chemical properties

General Information

Physical state:	Fluid
Colour:	According to product specification
Odour:	Characteristic
Odour threshold:	Not determined.
Melting point/freezing point:	Undetermined.
Boiling point or initial boiling point and boiling range:	100 °C
Flammability:	Not applicable.
Lower and upper explosion limit	
Lower:	2.6 Vol % (1,2-propylen-glycol)
Upper:	12.6 Vol % (1,2-propylen-glycol)
Flash point:	101 °C
Ignition temperature:	371 °C
Decomposition temperature:	Not determined.
pH at 20 °C:	6.5
Viscosity:	
Kinematic viscosity at 20 °C:	42 s (ISO 3 mm)
Dynamic:	Not determined.
Solubility water:	Fully miscible.
Partition coefficient n-octanol/water (log value):	Not determined.
Vapour pressure at 20 °C:	23 hPa
Density and/or relative density	

Safety Data Sheet

according to regulation (EU) No. 1907/2006

Printing date: 04/29/2026

pico-bimsclean

Page 5 of 8

Density at 20 °C: 1.02 g/cm³
Relative density: Not determined.
Vapour density: Not determined.

9.2. Other information

Appearance:

Form: Fluid

Important information on protection of health and environment, and on safety

Auto-ignition temperature:

Product is not selfigniting.

Explosive properties:

Product does not present an explosion hazard.

Solvent content:

Organic solvents: 15.0 %

Water: 84.0 %

Solids content: 0.0 %

Change in condition

Evaporation rate: Not determined.

Information with regard to physical hazard classes

Explosives: Void

Flammable gases: Void

Aerosols: Void

Oxidising gases: Void

Gases under pressure: Void

Flammable liquids: Void

Flammable solids: Void

Self-reactive substances and mixtures: Void

Pyrophoric liquids: Void

Pyrophoric solids: Void

Self-heating substances and mixtures: Void

Substances and mixtures, which emit flammable

gases in contact with water: Void

Oxidising liquids: Void

Oxidising solids: Void

Organic peroxides: Void

Corrosive to metals: Void

Desensitised explosives: Void

10. Stability and reactivity

10.1. Reactivity.

No further relevant information available.

10.2. Chemical stability

Thermal decomposition / conditions to be avoided:

No decomposition if used according to specifications.

10.3. Possibility of hazardous reactions:

No dangerous reactions known.

10.4. Conditions to avoid:

No further relevant information available.

10.5. Incompatible materials:

No further relevant information available.

10.6. Hazardous decomposition products:

No dangerous decomposition products known.

Safety Data Sheet

according to regulation (EU) No. 1907/2006

Printing date: 04/29/2026

pico-bimsclean

Page 6 of 8

11. Toxicological information *

11.1. Information on toxicological effects

Acute toxicity:

LD/LC50 values relevant for classification:

ATE (Acute Toxicity Estimates)

Oral LD50 13333 mg/kg (rat)

Skin corrosion/irritation: Based on available data, the classification criteria are not met.

Serious eye damage/irritation: Based on available data, the classification criteria are not met.

Respiratory or skin sensitisation: Based on available data, the classification criteria are not met.

CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction):

Germ cell mutagenicity: Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met.

Reproductive toxicity: Based on available data, the classification criteria are not met.

STOT-single exposure: Based on available data, the classification criteria are not met.

STOT-repeated exposure: Based on available data, the classification criteria are not met.

Aspiration hazard: Based on available data, the classification criteria are not met.

11.2. Information on other hazards

Endocrine disrupting properties: None of the ingredients is listed.

12. Ecological information *

12.1. Toxicity

Aquatic toxicity: No further relevant information available.

12.2. Persistence and degradability: No further relevant information available.

12.3. Bioaccumulative potential: No further relevant information available.

12.4. Mobility in soil: No further relevant information available.

12.5. Results of PBT and vPvB assessment

PBT: Not applicable.

vPvB: Not applicable.

12.6. Endocrine disrupting properties: For information on endocrine disrupting properties see section 11.

12.7. Other adverse effects: No further relevant information available.

Additional ecological information:

Safety Data Sheet

according to regulation (EU) No. 1907/2006

Printing date: 04/29/2026

pico-bimsclean

Page 7 of 8

General notes:

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

13. Disposal considerations *

13.1. Waste treatment methods

Recommendation:

Must not be disposed together with household garbage. Do not allow product to reach sewage system.

Uncleaned packaging

Recommendation:

Disposal must be made according to official regulations.

Recommended cleansing agents:

Water, if necessary together with cleansing agents.

14. Transport information *

14.1. UN-Number:

ADR, ADN, IMDG, IATA:

Void

14.2. UN proper shipping name:

ADR, ADN, IMDG, IATA:

Void

14.3. Transport hazard class(es):

ADR, IMDG, IATA

Class:

Void

Label:

-

ADN/R Class:

Void

14.4. Packing group:

ADR, IMDG, IATA:

Void

14.5. Environmental hazards:

Marine pollutant:

No

14.6. Special precautions for user

Hazard identification number (Kemler code):

-

EMS Number:

-

14.7. Maritime transport in bulk according to IMO instruments:

Not applicable.

Transport/Additional information:

Not dangerous according to the above specifications.

15. Regulatory information *

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

Labelling according to Regulation (EC)

No 1272/2008:

GHS label elements

Hazard pictograms:

Void

Signal word:

Void

Safety Data Sheet

according to regulation (EU) No. 1907/2006

Printing date: 04/29/2026

pico-bimsclean

Page 8 of 8

Hazard statements: Void

15.2. Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

16. Other information *

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

Recommended restriction of use: Product only for professional use

Abbreviations and acronyms

RID:	Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)
IATA-DGR:	Dangerous Goods Regulations by the „International Air Transport Association“ (IATA)
ICAO:	International Civil Aviation Organisation
ICAO-TI:	Technical Instructions by the „International Civil Aviation Organisation“ (ICAO)
ADR:	Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
IMDG:	International Maritime Code for Dangerous Goods
IATA:	International Air Transport Association
GHS:	Globally Harmonized System of Classification and Labelling of Chemicals
EINECS:	European Inventory of Existing Commercial Chemical Substances
ELINCS:	European List of Notified Chemical Substances
CAS:	Chemical Abstracts Service (division of the American Chemical Society)
LC50:	Lethal concentration, 50 percent
LD50:	Lethal dose, 50 percent
PBT:	Persistent, Bioaccumulative and Toxic
vPvB:	very Persistent and very Bioaccumulative

* Data compared to the previous version Rev.-No. 1.4 altered.