

# weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

weigoman

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### Identified Uses

PC29

Pharmaceuticals

### 1.3. Details of the supplier of the safety data sheet

#### Address:

Chemische Fabrik Dr. Weigert GmbH & Co. KG

Mühlenhagen 85

D-20539 Hamburg

Telephone no.

+49 40 789 60 0

Fax no.

+49 40 789 60 120

www.drweigert.com

#### E-mail address of person responsible for this SDS:

sida@drweigert.de

### 1.4. Emergency telephone number

Emergency telephone number: 112

## SECTION 2: Hazards identification \*\*\*

### 2.1. Classification of the substance or mixture

#### Classification (Regulation (EC) No. 1272/2008)

Classification (Regulation (EC) No. 1272/2008)

Flam. Liq. 2 H225

Eye Irrit. 2 H319

STOT SE 3 H336

\*

\*

\*

The product is classified and labelled in accordance with Regulation (EC) No 1272/2008  
For explanation of abbreviations see section 16.

### 2.2. Label elements

#### Labelling according to regulation (EC) No 1272/2008

#### Hazard pictograms



#### Signal word

Danger

# weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

## Hazard statements

H225 Highly flammable liquid and vapour.  
H319 Causes serious eye irritation.  
H336 May cause drowsiness or dizziness.

## Precautionary statements \*\*\*

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
P261 Avoid breathing dust/fume/gas/mist/vapours/spray.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.  
P304+P340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P310 Immediately call a POISON CENTER or doctor.

## Hazardous component(s) to be indicated on label (Regulation (EC) No. 1272/2008)

contains \*\*\* propan-1-ol; propan-2-ol

## 2.3. Other hazards

No special hazards have to be mentioned.  
The product contains no PBT substances. The product contains no vPvB substances. This product does not contain a substance that has endocrine disrupting properties with respect to human. The product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms.

## SECTION 3: Composition/information on ingredients \*\*\*

### 3.2. Mixtures

#### Hazardous ingredients \*\*\*

##### propan-2-ol

|                                                |                  |  |      |   |
|------------------------------------------------|------------------|--|------|---|
| CAS No.                                        | 67-63-0          |  |      |   |
| EINECS no.                                     | 200-661-7        |  |      |   |
| Registration no.                               | 01-2119457558-25 |  |      |   |
| Concentration                                  | >= 50            |  |      | % |
| Classification (Regulation (EC) No. 1272/2008) |                  |  |      |   |
|                                                | Flam. Liq. 2     |  | H225 |   |
|                                                | Eye Irrit. 2     |  | H319 |   |
|                                                | STOT SE 3        |  | H336 |   |

##### propan-1-ol

|                                                |                  |      |      |   |
|------------------------------------------------|------------------|------|------|---|
| CAS No.                                        | 71-23-8          |      |      |   |
| EINECS no.                                     | 200-746-9        |      |      |   |
| Registration no.                               | 01-2119486761-29 |      |      |   |
| Concentration                                  | >= 10            | < 25 |      | % |
| Classification (Regulation (EC) No. 1272/2008) |                  |      |      |   |
|                                                | Flam. Liq. 2     |      | H225 |   |
|                                                | Eye Dam. 1       |      | H318 |   |
|                                                | STOT SE 3        |      | H336 |   |

#### Further ingredients

##### glycerol

|               |           |      |  |   |
|---------------|-----------|------|--|---|
| CAS No.       | 56-81-5   |      |  |   |
| EINECS no.    | 200-289-5 |      |  |   |
| Concentration | >= 1      | < 10 |  | % |
| Advice: [3]   |           |      |  |   |

# weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

## Note

[3] Substance with occupational exposure limits

## Other information

Complete text of hazard statements in chapter 16

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

#### General information

Remove contaminated clothing immediately and dispose of safely. In case of persistent symptoms consult doctor.

#### After inhalation

Ensure supply of fresh air. Seek medical advice immediately.

#### After skin contact

Wash off immediately with soap and water. Seek medical advice immediately.

#### After eye contact

Separate eyelids, wash the eyes thoroughly with water (15 min.). Seek medical advice immediately.

#### After ingestion

Call in a physician immediately and show him the Safety Data Sheet.

#### Adhere to personal protective measures when giving first aid

First aider: Pay attention to self-protection!

### 4.2. Most important symptoms and effects, both acute and delayed

Until now no symptoms known so far.

### 4.3. Indication of any immediate medical attention and special treatment needed

#### Hints for the physician / hazards

In the case of swallowing with subsequent vomiting, aspiration of the lungs can occur which can lead to chemical pneumonia or asphyxiation.

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

#### Suitable extinguishing media

Dry powder, Foam, Carbon dioxide, Water spray jet

#### Non suitable extinguishing media

Full water jet

### 5.2. Special hazards arising from the substance or mixture

In case of combustion evolution of dangerous gases possible.

### 5.3. Advice for firefighters

#### Special protective equipment for fire-fighting

Do not inhale explosion and/or combustion gases. In case of combustion use a suitable breathing apparatus.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

Keep away sources of ignition. Ensure adequate ventilation. Use personal protective clothing.

### 6.2. Environmental precautions

# weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

Prevent spread over a wide area (e.g. by containment or oil barriers). Do not discharge into the drains/surface waters/groundwater. Do not discharge into the subsoil/soil.

## 6.3. Methods and material for containment and cleaning up

Pick up rest with suitable absorbent materials. Containers in which spilt substance has been collected must be adequately labelled. Dispose of as prescribed.

## 6.4. Reference to other sections

Refer to protective measures listed in Sections 7 and 8.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

#### Advice on safe handling

Provide good ventilation of working area (local exhaust ventilation if necessary). Avoid formation of aerosols. Keep container tightly closed.

#### Advice on protection against fire and explosion

Keep away from sources of heat and ignition. No smoking. Take action to prevent static discharges.

### 7.2. Conditions for safe storage, including any incompatibilities

#### Recommended storage temperature

|       |   |    |    |
|-------|---|----|----|
| Value | < | 25 | °C |
|-------|---|----|----|

#### Requirements for storage rooms and vessels

Store product in closed containers.

#### Hints on storage assembly

Do not store with strong oxidizing agents.

#### Storage classes

|                                     |   |                  |
|-------------------------------------|---|------------------|
| Storage class according to TRGS 510 | 3 | Flammable liquid |
|-------------------------------------|---|------------------|

#### Further information on storage conditions

Keep container tightly closed. Keep container in a cool, well-ventilated place. Protect from heat and direct sunlight.

### 7.3. Specific end use(s)

no data

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

#### Exposure limit values

##### propan-1-ol

|                                                    |      |                   |     |        |
|----------------------------------------------------|------|-------------------|-----|--------|
| List                                               | EH40 |                   |     |        |
| Type                                               | WEL  |                   |     |        |
| Value                                              | 500  | mg/m <sup>3</sup> | 200 | ppm(V) |
| Short term exposure limit                          | 625  | mg/m <sup>3</sup> | 250 | ppm(V) |
| Skin resorption / sensibilisation: Sk; Remarks: Sk |      |                   |     |        |

##### propan-2-ol

|                           |      |                   |     |        |
|---------------------------|------|-------------------|-----|--------|
| List                      | EH40 |                   |     |        |
| Type                      | WEL  |                   |     |        |
| Value                     | 999  | mg/m <sup>3</sup> | 400 | ppm(V) |
| Short term exposure limit | 1250 | mg/m <sup>3</sup> | 500 | ppm(V) |

##### Glycerol

|      |      |
|------|------|
| List | EH40 |
|------|------|

# weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

|       |     |                   |
|-------|-----|-------------------|
| Type  | WEL |                   |
| Value | 10  | mg/m <sup>3</sup> |

## Other information

There are not known any further control parameters.

## 8.2. Exposure controls

### General protective and hygiene measures

Do not smoke during work time. Do not inhale gases/vapours/aerosols. Do not eat or drink during work time. Take off immediately all contaminated clothing. Wash hands before breaks and after work.

### Respiratory protection

If workplace limits are exceeded, a respiratory protection approved for this particular job must be worn.

### Hand protection

Chemical resistant gloves

|                                          |                         |
|------------------------------------------|-------------------------|
| Use                                      | Permanent hand contact  |
| Appropriate Material                     | butyl                   |
| Material thickness                       | >= 0,7 mm               |
| Breakthrough time                        | > 480 min               |
| Appropriate Material                     | nitrile                 |
| Material thickness                       | > 0,4 mm                |
| Breakthrough time                        | > 120 min               |
| Use                                      | Short-term hand contact |
| Appropriate Material                     | nitrile                 |
| Material thickness                       | > 0,28 mm               |
| Hand protection must comply with EN 374. |                         |

### Eye protection

Safety glasses with side protection shield; Eye protection must comply with EN 166.

### Body protection

Clothing as usual in the chemical industry.

## SECTION 9: Physical and chemical properties

### 9.1. Information on basic physical and chemical properties

|                |                |
|----------------|----------------|
| Physical state | liquid         |
| Colour         | colourless     |
| Odour          | characteristic |

### Melting point

|         |                |
|---------|----------------|
| Remarks | not determined |
|---------|----------------|

### Freezing point

|         |                |
|---------|----------------|
| Remarks | not determined |
|---------|----------------|

### Boiling point or initial boiling point and boiling range

|       |    |    |
|-------|----|----|
| Value | 83 | °C |
|-------|----|----|

### Flammability

|            |                |
|------------|----------------|
| evaluation | not determined |
|------------|----------------|

### Upper and lower explosive limits

|                       |      |      |
|-----------------------|------|------|
| Lower explosion limit | 2,0  | %(V) |
| Upper explosion limit | 13,5 | %(V) |

### Flash point

|        |           |    |
|--------|-----------|----|
| Value  | 21        | °C |
| Method | DIN 51755 |    |

### Ignition temperature

## weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

Remarks not determined

### Decomposition temperature

Remarks  
Remarks not determined

### pH value

Value appr. 6

### Viscosity

Remarks not determined

### Solubility(ies)

Remarks not determined

### Partition coefficient n-octanol/water (log value)

Remarks not determined

### Vapour pressure

Value 43 hPa  
Temperature 20 °C

### Density and/or relative density

Value 0,85 g/cm<sup>3</sup>  
Temperature 20 °C

### Relative vapour density

Remarks not determined

## 9.2. Other information

### Odour threshold

Remarks not determined

### Evaporation rate (ether = 1) :

Remarks not determined

### Solubility in water

Remarks miscible in all proportions

### Explosive properties

evaluation not determined

### Oxidising properties

evaluation None known

### Other information

None known

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

No hazardous reactions when stored and handled according to prescribed instructions.

### 10.2. Chemical stability

No hazardous reactions known.

### 10.3. Possibility of hazardous reactions

No hazardous reactions known.

### 10.4. Conditions to avoid

No hazardous reactions known. Evolution of flammable mixtures possible in air when heated above flash point and/or during spraying or misting.

### 10.5. Incompatible materials

# weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

Reactions with strong oxidising agents.

## 10.6. Hazardous decomposition products

No hazardous decomposition products known.

## SECTION 11: Toxicological information

### 11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

#### Acute oral toxicity

|         |                                                                   |      |       |
|---------|-------------------------------------------------------------------|------|-------|
| ATE     | >                                                                 | 2000 | mg/kg |
| Remarks | Based on available data, the classification criteria are not met. |      |       |

#### Acute oral toxicity (Components)

##### propan-1-ol

|         |       |      |       |
|---------|-------|------|-------|
| Species | rat   |      |       |
| LD50    | appr. | 5400 | mg/kg |

##### propan-2-ol

|         |          |      |       |
|---------|----------|------|-------|
| Species | rat      |      |       |
| LD50    |          | 5840 | mg/kg |
| Method  | OECD 401 |      |       |

#### Acute dermal toxicity

|         |                                                                   |
|---------|-------------------------------------------------------------------|
| Remarks | Based on available data, the classification criteria are not met. |
|---------|-------------------------------------------------------------------|

#### Acute dermal toxicity (Components)

##### propan-1-ol

|                      |          |      |       |
|----------------------|----------|------|-------|
| Species              | rabbit   |      |       |
| LD50                 |          | 4032 | mg/kg |
| Duration of exposure | 24       | h    |       |
| Method               | OECD 402 |      |       |

##### propan-2-ol

|         |          |       |       |
|---------|----------|-------|-------|
| Species | rabbit   |       |       |
| LD50    |          | 13900 | mg/kg |
| Method  | OECD 402 |       |       |

#### Acute inhalational toxicity

|         |                                                                   |
|---------|-------------------------------------------------------------------|
| Remarks | Based on available data, the classification criteria are not met. |
|---------|-------------------------------------------------------------------|

#### Acute inhalative toxicity (Components)

##### propan-1-ol

|                      |          |      |      |
|----------------------|----------|------|------|
| Species              | rat      |      |      |
| LC50                 | >        | 33,8 | mg/l |
| Duration of exposure | 4        | h    |      |
| Administration/Form  | Vapors   |      |      |
| Method               | OECD 403 |      |      |

##### propan-2-ol

|                      |          |    |      |
|----------------------|----------|----|------|
| Species              | rat      |    |      |
| LC50                 | >        | 25 | mg/l |
| Duration of exposure | 6        | h  |      |
| Administration/Form  | Vapors   |    |      |
| Method               | OECD 403 |    |      |

#### Skin corrosion/irritation

|         |                                                                   |
|---------|-------------------------------------------------------------------|
| Remarks | Based on available data, the classification criteria are not met. |
|---------|-------------------------------------------------------------------|

#### Serious eye damage/irritation

|            |          |
|------------|----------|
| evaluation | irritant |
| Method     | OECD 405 |

# weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

Remarks  
Remarks

Test conducted with a similar formulation.  
The classification criteria are met.

## Sensitization

Remarks

Based on available data, the classification criteria are not met.

## Subacute, subchronic, chronic toxicity

Remarks

Based on available data, the classification criteria are not met.

## Mutagenicity

Remarks

Based on available data, the classification criteria are not met.

## Reproductive toxicity

Remarks

Based on available data, the classification criteria are not met.

## Carcinogenicity

Remarks

Based on available data, the classification criteria are not met.

## Specific Target Organ Toxicity (STOT)

Single exposure  
evaluation

May cause drowsiness or dizziness.

Remarks

The classification criteria are met.

Repeated exposure

Remarks

Based on available data, the classification criteria are not met.

## Aspiration hazard

Based on available data, the classification criteria are not met.

## 11.2 Information on other hazards

### Endocrine disrupting properties with respect to humans

The product does not contain a substance that has endocrine disrupting properties with respect to humans.

### Other information

There is no data available on the product apart from the information given in this subsection.

## SECTION 12: Ecological information

### 12.1. Toxicity

#### General information

not determined

#### Fish toxicity (Components)

##### propan-1-ol

Species

Fathead minnow (*Pimephales promelas*)

LC50

4555

mg/l

Duration of exposure

96

h

Method

OECD 203

##### propan-2-ol

Species

Fathead minnow (*Pimephales promelas*)

LC50

9640

mg/l

Duration of exposure

96

h

#### Daphnia toxicity (Components)

##### propan-1-ol

Species

Daphnia magna

EC50

3644

mg/l

Duration of exposure

48

h

Method

DIN 38412 / Part 11

# weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

## propan-2-ol

|                      |               |      |
|----------------------|---------------|------|
| Species              | Daphnia magna |      |
| LC50                 | appr. 10000   | mg/l |
| Duration of exposure | 48 h          |      |

## Algae toxicity (Components)

### propan-1-ol

|                      |                           |      |
|----------------------|---------------------------|------|
| Species              | Selenastrum capricornutum |      |
| EC50                 | 9170                      | mg/l |
| Duration of exposure | 48 h                      |      |

### propan-2-ol

|                      |                         |      |
|----------------------|-------------------------|------|
| Species              | Scenedesmus subspicatus |      |
| IC50                 | > 1000                  | mg/l |
| Duration of exposure | 72 h                    |      |

## Bacteria toxicity (Components)

### propan-1-ol

|                      |                  |      |
|----------------------|------------------|------|
| Species              | activated sludge |      |
| IC50                 | > 1000           | mg/l |
| Duration of exposure | 3 h              |      |
| Method               | OECD 209         |      |

### propan-2-ol

|         |                  |      |
|---------|------------------|------|
| Species | activated sludge |      |
| EC50    | > 100            | mg/l |

## 12.2. Persistence and degradability

### General information

There is no data available on the product apart from the information given in this subsection.

## 12.3. Bioaccumulative potential

### General information

There is no data available on the product apart from the information given in this subsection.

### Partition coefficient n-octanol/water (log value)

Remarks not determined

## 12.4. Mobility in soil

### General information

For this subsection there is no ecotoxicological data available on the product as such.

## 12.5. Results of PBT and vPvB assessment

### General information

Not applicable

### Results of PBT and vPvB assessment

The product contains no PBT substances  
The product contains no vPvB substances.

## 12.6 Endocrine disrupting properties

### Endocrine disrupting properties with respect to the environment

The product does not contain a substance that has endocrine disrupting properties with respect to non-target organisms.

## 12.7. Other adverse effects

### General information

not determined

# weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

## General information / ecology

Do not discharge product unmonitored into the environment.

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

#### Disposal recommendations for the product

EWC waste code 18 01 06\* chemicals consisting of or containing dangerous substances

EWC waste code 20 01 29\* detergents containing dangerous substances

The listed waste code numbers, according to the European Waste Catalogue (EWC), are to be understood as a recommendation. A final decision must be made in agreement with the regional waste disposal company.

#### Disposal recommendations for packaging

EWC waste code 15 01 02 plastic packaging

Completely emptied packagings can be given for recycling.

EWC waste code 15 01 10\* packaging containing residues of or contaminated by dangerous substances

Packaging that cannot be cleaned should be disposed off in agreement with the regional waste disposal company.

## SECTION 14: Transport information

|                                  | Land transport ADR/RID                                                              | Marine transport IMDG/GGVSee                                                         | Air transport ICAO/IATA                                                               |
|----------------------------------|-------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| Tunnel restriction code          | D/E                                                                                 |                                                                                      |                                                                                       |
| IMDG-Code segregation group      |                                                                                     | 0 Not applicable                                                                     |                                                                                       |
| 14.1. UN number or ID number     | 1987                                                                                | 1987                                                                                 | 1987                                                                                  |
| 14.2. UN proper shipping name    | ALCOHOLS, N.O.S. (propan-2-ol, propan-1-ol)                                         | ALCOHOLS, N.O.S. (propan-2-ol, propan-1-ol)                                          | ALCOHOLS, N.O.S. (propan-2-ol, propan-1-ol)                                           |
| 14.3. Transport hazard class(es) | 3                                                                                   | 3                                                                                    | 3                                                                                     |
| Label                            |  |  |  |
| 14.4. Packing group              | II                                                                                  | II                                                                                   | II                                                                                    |
| Special provision                | 640D                                                                                |                                                                                      |                                                                                       |
| Limited Quantity                 | 1 I                                                                                 | 1 I                                                                                  |                                                                                       |
| Transport category               | 2                                                                                   |                                                                                      |                                                                                       |
| 14.5. Environmental hazards      |                                                                                     | no                                                                                   |                                                                                       |

## Information for all modes of transport

### 14.6. Special precautions for user

# weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

See Sections 6 to 8

## Other information

### 14.7 Maritime transport in bulk according to IMO instruments

Not applicable

## SECTION 15: Regulatory information \*\*\*

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### Major-accident categories acc. 2012/18/EU \*\*\*

|          |     |                  |           |    |            |    |
|----------|-----|------------------|-----------|----|------------|----|
| Category | P5c | FLAMMABLE LIQUID | 5.000.000 | kg | 50.000.000 | kg |
|----------|-----|------------------|-----------|----|------------|----|

#### VOC \*\*\*

|          |       |   |       |     |
|----------|-------|---|-------|-----|
| VOC (EU) | 77,44 | % | 658,2 | g/l |
|----------|-------|---|-------|-----|

### 15.2. Chemical safety assessment

For this preparation a chemical safety assessment has not been carried out.

## SECTION 16: Other information

### Classification and procedure used to derive the classification for mixtures according to Regulation (EC) 1272/2008 [CLP]:

|                                                |      |
|------------------------------------------------|------|
| Classification (Regulation (EC) No. 1272/2008) |      |
| Flam. Liq. 2                                   | H225 |
| Eye Irrit. 2                                   | H319 |
| STOT SE 3                                      | H336 |

### Hazard statements listed in Chapter 2/3

|      |                                     |
|------|-------------------------------------|
| H225 | Highly flammable liquid and vapour. |
| H318 | Causes serious eye damage.          |
| H319 | Causes serious eye irritation.      |
| H336 | May cause drowsiness or dizziness.  |

### CLP categories listed in Chapter 2/3

|              |                                                              |
|--------------|--------------------------------------------------------------|
| Eye Dam. 1   | Serious eye damage, Category 1                               |
| Eye Irrit. 2 | Eye irritation, Category 2                                   |
| Flam. Liq. 2 | Flammable liquid, Category 2                                 |
| STOT SE 3    | Specific target organ toxicity - single exposure, Category 3 |

### Abbreviations

ADR: Accord européen relatif au transport international des marchandises Dangereuses par Route  
RID: Règlement concernant le transport international ferroviaire de marchandises dangereuses  
IMDG: International Maritime Code for Dangerous Goods  
ICAO: International Civil Aviation Organization  
IATA: International Air Transport Association  
MARPOL 73/78: International Convention for the Prevention of Pollution From Ships, 1973 as modified by the Protocol of 1978 (MARPOL: Marine Pollution)  
IBC: Intermediate Bulk Container  
CAS: Chemical Abstracts Service  
VOC: Volatile Organic Compound  
LD: Lethal dose  
LC: Lethal concentration  
PBT: Persistent, Bioaccumulative and Toxic  
vPvB: Very persistent and very bioaccumulative  
SVHC: Substances of very high concern  
EU: European Union

## weigoman

Version: 2 / GB

Replaces Version: 1 / GB

Date revised: 21.12.2022

Print date: 20.01.23

### Supplemental information

Relevant changes compared with the previous version of the safety data sheet are marked with: \*\*\*  
This information is based on our present state of knowledge. However, it should not constitute a  
guarantee for any specific product properties and shall not establish a legally valid relationship.