SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Occlu Spray Plus Spraydose grün
Article number: 554205

1.2 Relevant identified uses of the substance or mixture and uses advised against

1.2.1 Relevant uses

Occlusion spray

1.2.2 Uses advised against

None known.

1.3 Details of the supplier of the safety data sheet

Company

Hager & Werken GmbH & Co. KG
Ackerstr. 1
47269 Duisburg / GERMANY
Phone +49(0)203-99269-0
Fax +49 (0)203 29 92 83
Homepage www.hagerwerken.de
E-mail info@hagerwerken.de

Address enquiries to

Technical information info@hagerwerken.de

Safety Data Sheet sdb@chemiebuero.de (No dispatch of safety data sheets)

Safety data sheets are available from the supplier.

1.4 Emergency telephone number

Advisory body +49 (0) 551-19240 Giftinformationszentrum-Nord

SECTION 2: Hazards identification

2.1 Classification of the substance or mixture [REGULATION (EC) No 1272/2008]

Aerosol 1: H222 Extremely flammable aerosol. H229 Pressurised container: May burst if heated.
Aquatic Chronic 3: H412 Harmful to aquatic life with long lasting effects.

2.2 Label elements

The product is required to be labelled in accordance with regulation CLP.

Hazard pictograms



Signal word DANGER

Hazard statements

H222 Extremely flammable aerosol.
H229 Pressurised container: May burst if heated.
H412 Harmful to aquatic life with long lasting effects.

Precautionary statements

P210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.
P211 Do not spray on an open flame or other ignition source.
P251 Do not pierce or burn, even after use.
P410+P412 Protect from sunlight. Do not expose to temperatures exceeding 50 °C / 122°F.
P273 Avoid release to the environment.
P501 Dispose of contents/container in accordance with local/national regulation.

2.3 Other hazards

Environmental hazards Does not contain any PBT or vPvB substances.
Contains no ingredients with endocrine-disrupting properties.

Other hazards Further hazards were not determined with the current level of knowledge.

SECTION 3: Composition / Information on ingredients

3.1 Substances

not applicable

3.2 Mixtures

The product is a mixture.

Range [%]	Substance
30 - <70	Butane
	CAS: 106-97-8, EINECS/ELINCS: 203-448-7, EU-INDEX: 601-004-00-0, Reg-No.: 01-2119474691-32-XXXX
	GHS/CLP: Flam. Gas 1A: H220 - Press. Gas: H280
25 - <50	Propane
	CAS: 74-98-6, EINECS/ELINCS: 200-827-9, EU-INDEX: 601-003-00-5, Reg-No.: 01-2119486944-21-XXXX
	GHS/CLP: Flam. Gas 1A: H220 - Press. Gas: H280
2,5 - <5	Pentane
	CAS: 109-66-0, EINECS/ELINCS: 203-692-4, EU-INDEX: 601-006-00-1, Reg-No.: 01-2119459286-30-XXXX
	GHS/CLP: Flam. Liq. 1: H224 - Asp. Tox. 1: H304 - STOT SE 3: H336 - Aquatic Chronic 2: H411 - EUH066
1 - <5	iso-Butane
	CAS: 75-28-5, EINECS/ELINCS: 200-857-2, EU-INDEX: 601-004-00-0, Reg-No.: 01-2119485395-27-XXXX
	GHS/CLP: Flam. Gas 1A: H220 - Press. Gas: H280
<1,5	Isopentane
	CAS: 78-78-4, EINECS/ELINCS: 201-142-8, EU-INDEX: 601-006-00-1, Reg-No.: 01-2119475602-38-XXXX
	GHS/CLP: Flam. Liq. 1: H224 - Asp. Tox. 1: H304 - STOT SE 3: H336 - Aquatic Chronic 2: H411 - EUH066

Comment on component parts For full text of H-statements: see SECTION 16.

SECTION 4: First aid measures

4.1 Description of first aid measures

General information Change soaked clothing.

Inhalation Ensure supply of fresh air.
In the event of symptoms seek medical treatment.

Skin contact In case of contact with skin wash off immediately with soap and water.
Consult a doctor if skin irritation persists.

Eye contact In case of contact with eyes rinse thoroughly with plenty of water and seek medical advice.

Ingestion Do not induce vomiting.
In the event of symptoms seek medical treatment.

4.2 Most important symptoms and effects, both acute and delayed

Headache
Vertigo
Nausea, vomiting.
Irritant effects

4.3 Indication of any immediate medical attention and special treatment needed

treat symptomatically

SECTION 5: Fire-fighting measures

5.1 Extinguishing media

Suitable extinguishing media Carbon dioxide.
Dry powder.
Foam.

Extinguishing media that must not be used Water.

5.2 Special hazards arising from the substance or mixture

risk of formation of toxic pyrolysis products, carbon monoxide (CO), not combusted hydrocarbons
Bursting aerosols can be forcibly projected from a fire.

5.3 Advice for firefighters

Use self-contained breathing apparatus.

Fire residues and contaminated firefighting water must be disposed of in accordance within the local regulations.
Cool containers at risk with water spray jet.

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation.
Keep away from all sources of ignition.
Use personal protective equipment (protective gloves, safety glasses, protective clothing).

6.2 Environmental precautions

Do not discharge into the drains/surface waters/groundwater.

6.3 Methods and material for containment and cleaning up

Pick up with absorbent material (e.g. sand, sawdust, universal absorbent, diatomaceous earth).
Dispose of absorbed material in accordance within the regulations.

6.4 Reference to other sections

See SECTION 8+13

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Use only in well-ventilated areas.
Avoid spilling or spraying in enclosed areas.

Keep away from all sources of ignition - Refrain from smoking.
Take precautionary measures against static discharges.

Do not eat, drink, smoke or take drugs at work.
Wash hands before breaks and after work.
Use barrier skin cream.

7.2 Conditions for safe storage, including any incompatibilities

Keep only in original container.

Do not store together with oxidizing agents.

Keep in a cool place, heat causes increase in pressure and risk of bursting.
Protect from heat/overheating and from sun.

7.3 Specific end use(s)

See product use, SECTION 1.2

SECTION 8: Exposure controls / personal protection

8.1 Control parameters

Ingredients with occupational exposure limits to be monitored EU (2004/37/EG)

Substance / EC LIMIT VALUES
Pentane
CAS: 109-66-0, EINECS/ELINCS: 203-692-4, EU-INDEX: 601-006-00-1, Reg-No.: 01-2119459286-30-XXXX
Eight hours: 1000 ppm, 3000 mg/m³
Isopentane
CAS: 78-78-4, EINECS/ELINCS: 201-142-8, EU-INDEX: 601-006-00-1, Reg-No.: 01-2119475602-38-XXXX
Eight hours: 1000 ppm, 3000 mg/m³

DNEL

Substance
Butane, CAS: 106-97-8
There are no DNEL values established for the substance.
Propane, CAS: 74-98-6
There are no DNEL values established for the substance.
Pentane, CAS: 109-66-0
Industrial, inhalative, Long-term - systemic effects, 3000 mg/m³
Industrial, dermal, Long-term - systemic effects, 432 mg/kg bw/day
general population, dermal, Long-term - systemic effects, 214 mg/kg bw/day
general population, oral, Long-term - systemic effects, 214 mg/kg bw/day
general population, inhalative, Long-term - systemic effects, 643 mg/m³
iso-Butane, CAS: 75-28-5
There are no DNEL values established for the substance.
Isopentane, CAS: 78-78-4
Industrial, inhalative, Long-term - systemic effects, 3000 mg/m³
Industrial, dermal, Long-term - systemic effects, 432 mg/kg bw/day
general population, inhalative, Long-term - systemic effects, 643 mg/m³
general population, dermal, Long-term - systemic effects, 214 mg/kg bw/day
general population, oral, Long-term - systemic effects, 214 mg/kg bw/day

PNEC

Substance
Butane, CAS: 106-97-8
There are no PNEC values established for the substance.
Propane, CAS: 74-98-6
There are no PNEC values established for the substance.
iso-Butane, CAS: 75-28-5
There are no PNEC values established for the substance.
Isopentane, CAS: 78-78-4
There are no PNEC values established for the substance.

8.2 Exposure controls

Additional advice on system design	Ensure adequate ventilation on workstation. Measurement methods for taking workplace measurements must meet the performance requirements of DIN EN 482. For example, recommendations are given in the IFA's list of hazardous substances.
Eye protection	Safety glasses. (EN 166:2001)
Hand protection	0,7 mm Butyl rubber, >480 min (EN 374-1/-2/-3). The details concerned are recommendations. Please contact the glove supplier for further information.
Skin protection	Solvent-resistant protective clothing (EN 340)
Other	Personal protective equipment should be selected specifically for the working place, depending on concentration and quantity handled. The resistance of this equipment to chemicals should be ascertained with the respective supplier. Avoid contact with eyes and skin. Do not inhale aerosols.
Respiratory protection	In the event of occupational exposure limits being exceeded or of inadequate ventilation: wear appropriate respiratory protection. Short term: filter apparatus, filter AX (DIN EN 14387).
Thermal hazards	No information available.
Delimitation and monitoring of the environmental exposition	not determined

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

Physical state	solid
Form	aerosol
Color	green
Odor	characteristic
Odour threshold	not determined
pH-value	7 - 10
pH-value [1%]	not applicable
Boiling point or initial boiling point and boiling range [°C]	not applicable
Flash point [°C]	not determined
Flammability	not applicable
Lower explosion limit	not determined
Upper explosion limit	not determined
Oxidising properties	no
Vapour pressure/gas pressure [kPa]	270
Density [g/cm³]	0,79 - 0,89 (Liquid)
Relative density	not determined
Bulk density [kg/m³]	not applicable
Solubility in water	partially miscible
Solubility other solvents	No information available.
Partition coefficient n-octanol/water (log value)	not determined
Kinematic viscosity	not determined
Relative vapour density	not determined
Melting point [°C]	not determined
Auto-ignition temperature [°C]	not determined
Decomposition temperature [°C]	not applicable
Particle characteristics	not applicable

9.2 Other information

none

SECTION 10: Stability and reactivity

10.1 Reactivity

No dangerous reactions known if used as directed.

10.2 Chemical stability

Stable under normal ambient conditions (ambient temperature).

10.3 Possibility of hazardous reactions

Heat causes increase in pressure and risk of bursting.

10.4 Conditions to avoid

Strong heating.

10.5 Incompatible materials

Strong oxidizing agent.

10.6 Hazardous decomposition products

No hazardous decomposition products known.

SECTION 11: Toxicological information

11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

Acute oral toxicity Based on available data, the classification criteria are not met.

Substance
Pentane, CAS: 109-66-0
LD50, oral, Rat, 2000 - 5000 mg/kg (IUCLID)
Isopentane, CAS: 78-78-4
LD50, oral, Rat, >2000 mg/kg bw (OECD 401)

Acute dermal toxicity Based on available data, the classification criteria are not met.

Acute inhalational toxicity Based on available data, the classification criteria are not met.

Substance
Butane, CAS: 106-97-8
LC50, inhalative, Rat, 658 mg/L (IUCLID)
Propane, CAS: 74-98-6
LC50, inhalative, Rat, > 1443 mg/l (15 min) (Lit.)
Pentane, CAS: 109-66-0
LC50, inhalative, Rat, 25,3 mg/l (4h) (IUCLID)
iso-Butane, CAS: 75-28-5
LC50, inhalative, Rat, 1443 mg/L
Isopentane, CAS: 78-78-4
LC50, inhalativ (gas), Rat, > 25,3 mg/L, 4h

Serious eye damage/irritation Based on available data, the classification criteria are not met.

Substance
Butane, CAS: 106-97-8
Eye, non-irritating
Propane, CAS: 74-98-6
Eye, non-irritating
Pentane, CAS: 109-66-0
Eye, non-irritating
iso-Butane, CAS: 75-28-5
Eye, non-irritating
Isopentane, CAS: 78-78-4
no adverse effect observed
Eye, no adverse effect observed

Skin corrosion/irritation Based on available data, the classification criteria are not met.

Substance
Butane, CAS: 106-97-8
dermal, non-irritating
Propane, CAS: 74-98-6
dermal, non-irritating
Pentane, CAS: 109-66-0
dermal, non-irritating
iso-Butane, CAS: 75-28-5
dermal, non-irritating
Isopentane, CAS: 78-78-4
no adverse effect observed

Respiratory or skin sensitisation Based on available data, the classification criteria are not met.

Substance
Butane, CAS: 106-97-8

dermal, non-sensitizing
inhalative, non-sensitizing
Propane, CAS: 74-98-6
dermal, non-sensitizing
inhalative, non-sensitizing
Pentane, CAS: 109-66-0
dermal, non-sensitizing
iso-Butane, CAS: 75-28-5
dermal, non-sensitizing
inhalative, non-sensitizing
Isopentane, CAS: 78-78-4
dermal, no adverse effect observed

Specific target organ toxicity — single exposure Based on available data, the classification criteria are not met.

Substance
Butane, CAS: 106-97-8
inhalative, non-irritating
Propane, CAS: 74-98-6
inhalative, non-irritating
iso-Butane, CAS: 75-28-5
inhalative, non-irritating

Specific target organ toxicity — repeated exposure Based on available data, the classification criteria are not met.

Substance
Pentane, CAS: 109-66-0
NOAEC, inhalative, Rat, 20000 mg/m³, no adverse effect observed
Isopentane, CAS: 78-78-4
NOAEC, inhalative, Rat, 20000 mg/m³, no adverse effect observed

Mutagenicity Does not contain a relevant substance that meets the classification criteria.

Substance
Pentane, CAS: 109-66-0
in vitro, negativ
in vivo, negativ
Isopentane, CAS: 78-78-4
in vitro, negativ
in vivo, negativ

Reproduction toxicity Does not contain a relevant substance that meets the classification criteria.

- Fertility

Substance
Pentane, CAS: 109-66-0
NOAEL, oral, Rat, 1000 mg/kg bw/day, no adverse effect observed
NOAEC, inhalative, Rat, 20000 mg/m³, no adverse effect observed
Isopentane, CAS: 78-78-4
NOAEC, inhalative, Rat, 24080 mg/m³, no adverse effect observed

- Development

Substance
Pentane, CAS: 109-66-0
NOAEL, oral, Rat, 1000 mg/kg bw/day, no adverse effect observed
NOAEC, inhalative, Rat, 24080 mg/m³, no adverse effect observed

Isopentane, CAS: 78-78-4
NOAEL, oral, Rat, 1000 mg/kg bw/day, no adverse effect observed

Carcinogenicity Does not contain a relevant substance that meets the classification criteria.

Aspiration hazard Based on available data, the classification criteria are not met.

General remarks

Toxicological data of complete product are not available.
The toxicity data listed pertaining to the ingredients are intended for those working in the medicinal professions, experts for occupational health and safety and toxicologists. The toxicity data pertaining to the ingredients were supplied by the manufacturers of raw materials.

11.2 Information on other hazards

11.2.1 Endocrine disrupting properties Does not contain a relevant substance that meets the classification criteria.

11.2.2 Other information

SECTION 12: Ecological information

12.1 Toxicity

Substance
Butane, CAS: 106-97-8
LC50, (48h), Invertebrates, 14,22 - 69,43 mg/L
Pentane, CAS: 109-66-0
EL50, (48h), Daphnia magna, 48,11 mg/l
EL50, (72h), Algae, 20,33 mg/l
LL50, (96h), Fish, 27,55 mg/l
Isopentane, CAS: 78-78-4
EL50, (48h), Daphnia magna, 59,9 mg/L
EL50, (72h), Pseudokirchneriella subcapitata, 25,3 mg/L
LL50, (96h), Oncorhynchus mykiss, 34,3 mg/L
EL10, (21d), Daphnia magna, 11,5 mg/L
EL10, (60d), Oncorhynchus mykiss, 6,57 mg/L

12.2 Persistence and degradability

Behaviour in environment compartments not determined

Behaviour in sewage plant not determined

Biological degradability not determined

Substance
Isopentane, CAS: 78-78-4
The product is readily biodegradable.

12.3 Bioaccumulative potential

Accumulation in organisms is not expected.

12.4 Mobility in soil

No information available.

12.5 Results of PBT and vPvB assessment

Based on all available information not to be classified as PBT or vPvB respectively.

12.6 Endocrine disrupting properties

Does not contain a relevant substance that meets the classification criteria.

12.7 Other adverse effects

Ecological data of complete product are not available.

SECTION 13: Disposal considerations

13.1 Waste treatment methods

Waste material must be disposed of in accordance with the Directive on waste 2008/98/EC as well as other national and local regulations. It is not possible to determine a waste code for this product in accordance with the European Waste Catalogue (EWC) since it is only possible to classify it according to how it is used by the customer. The waste code is to be determined within the EU in liaison with the waste-disposal operator.

Product

Dispose of as hazardous waste.

Waste no. (recommended) 160504* gases in pressure containers (including halons) containing dangerous substances

Contaminated packaging

Uncontaminated packaging may be taken for recycling.

Waste no. (recommended) 150110* packaging containing residues of or contaminated by hazardous substances

SECTION 14: Transport information

14.1 UN number or ID number

Transport by land according to ADR/RID	1950
Inland navigation (ADN)	1950
Marine transport in accordance with IMDG	1950
Air transport in accordance with IATA	1950

14.2 UN proper shipping name

Transport by land according to ADR/RID	Aerosols
- Classification Code	5F
- Label	
- ADR LQ	1 l
- ADR 1.1.3.6 (8.6)	Transport category (tunnel restriction code) 2 (D)
Inland navigation (ADN)	Aerosols
- Classification Code	5F
- Label	
Marine transport in accordance with IMDG	Aerosols
- EMS	F-D, S-U
- Label	
- IMDG LQ	1 l
Air transport in accordance with IATA	Aerosols, flammable
- Label	

14.3 Transport hazard class(es)

Transport by land according to ADR/RID	2
Inland navigation (ADN)	2
Marine transport in accordance with IMDG	2.1
Air transport in accordance with IATA	2.1

14.4 Packing group

Transport by land according to ADR/RID	not applicable
Inland navigation (ADN)	not applicable
Marine transport in accordance with IMDG	not applicable
Air transport in accordance with IATA	not applicable

14.5 Environmental hazards

Transport by land according to ADR/RID	no
Inland navigation (ADN)	no
Marine transport in accordance with IMDG	no
Air transport in accordance with IATA	no

14.6 Special precautions for user

Relevant information under SECTION 6 to 8.

14.7 Maritime transport in bulk according to IMO instruments

not applicable

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

EEC-REGULATIONS	2008/98/EG (2000/532/EC); 2010/75/EU; 2004/42/EG; (EG) 648/2004; (EC) 1907/2006 (REACH); (EU) 1272/2008; 75/324/EWG ((EC) 2016/2037); (EU) 2020/878; (EU) 2016/131; (EU) 2024/573; (EU) 2019/1148; (EU) 2019/1021, (EU) 2023/707
- Comment on component parts	Substances of Very High Concern - SVHC: substances are not contained or are below 0.1%.
- Annex XIV (REACH)	According to Annex XIV of Regulation (EC) 1907/2006 (REACH) the product does not contain any substances ≥ 0.1% that are subject to authorisation.
- Annex XVII (REACH)	According to Annex XVII of Regulation (EC) 1907/2006 (REACH) the product contains ≥ 0.1% of substances with the following restrictions. 40 According to Annex XVII of Regulation (EC) 1907/2006 (REACH) the product is not subject to any restrictions.
TRANSPORT-REGULATIONS	ADR (2025); IMDG-Code (2025, 42. Amdt.); IATA-DGR (2025)
NATIONAL REGULATIONS (EU):	
- Observe employment restrictions for people	Observe employment restrictions for young people.
- VOC (2010/75/CE)	99 %

15.2 Chemical safety assessment

No information available.

SECTION 16: Other information

16.1 Hazard statements (SECTION 3)

EUH066 Repeated exposure may cause skin dryness or cracking.
H411 Toxic to aquatic life with long lasting effects.
H336 May cause drowsiness or dizziness.
H304 May be fatal if swallowed and enters airways.
H224 Extremely flammable liquid and vapour.
H280 Contains gas under pressure; may explode if heated.
H220 Extremely flammable gas.

16.2 Abbreviations and acronyms:

ADR = Accord européen relatif au transport international des marchandises Dangereuses par Route
RID = Règlement concernant le transport international ferroviaire de marchandises dangereuses
ADN = Accord européen relatif au transport international des marchandises dangereuses par voie de navigation intérieure
ATE = acute toxicity estimate
CAS = Chemical Abstracts Service
CLP = Classification, Labelling and Packaging
DMEL = Derived Minimum Effect Level
DNEL = Derived No Effect Level
EC50 = Median effective concentration
ECB = European Chemicals Bureau
EEC = European Economic Community
EINECS = European Inventory of Existing Commercial Chemical Substances
EL50 = Median effective loading
ELINCS = European List of Notified Chemical Substances
EmS = Emergency Schedules
GHS = Globally Harmonized System of Classification and Labelling of Chemicals
IATA = International Air Transport Association
IBC-Code = International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk
IC50 = Inhibition concentration, 50%
IMDG = International Maritime Code for Dangerous Goods
IUCLID = International Uniform ChemicaL Information Database
IVIS = In vitro irritation score
LC50 = Lethal concentration, 50%
LD50 = Median lethal dose
LC0 = lethal concentration, 0%
LOAEL = lowest-observed-adverse-effect level
LL50 = Median lethal loading
LQ = Limited Quantities
MARPOL = International Convention for the Prevention of Marine Pollution from Ships
NOAEL = No Observed Adverse Effect Level
NOEC = No Observed Effect Concentration
PBT = Persistent, Bioaccumulative and Toxic substance
PNEC = Predicted No-Effect Concentration
REACH = Registration, Evaluation, Authorisation and Restriction of Chemicals
STP = Sewage Treatment Plant
TLV®/TWA = Threshold limit value – time-weighted average
TLV®STEL = Threshold limit value – short-time exposure limit
VOC = Volatile Organic Compounds
vPvB = very Persistent and very Bioaccumulative

16.3 Other information

This document does not comply with Regulation (EC) No 1907/2006, article 31 (5) and may be used for internal purposes only.

Classification procedure

Aerosol 1: H222 Extremely flammable aerosol. (Bridging principle "Aerosols") H229 Pressurised container: May burst if heated. (Bridging principle "Aerosols")
Aquatic Chronic 3: H412 Harmful to aquatic life with long lasting effects. ()

Modified position

1.1, 1.3, 3.2, 8.1, 9.1, 11.1, 11.2, 12.6, 12.7, 15.1, 16.2, 16.3

Copyright: Chemiebüro®